



Retail Packaged Aluminum Crutches with Accessories

HCPCS: E0114,E0114

- Lighter than wood, the standard aluminum crutch provides stability and durability
- Double extruded center tube provides additional strength to weight-bearing area
- Pushpin adjustment makes it easy-to-adjust crutch height in 1" increments
- Non-skid, jumbo-size vinyl contoured tips provide excellent traction
- Easy wing nut adjustments for hand grips
- Comfortable, durable accessories are latex free
- Weight Capacity 350 lbs on aluminum crutch (excluding pediatric)

PRODUCTS

Description

RTL10400	Crutch, Adult,1PR,8/CS
RTL10402	Crutch,Tall,1PR,8/CS

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Weight Capacity	350 lbs
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]