

Universal Aluminum Crutch with Accessories

HCPCS: E0114



FIGURE A



FIGURE B



FIGURE C




FIGURE D



FIGURE E

- Universal adjustment for use by a child, adult or tall adult (Figure A)
- Euro-Style clips and pushpin adjustment make it easy-to-adjust height of hand grip (Figure B)
- Lighter than wood, the aluminum crutch provides stability and durability
- Crutch disassembles for convenient storage (Figure C)
- Non-skid, jumbo-size vinyl contoured tips provide excellent traction
- Comfortable, durable underarm pad and hand grip

PRODUCTS

Description

RTL10433	Universal Quick Adjust Crutch,1PR,8/CS
----------	--

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Weight Capacity	300 lb
Height	38.5" - 62.5"