

Premium Plastic Raised, Regular/Elongated Toilet Seat, with Lock

HCPCS: E0244



FIGURE A



FIGURE B

- Designed for individuals who have difficulty sitting down or standing up from the toilet
- Wide opening in front and back for personal hygiene
- Lightweight and portable
- Standard locking mechanism ensures safety (Figure A)
- Fits most elongated toilets
- No tools required for installation
- Easy-to-clean
- Made using recycled materials

PRODUCTS

Description

12014	Raised Toilet St. w/o Arms,1/cs, RTL
-------	--------------------------------------

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Seat Dimensions	17" (W) x 16.5" (D) x 4.5" (H)
Weight Capacity	300 lb
Carton Shipping Weight	6 lb
Weight	4 lb