

## Scout Compact 3-Wheel Scooter

HCPCS: K0800-1,K0800



FIGURE A



FIGURE B



FIGURE C



FIGURE D




FIGURE E



FIGURE F

The Scout transportable scooter combines easy and convenient mobility with superior value. The Scout is compact, lightweight and easy to disassemble-perfect for taking wherever you want to go, either around the neighborhood or far from home. It can be disassembled quickly and easily into five easy-to-handle pieces for loading into a car, truck, or van and taking shopping or on vacation. The three-wheel version of the Scout allows for nimbly negotiating indoor obstacles. The Scout also offers great features like a quick connect battery pack and a delta-style tiller to make operating the scooter easier for those with limited strength or dexterity. This model comes with color-through panels to hide scuffs and scratches and black wheel rims.

- Quick and easy disassembly
- Interchangeable color panels hide scuffs and scratches
- Armrests are padded and adjustable
- Quick connect batteries
- Ergonomic throttle control
- Height-adjustable swivel seat with fold-down backrest
- Flat-free, non-marking tires
- Large, plastic carry basket
- Easy-to-adjust tiller
- Anti-tip wheels
- Delta tiller
- Easy freewheel operation

| PRODUCTS           | Description                              |
|--------------------|------------------------------------------|
| SFSCOUT3           | Scout 3-Wheel Scooter, Red/Blue, 1/ea    |
| SPECIFICATIONS     | Description                              |
| Warranty           | 24 Month Limited Warranty on Electronics |
| Product Weight     | 85 lb                                    |
| Product Dimensions | 42.5" (L) x 20.5" (W) x 34" (H)          |
| Seat Dimensions    | 16" (D) x 16.5" (W)                      |
| Weight Capacity    | 300 lb                                   |
| Top Speed          | 4.25 mph                                 |
| Batteries (pair)   | 12V x 12AH                               |
| Battery Charger    | 2A Offboard                              |
| Climbing Angle     | 6°                                       |
| Color              | Red & Blue                               |
| Controller         | Dynamic 45A                              |
| Ground Clearance   | 2.5"                                     |

| SPECIFICATIONS        | Description                        |
|-----------------------|------------------------------------|
| HCPCS                 | K0800                              |
| Maximum Range         | 9 miles                            |
| Motor                 | 24V x 270W x 4700 RPM              |
| Seat Type             | Stadium                            |
| Turning Radius        | 45.5"                              |
| Warranty on Batteries | 12 Month Warranty on Batteries     |
| Warranty on Frame     | Limited Lifetime Warranty on Frame |
| Wheel Size            | Front: 8" x 2"   Rear: 8" x 2"     |
| ACCESSORIES           | Description                        |
| SF8030                | Crutch Holder Set, SType, 1/ea, EA |