

Silver Sport 1 Wheelchair

HCPCS: K0001-1,K0001-1



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Easy to clean powder coated silver vein steel frame
- Padded fixed armrests
- Nylon upholstery is durable, lightweight, attractive and easy-to-clean
- polyurethane tires mounted on composite wheels provide durability and low maintenance
- Comes with push-to-lock wheel locks
- Carry pocket on backrest provides additional convenience
- Comes with swing-away footrests
- Folding back-canes for transport and storage
- Weight Capacity: 250 lb

PRODUCTS	Description
SSP118FA-SF	SlvrSprt-1,18"SlvrVeinFxdArm,SwgFt,1/cs

SPECIFICATIONS	Description
Warranty	Limited 3 Year
Product Dimensions	Open w/ Riggings: 42" (L) x 24" (W) x 36" (H) Folded without Riggings: 32" (L) x 11" (W) x 22" (H)
Seat Dimensions	18" (W) x 16" (D)
Weight Capacity	250 lb
Seat To Floor Height	19.5"
Carton Shipping Weight	42.8 lb
Weight	37.6 lb
Back Height	16"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov]

ACCESSORIES	Description
STDS831	Heel Loops, 1pr, BX
STDS1040S	Cup Holder, Swivel,1BX,6/CS, BX
STDS833N	Rear Anti Tipper w.wheels SSP&CS,1/pr, PR
STDS6008-1	Carry Bag,W/C,Oxygen,1/ea, EA
LK3JELR	ELR Set,Black,1Pair eq EA, PR
STDS818	Anti-Tippers,rear,SentraRecliner,1pr/bx, BX
STDS850	Seat Belt, Clasp Type,1/ea, EA
STDS851	Seat Belt, Velcro Close,1/ea, EA
STDS855	Seat Belt,60"" Auto-type,1/ea, EA

ACCESSORIES	Description
STDS820	Telescoping IV Pole,W/C, 1/bx, BX
STDS856	Seat Belt,60"Velcro-type,1/ea, EA
STDS806	Anti-Fold/Theft Device for w/c,1/ea, BX
STDS834	OverheadAnti- Theft Bar,Single,1/ea, EA
STDS801GT	WheelLockExtension,ViperGT,SSP2,1/ea, EA
WASR	Limb Support, Right,1/cs, CS
WASL	Limb Support,Left,1/cs, CS
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
STDELR-AL	ArticulatingElev.Legrest, Std,1pr/cs, CS
STDSNUT	Tinnerman Nut,1pr/bx, EA
14300	16x17.5 Gen.use adj back cush 16-21in, EA