

Toilet Safety Frame with Padded Arms

HCPCS:



FIGURE A



FIGURE B



FIGURE C

- Anodized aluminum is sturdy and lightweight
- Padded armrests provide additional comfort and support
- Comes with height-adjustable (1" increments) and width-adjustable arms (2 positions) (Figure A)
- Powder-coated aluminum bracket easily attaches frame to toilet (Figure B)
- Bracket secures safety frame in place

PRODUCTS	Description
RTL12000	Toilet Safety Frame,1EA,4/CS
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	300 lb
Carton Shipping Weight	20 lb
Weight	4 lb
Width Between Arms	19" (smallest) 23.50" (largest)