

Raised Toilet Seat with/without Lid

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Designed for individuals who have difficulty sitting down or standing up from the toilet
- Heavy-duty molded plastic construction
- Locking device with larger, heavy-duty "worm screw" and locking plate for a safe, secure fit (Figure A)
- Lightweight and portable
- Fits Standard Toilets
- No tools required for installation
- Easy-to-clean
- Made using recycled materials

PRODUCTS	Description
12062	Raised ToiletSeat,2"w/oLid,1/cs, RTL
12063	Raised ToiletSeat,2in WithLid,1/cs, RTL
12065	RaisedToiletSeat,4"W/Lid, 1/bx, RTL
12066	RaisedToiletSeat,6"w/o Lid, 1/cs
RTL12064	Raised Toilet Seat,4" Without Lid,1 EA

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	14" (W) x 16" (D)
Weight Capacity	400 lb

ACCESSORIES	Description
12062-C	Clamps, 12062 Series, 12L Ser.#, 1PR/BX, BX