

Trimline CPAP Tubing

HCPCS: A7037,A7037




FIGURE A

- Premier 6, trim (15mm diameter) CPAP tubing
- 15 mm diameter significantly reduces noise from PAP machine improving the patient therapy experience
- 15 mm diameter also significantly reduces drag on mask, reducing leaks and improving patient comfort and compliance
- 22mm cuff connections are over-molded onto the tubing
- Retail packaged available

PRODUCTS

Description

CPAPTUB6TL

CPAPTube 6' Grey Trimline,1/ea,12/CV