

Wallaby Pediatric Wheelchair

HCPCS:



FIGURE A



- Standard Transit Tie Down Brackets
- WC19 Compliant when used with WB8021 for users between 50 and 80 lbs
- Height adjustable push handles
- Flip back padded desk arms
- Pelvic belt standard
- Push to lock wheel locks
- 22" Rear wheels with 6" front casters
- Swingaway footrests with tool free height adjustment and calf strap

PRODUCTS	Description
WB1200-2GFR	Wallaby, 12", Fire Truck Red
WB1400-2GJB	Wallaby, 14", Jet Fighter Blue
SPECIFICATIONS	Description
Product Weight	33 lb (without footrests)
Seat Dimensions	12"W x 12"D 14"W x 14"D
Handle Height	33" - 39"
Weight Capacity	150 lb
Seat To Floor Height	19"
Overall Width	20" 22"
Back Height	14"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]
ACCESSORIES	Description
WB 8021	Headrest Extension, Wallaby, 1ea, EA
WB ELR	Elevating Leg Rests, Wallaby, 1pr, PR
STDS833N	Rear Anti Tipper w.wheels SSP&CS,1/pr, PR