

EXCURSION™ ROLLATORS



Two smart choices offering premium design and functionality

The M+D® by Drive® Excursion™ Series delivers premium comfort to help you stand out in today's mobility market. Both models offer smooth rolling, comfortable support, and the easy folding users want, while giving retailers MAP protection, strong display appeal, and fewer assembly challenges.

With two configurations—the Excursion X20 and the higher-capacity Excursion XWD—the series lets providers meet mobility needs for a broad range of users with solutions that look great, perform reliably, and help drive long-term customer satisfaction.

Premium Comfort & Support

Offers comfortable navigation and rest with padded seats and full, contoured backrests and comfortable brake access with ergonomic handgrips.

Smooth, Confident Mobility

Provides a stable, smooth roll indoors and outdoors thanks to wide, 8" treaded wheels positioned outside the handles, plus easy-to-activate hand brakes with parking brake function and curb climber.

Compact Portability

Easy to fold with just the pull of a strap and weighs just 19.4 lbs. (X20) / 20.5 lbs. (XWD) so users can pop it into trunks and storage spaces with minimal effort.



Product Specs

Seat Dimensions
Seat to Floor Height
Handle Height
Weight Capacity
Product Weight
Frame
Wheels
Warranty
Accessories

X20

17.5" W x 10" D
20"
30.5" – 37"
300 lb
19.4 lb
Side-fold
8" wide, treaded
Limited Lifetime
Optional slow brakes item # 730-942

XWD

19.5" W x 12" D
20"
30.5" – 37"
400 lb
20.5 lb
Side-fold
8" wide, treaded
Limited Lifetime

Item #	Description	UPC
MD700920BK	Excursion X20, Onyx	822383030500
MD700920BL	Excursion X20, Sapphire	822383030517
MD700920PR	Excursion X20, Lavender	822383030524
MD700929BL	Excursion XWD, Sapphire	822383030531