

TRIUMPH
escape

**WHERE
COMFORT**
MEETS
COMPACT

Escape Petite



Anatomical handles with
easy-to-apply, hidden
cable braking system

Ergonomic padded
backrest with
adjustable depth

Padded seat

15.5"

A compact folding system
activated by the lift of a
release handle


Strong
aluminum
T-frame

Large 8-inch wheels for
smoother navigation

Width between handles: 15.5"
Seat heights: 19"

Crafted for our petite
clients 5'4" and under

Easy to adjust handle height,
pressure-reducing anatomical
handgrips



ACCESSORIES

Sold separately



Cane Holder



Cup Holder



Phone Holder



Bamboo Tray

SPECIFICATIONS

Colour	Model Numbers
Charcoal	500-15191
Description	
Seat height	19" (48 cm)
Handle height range	28 - 36" (71 - 91.5 cm)
Distance between handles	15.5" (39 cm)
Distance between wheels	15.75" (40 cm)
Seat size	14.5 x 9.5" (47 x 24 cm)
Overall width	22.5" (63.5 cm)
Overall length	25.6" (65 cm)
Depth folded	11" (28 cm)
Wheel diameter	8" (200 cm)
Weight capacity of rollator	300 lb (136 kg)
Weight capacity of shopping bag	11 lb (5 kg)
Rollator Weight	15.4 lb (7 kg)

