

Specifications

Colour	Model Numbers		
Charcoal	500-10191	500-10211	500-10241
Red	500-10195	500-10215	500-10245
Blue	500-10192	500-10212	500-10242
Description	Super Low	Low	Standard
Seat height	19" (48 cm)	21" (54 cm)	24" (61 cm)
Handle height range	28 - 36" (71 - 91.5 cm)	29.5 - 37.5" (75 - 95.3 cm)	32 - 40" (81.3 - 101.6 cm)
Distance between handles	19" (45.7 cm)	19" (45.7 cm)	19" (45.7 cm)
Distance between wheels	19" (48.25 cm)	19" (48.25 cm)	19" (48.25 cm)
Seat size	18.5 x 9.5" (47 x 24 cm)	18.5 x 9.5" (47 x 24 cm)	18.5 x 9.5" (47 x 24 cm)
Overall width	25" (63.5 cm)	25" (63.5 cm)	25" (63.5 cm)
Overall length	25.6" (65 cm)	25.6" (65 cm)	25.6" (65 cm)
Depth folded	11" (28 cm)	11" (28 cm)	11" (28 cm)
Wheel diameter	8" (200 mm)	8" (200 mm)	8" (200 mm)
Weight capacity of rollator	275 lb (125 kg)	275 lb (125 kg)	275 lb (125 kg)
Weight capacity of shopping bag	11 lb (5 kg)	11 lb (5 kg)	11 lb (5 kg)
Rollator weight	15.4 lb (7 kg)	15.6 lb (7.1 kg)	15.8 lb (7.2 kg)




Limited lifetime warranty
triumphmobility.com

Authorized Dealer

Contact us

Triumph Mobility
151 Bentley Ave, Unit 4
Nepean, Ontario, K2E 6T7
Tel: +1 (855) 546 0711
Fax: +1 (866) 225 3059
Email: info@triumphmobility.com
www.triumphmobility.com

[f](#) [@](#) [in](#) [@](#) [d](#) [X](#) **TriumphMobility**



SCAN ME!

©2023 Triumph Mobility Inc. All rights reserved. Triumph™ is the trademark of Triumph Mobility.

Printed in Canada

TRIUMPH
escape





Freedom to go anywhere you want to go

The Escape rollator features adjustable handle height with the pressure-reducing anatomical handgrips to maintain your hand strength and the brakes are the easiest-to-apply hand brakes available on any rollator.

Features

-  Largest handle height range in the market
-  Adjustable padded seat & backrest
-  3 seat heights (19", 21", 24")
3 colour options (Red, Blue, Grey)
-  Strong aluminum frame
-  Gentle touch brake system
-  Spring assisted easy folding
-  Anti-shock 8" wheels
-  Curb climber



Quality Accessories

Sold Separately



Cup Holder



Phone Holder



Tray



Cane Holder
(cane not included)