



Product Catalog

STYLE **FUNCTION** COMFORT

OUR COMPANY

Flight Dental Systems is a Canadian owned and operated dental equipment manufacturer with sales in over 25 countries globally. It is one of the fastest growing dental equipment manufacturers in the industry and its philosophy has always been to produce an innovative product with exceptional design, reliability and quality, at an affordable price. This philosophy is evident throughout the Flight product line, as each item is meticulously designed and engineered to suit the needs and requirements of every practitioner. In addition, Flight Dental Systems continues to improve its products through extensive R&D to provide the best possible product for the best possible price to it's customer.

A6 PATIENT CHAIR

STANDARD FEATURES



Linak Motor
From Denmark



Memory Foam



Hydraulic System
Made In Japan



ERGONOMIC DESIGN

Thin Backrest and Seatrest allows better access to the oral cavity while maintaining patient comfort with the memory foam support. The synchronized compensating backrest provides optimal patient and operator comfort during any procedure.

DESIGNED & ENGINEERED IN CANADA

Designed and Engineered in Canada according to the highest ISO 13485:2016 Standards. Perfect form and function to ensure the operator achieves the best ergonomics to achieve the best results while providing the patient the most comfortable stress-free experience.



FUNCTIONAL DESIGN

Dual Drop Down Armrest allows for easy Patient Access and exceptional support when the patient is in the supine position.



- ✓ A6 Hybrid Patient Chair (Hydraulic and Electromechanical)
- ✓ Hydraulic Base and Electromechanical Back Rest
- ✓ Seat Height 14.6" to 27.4" (370mm to 697mm)
- ✓ Backrest Reclining 4 degree to 67 degree
- ✓ Double Articulating Headrest
- ✓ Electronic Deluxe Foot Control: 2 Programmable Positions,
- ✓ Auto Return and Last Position
- ✓ Input Voltage (110v or 220v)



STOOLS



DELUXE DOCTOR STOOL

Dimensions:

| | |
|------------------------------------|--------------------|
| Backrest Vertical Adjustment Range | 0~ 4.7"(12cm) |
| Seat Width | 17"(43.18cm) |
| Seat Height Ranges | 20"-25.5"(51-68cm) |
| 5 Leg Base Diameter | 24"(61cm) |



DELUXE ASSISTANT STOOL

Dimensions:

| | |
|--------------------------------------|--------------------|
| Abdominal Vertical Adjustment Range: | 0~3.1"(8cm) |
| Seat Width | 18"(45.7cm) |
| Seat Height Ranges | 22"-27.5"(56-70cm) |
| 5 Leg Base Diameter | 24"(61cm) |



SADDLE STOOL

Dimensions:

| | |
|---------------|--------------------------|
| Seat Height | 21.6" to 29.5" (55-75cm) |
| Base Diameter | 18.9" (48cm) |
| Seat Depth | 14" (35.5cm) |
| Seat Width | 17.1" (45.3cm) |



SADDLE STOOL WITH BACKREST

** AVAILABLE IN BLACK ONLY**

Dimensions:

| | |
|-------------------------------------|----------------------|
| Abdominal Vertical Adjustment Range | 0~ 3.1"(8cm) |
| Seat Width | 16.1"(41cm) |
| Seat Height Ranges | 24.5"-31.5"(62-80cm) |
| 5 Leg Base Diameter | 24"(61cm) |



SPLIT-SADDLE STOOL

** AVAILABLE IN BLACK ONLY**

Dimensions:

| | |
|---------------|--------------------------|
| Seat Height | 21.6" to 29.5" (55-75cm) |
| Base Diameter | 18.9" (48cm) |
| Seat Depth | 14" (35.5cm) |
| Seat Width | 17.1" (45.3cm) |

SADDLE STOOL COLOR OPTIONS



UPHOLSTERY STYLES



SEAMLESS



STITCHED

COLOURS

Standard Plush Seamless and Sewn Upholstery



*Custom Colours Available for Both Ultraleather and Plush – Upcharge

FLIGHT DENTAL SYETEMS - First Class has Never been so Affordable!



Flight Dental Systems Canada Head Office
21 Kenview Blvd, Unit 9, Brampton, ON, L6T 5G7, CANADA
Tel: 905-799-0517 Fax: 905-793-2406

Flight Dental Systems USA
3103 Farrell Road, Houston, Texas, 77073, USA
Tel: 346-398-4817

www.flightdentalsystems.com
Toll Free: 1-866-799-0517