

CURING LIGHTS

XLITE 2neo

Features

Cures in 3 Seconds	Large Lens 12 mm
Convenient	Powerful

Technical Specifications

Light Source: High Powered Blue Violet LED

Wavelength Range: 385-515nm

Light Intensity: >2,300mW/cm2

Battery: Lithium Ion Battery
Continuous 120 Minutes of
Running Time on Full Battery

Input: 100-240V AC 50/60Hz

Output: 5V DC

Charging Time: 90 minutes

Large
Lens
12 mm

ULTRASONIC SCALERS

MAGNETOSTRICTIVE SCALERS



FMS-2530

BI-2530K

Features

- Dual Frequency Ultrasonic System: Accepts both 25 and 30K Inserts
- Automatically adjusts to 25 or 30k
- Turbo Function controlled by a 2-position foot switch to provide extra power to remove heavy calculus
- Detachable Autoclavable Handpiece Sleeve
- Individual Water Control
- Power Control

Operates with our 30Khz Inserts

Compatible with Dentsply Cavitron Units



TM0004-192
IF-5 #3 Beavertail 30k



TM0004-462
IF-100 #100 Plastic slim 30k internal



TM0004-362
IF-50 #10 Plastic universal 30k internal



TM0005-072
FS Triple Bend Internal Flow Insert #22BW 30K

BUILT- IN PIEZO SCALER



N6LED-1902

N6LED-1902

- Easy & Compact to Install in Any Delivery Unit
- Piezo Technology
- LED Illumination
- Microcomputer control and Automatic Frequency Tracking
- Compatible with EMS Tips