



Product Catalog

STYLE **FUNCTION** COMFORT

OUR COMPANY

Flight Dental Systems is a Canadian owned and operated dental equipment manufacturer with sales in over 25 countries globally. It is one of the fastest growing dental equipment manufacturers in the industry and its philosophy has always been to produce an innovative product with exceptional design, reliability and quality, at an affordable price. This philosophy is evident throughout the Flight product line, as each item is meticulously designed and engineered to suit the needs and requirements of every practitioner. In addition, Flight Dental Systems continues to improve its products through extensive R&D to provide the best possible product for the best possible price to it's customer.

DELIVERY UNITS

REAR DELIVERY

TRADITIONAL DELIVERY UNIT



TRAD-2002



TRAD-2001

CONTINENTAL DELIVERY UNIT



CONT-3002

PANEL MOUNT DELIVERY



RD-2102/VAC

FEATURES

- ✓ Easy to Install Cabinet Panel Mount Delivery
- ✓ Includes 3 HP Control Block + Syringe
- ✓ Optional Vacuum Package
- ✓ Clean Water Bottle System Included

UNDER CABINET MOUNT DELIVERY

RD-3100
3 HP Under Cabinet Mount



RD-3150
3 HP Under Cabinet Mount with Vacuum (HVE, SE)



AI-105



MOBILE CART DELIVERY

FEATURES

- ✓ Optional Vacuum with HVE + SE (MC-1300FVAC/MC-1302FVAC)
- ✓ Solid Frame with Adjustable Height
- ✓ Dual Casters for all surfaces
- ✓ Water Bottle System Included



MC-1300F



MC-1302F

OPTIONS

- ✓ Water Bottle System (110-032)
- ✓ Fiber Optic System (FO-101KIT)

- ✓ Under Cabinet Mount Vacuum Arm (AI-105)
- ✓ Extra HVE (XHVE-102)

MK DENT HIGHSPEED HANDPIECE



HB23K - MK Highspeed

Single Spray, No Light, Mini Head

HB23KL - MK Dent Highspeed

Single Spray, With Light, Mini Head

HE22KL - MK Dent Highspeed

Triple Spray, with Light, Mini Head

HE21KL - MK Dent Highspeed

Quad Spray, with light, Standard Head



Triple Spray
Quick & Efficient
Cooling

**Fiber Optic
Illumination**
Clear & Bright
Illumination


Push Button Chuck
Easy to Release
and Engage Bur

**German
Ceramic
Bearings**
Quality
Precision


**Quiet and
Comfortable
Operation**
Low Noise
and Powerful

HIGH SPEED HANDPIECES

COUPLERS



SELF-GENERATED LED 45° SURGICAL HANDPIECES



MODEL: F-160

FEATURES

- ✓ 350,000 RPM
- ✓ Single Water Spray
- ✓ Rear Exhaust Air
- ✓ 25,000 Lux LED Light – Self Generating
- ✓ 4 Hole Midwest Connection
- ✓ 21W Torque

SLOW SPEED HANDPIECE COLLECTION



ONE PIECE HANDPIECES



AIR MOTOR



AM0004
MK DENT 20K AIR MOTOR
NO WATER

EXTERNAL WATER - NSK COMPATIBLE



INTERNAL WATER SPRAY - KAVO COMPATIBLE



ELECTRIC MOTORS

LIGHTWEIGHT, POWERFUL AND QUIET

Flight Electric Motors systems helps to improve efficiency, reduce aerosols while providing a more comfortable patient experience. The Flight Electric system combines the state of the art light weight micromotors with a LED system and quiet and vibration free electric solution.



EM-9001
Flight Tabletop Electric Motor,
3.4Ncm Torque, 1 Motor



EMB-9002
Flight Built in Electric Motor System,
3.4Ncm Torque, 1 motor



MX2 LED Micro-Series

Type	Brushless, with internal spray
Max Torque	3.5 N
Speed Range	100-40,000 rpm
Weight	94 gr

Bien Air
Dental

ELECTRIC CONTRA ANGLE ATTACHMENTS

MK-dent
GERMANY



LE11 MK Dent

▼ **LE11 MK Dent**
1:1 Contra Angle,
CA Burs



LE15L MK Dent

▼ **LE15L MK Dent**
1:5 Contra Angle
with Light, FG Burs

PREMIUM OIL LUBRICATION

▼ **LU1011**
MK Dent Premium Synthetic Oil
- Spray Can

▼ **LU1022**
MK Dent Premium Synthetic Oil
- Spray Can Compatible
with Kavo Quattrocare,
Scican Statmatic

▼ **O-3205**
Flight Premium Maintenance
Oil (1000 ml)



LU1011



LU1022



O-3205

ADAPTERS



▼ **PN: LT1014**
FG Spindle



▼ **PN: LT1013**
Slow Speed Attachments



▼ **PN: LT1012**
Fits Kavo High Speed

AUTOMATIC HANDPIECE MAINTENANCE SYSTEM

ACCOMMODATES ALL BRANDS OF HANDPIECES

All brands of E-Type low speed handpieces connect directly to the coupling supplied. The highspeed handpiece connects directly to a 4-hole connection or a quick disconnect coupler can be installed.



Model: HPL-3001

1 Low Speed Connection and
2 Highspeed Connections

Model: HPL-3002

2 Low Speed Connections
and 1 Highspeed Connection

FEATURES

- ✓ 2 High speed Connections
- ✓ 1 Low speed Connection
- ✓ Easy Fill Oil Reservoir
- ✓ Clean and mess-free oil sponge
- ✓ Compressed Air Input with Easy connection
- ✓ Individual programmed Cycles that provides thorough cleaning and lubrication
- ✓ Easily detachable door for easy cleaning



EASY TO CLEAN

The front door can be completely removed to allow fast, easy cleaning of the system



OIL RESERVOIR

The oil reservoir has a 500ml capacity, reducing the need for frequent refilling. The clean easy to use fill tank also has an oil level indicator to allow the user to see when to add more oil.



PROGRAMMED CYCLES

S: Short Cycle for High Speed Handpieces
(Lubricating and Cleaning Time: 35 seconds)

L: Long Cycle for Low Speed Handpieces
(Lubricating and Cleaning Time: 45 seconds)

E: Extra Long Cycle for Speed Reduction
Handpieces (Lubricating and Cleaning Time:
50 seconds)

Technical Specifications

Rated Voltage	AC120V +/-10% 60Hz
Input Voltage	120V: 19VA
Air Pressure	0.35-0.6MPa - (50-85psi)
Dimensions	W 280mm x D275mm x H 360mm W 11" x D10.8" x H 14.2"
Weight	6.8kg or 15lbs



Flight Dental Systems Canada Head Office
21 Kenview Blvd, Unit 9, Brampton, ON, L6T 5G7, CANADA
Tel: 905-799-0517 Fax: 905-793-2406

Flight Dental Systems USA
3103 Farrell Road, Houston, Texas, 77073, USA
Tel: 346-398-4817

www.flightdentalsystems.com
Toll Free: 1-866-799-0517