

GAMMEX™ Non-Latex PI Green

S U R G I C A L G L O V E S

Soft, skin-friendly polyisoprene surgical gloves, for allergy protection and superb comfort

- **Reliable allergy protection:** Made from synthetic polymers, using these sterile polyisoprene (PI) gloves eliminates Type I latex allergy risks while ZMBT, DPG and CPC* are not used to protect wearers from these known allergens
- **Superior comfort and dexterity:** Their soft yet strong formulation ensures wearers benefit from exceptional fit and feel
- **Excellent grip and sensitivity:** A smooth surface with a micro-textured finish offers the wearer secure handling and instrument grip during the most rigorous procedures
- **Effortless donning and doffing:** Their unique internal coating facilitates both dry and damp hand donning, either when first donned or during intra-operative changes
- **Chemo-tested:** Tested to ASTM D6978, these GAMMEX® gloves are safe to use with a wide range of chemotherapy drugs**

**Free from natural rubber latex proteins and do not contain Zinc 2-mercaptobenzothiazole (ZMBT), Diphenylguanidine (DPG) and cetylpyridinium chloride (CPC)*

***Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations*

Recommended For

- All surgeries, especially high risk procedures in orthopaedics, trauma, and cardiac
- Tested for use with chemotherapy drugs



GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.



GAMMEX™ Non-Latex PI Green

SURGICAL GLOVES

PRODUCT INFORMATION

Material	PI (Polyisoprene)
Powder Content	Powder-Free
Color	Green
Thickness	Standard
Cuff Style	Beaded
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Polymer Coated
Allergy Prevention	Latex (Type I)
Grip Level	Moderate
Double Gloving Recommendation	Underglove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	ASTM D3577, ASTM D6978, ASTM D7160, ASTM F1670, ASTM F739
Quality / Environmental Standards	ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Zinc Dialkyldithiocarbamate

ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	20685255	20685260	20685265	20685270	20685275	20685280	20685285	20685290

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 844 494 7854
Email: insidesalesus@ansell.com
Sandel US Tel: 866 764 3327
Sandel email: sandel-cs@ansell.com

Ansell Commercial Mexico
Tel: +52 442 296 2050
Email: latam.medical@ansell.com
Ansell Brasil Ltda
+55 (11) 3356-3100
Email: latam.medical@ansell.com
Ansell Colombia SAS
Tel: +57 1 288 3247
Email: latam.medical@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu
Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

India
Tel: +91 22 38124500
Email: info.india@ansell.com
Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com
China
Tel: +86(21) 3827 5005
Email: infochina@ansell.com

ansell.com

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking.
© 2025 Ansell Limited. All Rights Reserved.