

GAMMEX™ Non-Latex PI

S U R G I C A L G L O V E S

Soft, skin-friendly polyisoprene surgical gloves, for optimal comfort and allergy protection

- **Optimal allergy protection:** Being non-latex, these polyisoprene surgical gloves eliminate the risk of a latex allergy while minimizing sensitization to a Type IV chemical allergy as chemical accelerators known to cause skin allergies have been eliminated from its skin-friendly formulation
- **Superior comfort and dexterity:** Ansell's proprietary soft yet strong formulation delivers exceptional fit and feel
- **Excellent grip and sensitivity:** The smooth surface of these surgical gloves has a micro-textured finish to ensure confident handling and instrument grip during the most rigorous procedures
- **Effortless donning and doffing:** A unique internal coating aids dry and damp hand donning, either when first donned or during intra-operative changes
- **Chemo-tested:** Tested to ASTM D6978, these GAMMEX® sterile gloves are safe to use with a wide range of chemotherapy drugs**



**Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

Recommended For

- All surgeries, especially high risk procedures in orthopaedics, trauma, and cardiac
- Extra comfort during long, sustained procedures
- Protection from chemical accelerators such as thiurams, carbamates and mercaptos known to cause chemical Type IV allergies
- Having passed the Mouse Embryo Assay (MEA) test, may be used for In Vitro Fertilization (IVF) procedures

FEATURED TECHNOLOGIES



GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.



GAMMEX™ Non-Latex PI

SURGICAL GLOVES

PRODUCT INFORMATION

Material	PI (Polyisoprene)
Powder Content	Powder-Free
Color	White
Thickness	Standard
Cuff Style	Beaded with SUREFIT™ Technology
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Polymer Coated with DERMASHIELD™ Technology
Allergy Prevention	Latex (Type I)
Grip Level	Moderate
Double Gloving Recommendation	Outer glove or underglove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	AS/NZS 4179, ASTM D3577, ASTM D6978, ASTM D7160, ASTM F739, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282, JIS T9107
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Xanthates, Zinc Dialkyldithiocarbamate

ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	20685755	20685760	20685765	20685770	20685775	20685780	20685785	20685790

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 844 494 7854
Email: insidesalesus@ansell.com
Sandel US Tel: 866 764 3327
Sandel email: sandel-cs@ansell.com

Ansell Commercial Mexico
Tel: +52 442 296 2050
Email: latam.medical@ansell.com
Ansell Brasil Ltda
+55 (11) 3356-3100
Email: latam.medical@ansell.com
Ansell Colombia SAS
Tel: +57 1 288 3247
Email: latam.medical@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu
Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

India
Tel: +91 22 38124500
Email: info.india@ansell.com
Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com
China
Tel: +86(21) 3827 5005
Email: infochina@ansell.com

ansell.com

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking.
© 2025 Ansell Limited. All Rights Reserved.