

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 1/5/2026

Axiom Back Width	Fits Back Frame Widths	
	Inner Slots	Outer Slots
12" (30.5cm)	11" (28.0cm)	12" - 14" (30.5cm - 35.5cm)
14" (35.5cm)	12" - 13" (30.5cm - 33.0cm)	14" - 16" (35.5cm - 40.5cm)
16" (40.5cm)	14" - 15" (35.5cm - 38.0cm)	16" - 18" (40.5cm - 45.5cm)
18" (45.5cm)	16" - 17" (40.5cm - 43.0cm)	18" - 20" (45.5cm - 50.5cm)
20" (51.0cm)	18" - 19" (45.5cm - 48.5cm)	20" - 22" (51cm - 56cm)
22" (56.0cm)	20" - 21" (51.0cm - 53.5cm)	22" - 24" (56.0cm - 61.0cm)



Axiom Custom Back

Axiom Custom Back (E2617)


\$1499

Complex Axiom Custom Back (E2617)

\$2402

Axiom Custom Back Shape

Hardware is included with back. Selection must be made below. Back selections 12H or less are not available with Accessory Lateral supports, Headrest Mount or Anterior Support mounts. Backs taller than 12H include headrest and Anterior Trunk Support mounting locations. See Axiom Back Sizing Guide for additional standard back dimensions not listed.



Dimensions E, F & G can be found in Axiom Back Sizing Guide

Retail	N/C		N/C	N/C
Item	AXCPB	AXCPB	AXCDL	
Description	2.25" Lateral Depth (H)		3.75" Lateral Depth (H)	6" Lateral Depth (H)

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1077.5 Pricing and specifications subject to change without notice.

Axiom Custom Back Sizes

For sizes not listed see Complex Custom Back Dimensions below. Height measurement is for metal shell only. Finished Back is 2" taller as foam cushion extends 1" over all edges.

		Wheelchair Frame Width							
		12"	14"	16"	18"	20"	22"		
Retail		N/C	N/C	N/C	N/C	N/C	\$163		
Axiom Back Height	Low	12"	N/C	12W X 12H	14W x 12H	16W X 12H	18W x 12H	20W x 12H	
		14"	N/C	12W X 14H	14W X 14H	16W X 14H	18W x 14H	20W x 14H	22W x 14H
	Mid	16"	N/C	12W X 16H	14W X 16H	16W X 16H	18W x 16H	20W x 16H	
		18"	N/C	12W X 18H	14W X 18H	16W X 18H	18W x 18H	20W x 18H	22W x 18H
	Tall	20"	\$163		14W X 20H	16W x 20H	18W x 20H	20W x 20H	
		22"	\$163			16W x 22H	18W x 22H	20W x 22H	22W x 22H

Axiom Complex Custom Back Dimensions


Back Width Range from 12W to 24W available in 1" increments. Back Height Range (D) from 8H to 24H available in 1" increments. Lateral Height (BL/BR) is the dimension from the top edge to bottom edge of the lateral support measured at centerline. Lateral Height minimum height is 6" on PB, PL, & DL. Top of Lateral Location (CL/CR) is measured from the bottom edge of the backrest to the top edge of the lateral support at centerline. Dimensions are for metal shell only. Foam cushion extends 1" over all edges, so finished dimensions will be 1" greater in all directions.

Frame Width ^{1,2}	Height (D) ^{1,2}	Left Lateral Height (BL) ¹	Left Lateral Top Location (CL) ¹	Right Lateral Height (BR) ¹	Right Lateral Top Location (CR) ¹

¹ All inputs required in this section. Inputs must be in 1" increments only
² \$160 Upcharge for for Widths 22" and greater and/or Heights 20" and greater

Mounting Hardware

- DEX Hardware N/C
Removable adjustable hardware. 300lbs Weight Capacity.
- BTX Hardware N/C
Removeable adjustable hardware. Not recommended for or installed on Ultralight Rigid Frame Wheelchairs. Required on 23" or 24" widths. Not intended for use on frames wider than 25". 400lbs Weight Capacity.
- FXT Hardware N/C
Non-removable hardware. Not available with Lateral Supports. 300lbs Weight Capacity.



Axiom Custom Back Foam Selections

1/4" Closed Cell Foam is used as a base level under all level 1 foam. Max foam thickness is 4". Only one foam selection per level. You must select at least one foam level.

Item Number	Description	Level 1	Level 2	Level 3	Level 4	Retail	Notes
BCKFOM700 ^{1,2}	1-1/4" 2.5# Density HR Foam	N/C				\$99	
BCKFOM015	1/2" HR1					\$47	
BCKFOM011	1" HR1					\$80	
BCKFOM025	1/2" HR2					\$47	
BCKFOM021	1" HR2					\$80	
BCKFOM035	1/2" HR3					\$47	
BCKFOM031	1" HR3					\$80	
BCKFOM105	1/2" X-Soft SunMate					\$80	
BCKFOM101	1" X-Soft SunMate					\$136	
BCKFOM115	1/2" Soft SunMate					\$80	
BCKFOM111	1" Soft SunMate					\$136	
BCKFOM125	1/2" Medium SunMate					\$80	
BCKFOM121	1" Medium SunMate					\$136	
BCKFOM205	1/2" X-Soft T-Foam					\$107	
BCKFOM201	1" X-Soft T-Foam					\$184	
BCKFOM215	1/2" Soft T-Foam					\$107	
BCKFOM211	1" Soft T-Foam					\$184	
BCKFOM225	1/2" Medium T-Foam					\$107	
BCKFOM221	1" Medium T-Foam					\$184	

¹ This foam in level 1 only would be equivalent to a standard Axiom Back
² Max 3 levels of 2.5# Density HR foam can be selected

Custom Axiom Back Modifications

Available only on Axiom Custom Back

Item	Description	Retail	HCPCS
BCKBMD501	FIP Package ¹	\$685	E2617
BCKACC601	1 Unit of FIPS Solution	\$409	E2617/K0108
BCKACC701	Cube Cut Foam Modification ²	\$286	E2617

¹ Package includes 4" deep cover, 1" Medium Sunmate over 1/4" Closed cell base. Includes loose piece of 1/4" Soft Sunmate to use between molded foam and face of cover. Foam Chemical kit is included

² Foam cushion insert is cut into 2" x 2" cubes for increased surface tension reduction and pressure relief. Foam is cut to a depth of approximately 1/4" above the closed cell foam bottom base layer. Foam will be cut through all levels of foam selected above. Not available with Foam-in-Place (FIPs) Modification.

Axiom Custom Back Cover Options

Custom Back Cover Material Options ¹	Axiom Fabric Cover	TekStretch Wipeable Cover			
Custom Back Embroidery Options ²	Block	Script	Name		Color:
Military Service Branch Decal ³	Air Force	Army	Navy	Marines	Coast Guard
Extra Cover Material Options ¹	Axiom Fabric Cover	TekStretch Wipeable Cover			
Extra Cover Embroidery Options ²	Block	Script	Name		Color:
Extra Cover Military Service Branch Decal ³	Air Force	Army	Navy	Marines	Coast Guard

¹ Cover material choice is for contact fabric only. (Required selection)

² 10 Character maximum, including spaces. No designs available on Axiom Custom Back. Embroidery not available with 6" Deep Lateral Back (AXCDL)

³ Decal is affixed to rear shell of backrest

Lateral Trunk Supports for AXCPB Backs



Fixed Modular Mounting Hardware

Swing Away Mounting Hardware

Adjustable Link, Swing Away, Fixed Pad Mount Hardware

Adjustable Link, Swing Away, Swivel Pad Mount Hardware

Left	Right	Description	HCPCS	Retail
PBLFM3x4L	PBLFM3x4R	3 x 4 Lateral Pads, Flat with Fixed Modular Mounting Hardware	E0956	\$160
PBLFM4x5L	PBLFM4x5R	4 x 5 Lateral Pads, Flat with Fixed Modular Mounting Hardware	E0956	\$160
PBLS3x4L	PBLS3x4R	3 x 4 Lateral Pads, Flat with S/A Mounting Hardware	E0956, E1034	\$340
PBLS4x5L	PBLS4x5R	4 x 5 Lateral Pads, Flat with S/A Mounting Hardware	E0956, E1034	\$340
PBLSC3x4L	PBLSC3x4R	3 x 4 Lateral Pads, Curved with S/A Mounting Hardware	E0956, E1034	\$340
PBLSC4x5L	PBLSC4x5R	4 x 5 Lateral Pads, Curved with S/A Mounting Hardware	E0956, E1034	\$340
PBLALSF3X4PL	PBLALSF3X4PR	3 x 4 Curved Lateral Pads with Adjustable Link, Swing Away, Fixed Pad Mount Hardware	E0956, E1034	\$453
PBLALSF4X5PL	PBLALSF4X5PR	4 x 5 Curved Lateral Pads with Adjustable Link, Swing Away, Fixed Pad Mount Hardware	E0956, E1034	\$453
PBLALSS3X4PL	PBLALSS3X4PR	3 x 4 Curved Lateral Pads with Adjustable Link, Swing Away, Swivel Pad Mount Hardware	E0956, E1034	\$503
PBLALSS4X5PL	PBLALSS4X5PR	4 x 5 Curved Lateral Pads with Adjustable Link, Swing Away, Swivel Pad Mount Hardware	E0956, E1034	\$503

Pair	Description	HCPCS	Retail
15534P	Omnilink Hardware Cover 5" (Pair)	K0108	\$58

Available only with Adjustable Link Hardware

Lateral Trunk Supports for AXCP Backs



Fixed Modular Mounting Hardware



Swing Away Mounting Hardware



Adjustable Link, Swing Away, Fixed Pad Mount Hardware



Adjustable Link, Swing Away, Swivel Pad Mount Hardware

Left	Right	Description	HCPCS	Retail
PLLFM3x4L	PLLFM3x4R	3 x 4 Lateral Pads, Flat with Fixed Modular Mounting Hardware	E0956	\$160
PLLFM4x5L	PLLFM4x5R	4 x 5 Lateral Pads, Flat with Fixed Modular Mounting Hardware	E0956	\$160
PLLS3x4L	PLLS3x4R	3 x 4 Lateral Pads, Flat with S/A Mounting Hardware	E0956, E1034	\$340
PLLS4x5L	PLLS4x5R	4 x 5 Lateral Pads, Flat with S/A Mounting Hardware	E0956, E1034	\$340
PLLSC3x4L	PLLSC3x4R	3 x 4 Lateral Pads, Curved with S/A Mounting Hardware	E0956, E1034	\$340
PLLSC4x5L	PLLSC4x5R	4 x 5 Lateral Pads, Curved with S/A Mounting Hardware	E0956, E1034	\$340
PLLALSF3X4PL	PLLALSF3X4PR	3 x 4 Curved Lateral Pads with Adjustable Link, Swing Away, Fixed Pad Mount Hardware	E0956, E1034	\$453
PLLALSF4X5PL	PLLALSF4X5PR	4 x 5 Curved Lateral Pads with Adjustable Link, Swing Away, Fixed Pad Mount Hardware	E0956, E1034	\$453
PLLALSS3X4PL	PLLALSS3X4PR	3 x 4 Curved Lateral Pads with Adjustable Link, Swing Away, Swivel Pad Mount Hardware	E0956, E1034	\$503
PLLALSS4X5PL	PLLALSS4X5PR	4 x 5 Curved Lateral Pads with Adjustable Link, Swing Away, Swivel Pad Mount Hardware	E0956, E1034	\$503

Pair	Description	HCPCS	Retail
15534P	Omnilink Hardware Cover 5" (Pair)	K0108	\$58

Available only with Adjustable Link Hardware

Left Lateral Offset Adapter

Available with both AXCPB and AXCP Backs

- 1" Offset Adapter for Fixed Modular and S/A Mounting Hardware Left \$99
- 2" Offset Adapter for Fixed Modular and S/A Mounting Hardware Left \$99

Right Lateral Offset Adapter

Available with both AXCPB and AXCP Backs

- 1" Offset Adapter for Fixed Modular and S/A Mounting Hardware Right \$99
- 2" Offset Adapter for Fixed Modular and S/A Mounting Hardware Right \$99

Painted Back Shell Modificaton

Gloss Finish

If no selection is made, standard Shadow Black will be used. Touch-up paint is available for purchase through the parts manual.

Black	N/C
Ruby Red	N/C
Radiant Orange	N/C
Yellow	N/C
Electric Green	N/C
Rich Teal	N/C
Glacier Blue	N/C
Blue Steel	N/C
Deep Sapphire Blue	N/C
Purple Plum	N/C
Opalescent Violet	N/C
Dark Pink	N/C
Light Pink	N/C
Champagne	N/C

Matte Finish

If no selection is made, standard Shadow Black will be used. Touch-up paint is available for purchase through the parts manual.

Iridescent Navy	N/C
Bright Red	N/C
Tactical Green	N/C
Dry Earth Tan	N/C
Platinum	N/C
White	N/C

Headrests

Headrest Pads

Selection of a Headrest Pad requires a selection from Headrest Mounting Hardware section below



Contour Positioning



Tri-Pad Adjustable



Comfort Soft



TriSoft



Comfort Soft Adjustable

Headrest Pads			
Description	Item Code	HPCPS	Retail
Contour Positioning Headrest Pad, Small, 8W x 6.25H (E0955)	116919	E0955	\$190
Contour Positioning Headrest Pad, Medium, 9W x 6.5H (E0955)	116921	E0955	\$190
Contour Positioning Headrest Pad, Large, 10W x 7.5H (E0955)	116923	E0955	\$190
Tri-Pad Adjustable Headrest, Small (E0955)	116925	E0955	\$460
Tri-Pad Adjustable Headrest, Medium (E0955)	116927	E0955	\$460
Tri-Pad Adjustable Headrest, Large (E0955)	116929	E0955	\$460
Comfort Soft Headrest, 8" (E0955)	116931	E0955	\$204
Comfort Soft Headrest, 10" (E0955)	116933	E0955	\$204
Comfort Soft Headrest, 14" (E0955)	116935	E0955	\$204
Comfort Soft Headrest, 19" (E0955)	116937	E0955	\$235
TriSoft Headrest, 10" (E0955)	116939	E0955	\$204
TriSoft Headrest, 14" (E0955)	116941	E0955	\$204
Comfort Soft Adjustable Headrest, 10" (E0955)	116943	E0955	\$383
Comfort Soft Adjustable Headrest, 14" (E0955)	116945	E0955	\$383
Comfort Soft Adjustable Headrest, 19" (E0955)	116947	E0955	\$477

See Axiom Hardware Sizing Guide for full item sizing

Headrest Hardware

Selection of Headrest Hardware requires a selection from Headrest Pad section above



Std. Straight Removable Adjustable



Multi-Axis Removable Adjustable



Omnilink Multi-Axis Adjustable



Posalinc Removable Adjustable



Posalinc Flip-Back Adjustable

Headrest Hardware			
Description	Item Code	HPCPS	Retail
Standard Straight Removable Adjustable Hardware (E1033)	120038	E1033	\$198
Multi-Axis Removable Adjustable Hardware (E1033)	120039	E1033	\$295
Omnilink Multi-Axis Adjustable Hardware (E1033)	125494	E1033	\$417
POSAlinc Removable Adjustable Hardware (E1033)	120040	E1033	\$359
POSAlinc Removable Adjustable Hardware, Tall (E1033)	123528	E1033	\$396
POSAlinc Flip-Back Adjustable Hardware (E1033)	123530	E1033	\$430
POSAlinc Flip-Back Adjustable Hardware, Tall (E1033)	123531	E1033	\$474

Universal Headrest Mounting Plate (K0108)
 Not required with Axiom Headrests

\$54

Anterior Trunk Supports



2-Piece
Hook & Loop



1-Piece
with Buckle



Dynaform Standard,
Dynamic



Dynaform Slim Cut,
Dynamic



Airlogic Standard,
Structure



Airlogic Slim Cut,
Dynamic

Description	Item Code	HCPCS	Retail
Dynaform Chest Support: 2-Piece with Hook&Loop, 3" High Pad, Small	14941	E0960	\$146
Dynaform Chest Support: 2-Piece with Hook&Loop, 4" High Pad, Medium	14942	E0960	\$146
Dynaform Chest Support: 2-Piece with Hook&Loop, 5" High Pad, Large	14943	E0960	\$146
Dynaform Chest Support: 1-Piece with Buckle, 3" Height x 18" Width, Small	14924	E0960	\$146
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 20" Width, Medium	14927	E0960	\$146
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 22" Width, Large	14928	E0960	\$146
Dynaform Anterior Trunk Support: Standard, Dynamic, X-Small	14851	E0960	\$179
Dynaform Anterior Trunk Support: Standard, Dynamic, Small	14852	E0960	\$179
Dynaform Anterior Trunk Support: Standard, Dynamic, Medium	14853	E0960	\$179
Dynaform Anterior Trunk Support: Standard, Dynamic, Large	14854	E0960	\$179
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, X-Small	14856	E0960	\$179
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Small	14857	E0960	\$179
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Medium	14858	E0960	\$179
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Large	14859	E0960	\$179
AirLogic Anterior Trunk Support: Standard, Structure, Small	14282	E0960	\$203
AirLogic Anterior Trunk Support: Standard, Structure, Medium	14283	E0960	\$203
AirLogic Anterior Trunk Support: Standard, Structure, Large	14284	E0960	\$203
AirLogic Anterior Trunk Support: Standard, Structure, X-Large	14285	E0960	\$203
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Small	14257	E0960	\$203
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Medium	14258	E0960	\$203
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Large	14259	E0960	\$203
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, X-Large	14260	E0960	\$203

Accessories

Adjustable Shoulder Guides

\$101

Lumbar Pad

\$52

Lumbar Pad available for 12w to 22W Back width only. For odd widths a size up in width will be provided.