



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



Company: \_\_\_\_\_ Ship To: \_\_\_\_\_  
 Account Number: \_\_\_\_\_ Date: \_\_\_\_\_ Address: \_\_\_\_\_  
 Quote: PO Number: \_\_\_\_\_ City: \_\_\_\_\_ State: \_\_\_\_\_  
 ATP: \_\_\_\_\_ Zip: \_\_\_\_\_ Phone: \_\_\_\_\_  
 Marked For: \_\_\_\_\_ Email/Fax: \_\_\_\_\_  
 Associated Quote/Order#: \_\_\_\_\_ Effective: 1/5/2026

Flip - Rigid Base - Transit Included (E1233, E1161)	_____	\$3398
Flip - Rigid Base - Transit Not Included (E1233, E1161)	_____	\$3142
Flip - Folding Base - Transit Included (E1234, E1161)	_____	\$3563
Flip - Folding Base - Transit Not Included (E1234, E1161)	_____	\$3307

165 lb weight capacity Tilt range: 0° - 45° Maximum Frame Width of 18"

## Frame

*Growth Capabilities: Minimum 2" of seat width and minimum 3" of seat depth*

Frame Width									
Frame Width									
10"	11"	12"	13"	14"	15"	16"	17"	18"	

Frame Depth						
Select Frame Depth within desired Frame Depth Range. Depth cannot exceed largest orderable size.						
Frame Depth Range	Frame Depth					
Short (12" - 16")	12" N/C	13" N/C	14" N/C	15" N/C	16" N/C	
Medium (14" - 18")	14" N/C	15" N/C	16" N/C	17" N/C	18" N/C	
Long (16" - 20")	16" N/C	17" N/C	18" N/C	19" N/C	20" (E2203) \$672	

Transit Option	
Transit System (E1022)	\$256
<i>No charge when selected with Base Options with Transit Included</i>	
Omit Transit	N/C

Offset Hanger Option	
Offset Hanger	\$37
<i>Available only with Rigid Base</i>	

Tilt Mechanism	
Hand Tilt	N/C
Foot Tilt	N/C

Seat Pan	
Growing Solid Seat Pan (E2231)	\$241
Omit Seat Pan	N/C



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



## Backrest

Back Post / Height Range / Back Height		Adjustable Height Back Post Height			
Stroller Handle Back Post (Short 20")	N/C	18" Back Height	19" Back Height	20" Back Height	21" Back Height
Stroller Handle Back Post (Tall 24")	N/C	22" Back Height	23" Back Height	24" Back Height	25" Back Height
Fixed Height Back with Adjustable Handle (Short 22")	N/C	26" Back Height	27" Back Height		N/C
<i>Not available with Dynamic Back</i>					
Fixed Height Back with Adjustable Handle (Tall 26")	N/C				
<i>Not available with Dynamic Back</i>					
Reclining Back with Adjustable Handle (Short 22") (Non-Folding) (E1225)	\$962				
<i>Not available with Dynamic Back or Rigidizer Bar. Available only with Attendant Foot Lock. Minimum Seat Height of 15" with 12" Mag wheel and 16" with 16" Mag wheel.</i>					
Reclining Back with Adjustable Handle (Tall 26") (Non-Folding) (E1225)	\$962				
<i>Not available with Dynamic Back or Rigidizer Bar. Available only with Attendant Foot Lock. Minimum Seat Height of 15" with 12" Mag wheel and 16" with 16" Mag wheel.</i>					
Adjustable Height Back with Adjustable Handle (Short 18" - 25")	N/C				
<i>Not available with Dynamic Back</i>					
Adjustable Height Back with Adjustable Handle (Tall 20" - 27")	N/C				
<i>20"-22" Heights not available with Dynamic Back</i>					

Stroller Handle Back Post Option	
<i>Available only with Stroller Handle Back Post</i>	
Removable Stroller Handle Extension	\$101

Other Backrest Options	
Non Folding Back	N/C
<i>Standard with Dynamic Back</i>	
Dynamic Back (Non-Folding) (E2398)	\$1451
<i>Not available with Flip Up Armrests. Available only with Non Folding Back.</i>	
Adjustable Height Rigidizer Bar	\$86

## Footrest

Footrest Hanger Type					
Hanger	70°	80°	90°		Retail
Extension Tube Mount		1			N/C
Front Mount		1	1		N/C
Pro Elevating Leg Rests E0990 <sup>2</sup>					\$354
Pediatric Pro Elevating Leg Rests E0990 <sup>3</sup>					\$332
Omit Hangers					(\$20)

<sup>1</sup> Not available with 8" casters  
<sup>2</sup> Available only with Classic hanger release. Not available with Frame Widths less than 14" without Offset Hangers or Frame Widths less than 12" with Offset Hangers.  
<sup>3</sup> Available only with Classic Hanger Release

Pro ELR Cover	
<i>Available only on Pro ELR</i>	
Pro ELR Gel Pad Cover (PAIR)	\$468

Hanger Extension Tubes	
<i>Available only with Front Mount Hangers</i>	
3" Extension	\$59
3" Extension (Provided as Parts)	\$59



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



## Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements based on 70° Hanger Bend.

Footrest	Extension	Seat to Footrest Length											
		10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	15"	15.5"	18"	
Extension Tube Hangers (70°,80°)	Short	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	15"	15.5"		
	Medium	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"				
	Long	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"			
Front Mount Hangers <sup>1</sup>	Short Front Mount Hanger (70°,80°,90°)	Standard Selections									Requires 3" Extension		
		4" <sup>4</sup>	4.5" <sup>4</sup>	5"	5.5"	6"					8.5"	9"	9.5"
	Medium Front Mount Hanger (70°,80°,90°)	6.5"	7"	7.5"	8"						10"	10.5"	11"
		4" <sup>4</sup>	4.5" <sup>4</sup>	5"	5.5"	6"	6.5"	7"			10.5"	11"	11.5"
	Long Front Mount Hanger (70°,80°)	7.5"	8"	8.5"	9"	9.5"	10"				12"	12.5"	13"
		5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"		12.5"	13"	13.5"
		9"	9.5"	10"	10.5"	11"	11.5"	12"		14"	14.5"	15"	
Pro ELR Low Pivot Setting	Pro ELR Low Pivot Setting Short (10.5" - 16")												
	Pro ELR Low Pivot Setting Long (13" - 19.5")												
Pro ELR High Pivot Setting	Pro ELR High Pivot Setting Short (8.5" - 14.5")												
	Pro ELR High Pivot Setting Long (11.5" - 18")												
Peds Pro ELR	Low Pivot Setting	6" <sup>3</sup>	6.5" <sup>3</sup>	7" <sup>3</sup>	7.5"	8"	8.5"	9" <sup>2</sup>	9.5" <sup>2</sup>	10" <sup>2</sup>	10.5" <sup>2</sup>	11" <sup>2</sup>	11.5" <sup>2</sup>
	High Pivot Setting	6" <sup>3</sup>	6.5" <sup>3</sup>	7" <sup>3</sup>	7.5" <sup>2</sup>	8" <sup>2</sup>	8.5" <sup>2</sup>	9" <sup>2</sup>	9.5" <sup>2</sup>	10" <sup>2</sup>			

<sup>1</sup> See Minimum Front Seat Height chart below

<sup>2</sup> Comes standard with 3" Extension

<sup>3</sup> Not available on chairs 14" wide or more

<sup>4</sup> Available only with 90° Front Mount Hangers

Hanger Length	Short		Medium		Long	
	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13"	13"	13"	13.5"	14.5"	15.5"
Minimum Front Seat Height with 3" Extension	13.5"	14"	15.5"	16"	17"	18"

## Hanger Release

If no selection is made 4-Way - In-Line will be used

Classic	N/C
4-Way - In-Line	N/C
Not available with Pro ELR or Peds Pro ELR	
4-Way - Out	N/C
Not available with Pro ELR or Peds Pro ELR	

## Heel Loops

Heel Loops (E0951)	\$43
Not available with Frame Widths less than 14" without Offset Hangers or Frame Widths less than 12" with Offset Hangers	
Adjustable Heel Loops (E0951)	\$59
Omit Heel Loops	N/C

## Calf Straps

Velcro Adjustable (K0038)	\$48
Padded Velcro Adjustable (K0038)	\$75
Padded Airmesh Calf Strap (K0038)	\$134

## Footplates

Composite	N/C
Not available with Front Mount Hanger or Frame Widths less than 14"	
Composite Angle Adjustable (K0040)	\$203
Not available with Front Mount Hanger or 10"-13" Frame Widths	
Aluminum Angle Adjustable (K0040)	\$219
Aluminum Locking Multi-Angle Adjustable (K0108)	\$332
Available only with Front Mount Hanger	

## Residual Limb Support

	Limb Support	Price	Omit	Cover		Pad Size				
				Smooth Cover Left	Air Mesh Cover Left	5X8 Pad Left Side	7X10 Pad Left Side	9x10 Pad Left Side	9x14 Pad Left Side	11x10 Pad Left Side
Left	Left Limb Support (E1020)	\$449	Omit Left Hanger	Smooth Cover Left	Air Mesh Cover Left	5X8 Pad Left Side	7X10 Pad Left Side	9x10 Pad Left Side	9x14 Pad Left Side	11x10 Pad Left Side
Right	Right Limb Support (E1020)	\$449	Omit Right Hanger	Smooth Cover Right	Air Mesh Cover Right	5X8 Pad Right Side	7X10 Pad Right Side	9x10 Pad Right Side	9x14 Pad Right Side	11x10 Pad Right Side

## Seat Height, Casters, and Rear Wheels

Select Seat Height corresponding to desired Rear Wheel and Caster Size

Rear Wheel Size	Caster Size				
	5"	6"	6" x 2"	7"	8"
12" with Pneumatic or Pneumatic with Airless Insert Tire	Seat Height				
	13.5"				
	14"	14" <sup>4</sup>	14" <sup>4</sup>		
	14.5"	14.5"	14.5"		
	15"	15"	15"	15" <sup>4</sup>	
	15.5"	15.5"	15.5"	15.5"	
	16"	16"	16"	16"	
	16.5"	16.5"	16.5"	16.5"	16.5"
	17"	17"	17"	17" <sup>3</sup>	17" <sup>3</sup>
17.5"	17.5"	17.5"	17.5"	17.5"	
12" with Full Poly Tire	13.5"				
	14"				
	14.5"	14.5"	14.5"		
	15"	15"	15"		
	15.5"	15.5"	15.5"	15.5"	
	16"	16"	16"	16"	16" <sup>4</sup>
	16.5"	16.5"	16.5"	16.5"	16.5"
	17"	17"	17"	17" <sup>3</sup>	17" <sup>3</sup>
12" with Drum Brake <sup>13</sup>	13.5"				
	14"	14" <sup>4</sup>	14" <sup>4</sup>		
	14.5"	14.5"	14.5"		
	15"	15"	15"	15" <sup>4</sup>	
	15.5"	15.5"	15.5"	15.5"	
	16"	16"	16"	16"	16" <sup>4</sup>
	16.5"	16.5"	16.5"	16.5"	16.5"
	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"
16" with Pneumatic, Pneumatic with Airless Insert or Full Poly Tire	15"	15"	15"	15" <sup>4</sup>	
	15.5"	15.5"	15.5"	15.5"	
	16"	16"	16"	16"	
	16.5"	16.5"	16.5"	16.5"	16.5"
	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"
	18.5" <sup>4</sup>	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"
19.5" <sup>2</sup>	19.5" <sup>4</sup>	19.5" <sup>4</sup>	19.5" <sup>4</sup>		
16" with Low Poly Tire	14"	14" <sup>4</sup>	14" <sup>4</sup>		
	14.5"	14.5"	14.5"		
	15"	15"	15"		
	15.5"	15.5"	15.5"	15.5"	
	16"	16"	16"	16"	
	16.5"	16.5"	16.5"	16.5"	16.5"
	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"

1 Self Propel Axle Plate selections limit tilt range  
 2 Not available with 12" Frame Depth  
 3 Not available with 12-16" Frame Depth  
 4 Not available with 17-20" Frame Depth

Rear Wheel Size	Caster Size				
	5"	6"	6" x 2"	7"	8"
16" with Drum Brake	Seat Height				
	15"	15"	15"		
	15.5"	15.5"	15.5"	15.5"	
	16"	16"	16"	16"	16"
	16.5"	16.5"	16.5"	16.5"	16.5"
	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"
	18.5" <sup>4</sup>	18.5"	18.5" <sup>4</sup>	18.5" <sup>4</sup>	18.5"
	19"	19"	19"	19"	19"
20" with Pneumatic, Pneumatic with Airless Insert, Full Poly Tire Or Drum Brake	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"
	19.5"	19.5"	19.5"	19.5"	19.5"
	20" <sup>4</sup>	20" <sup>4</sup>	20" <sup>4</sup>	20"	20"
20" Self Propel <sup>1</sup> or Self Propel with Drum Brake	15.5" <sup>10</sup>	15.5" <sup>10</sup>	15.5" <sup>10</sup>	15.5" <sup>10</sup>	
	16" <sup>10</sup>	16" <sup>10</sup>	16" <sup>10</sup>	16" <sup>10</sup>	16" <sup>10</sup>
	16.5" <sup>10</sup>	16.5" <sup>10</sup>	16.5" <sup>10</sup>	16.5" <sup>10</sup>	16.5" <sup>10</sup>
	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"
	18.5" <sup>4</sup>	18.5"	18.5"	18.5" <sup>4</sup>	18.5"
22" Self Propel or Self <sup>1</sup> Propel with Drum Brake	16.5" <sup>10</sup>	16.5" <sup>10</sup>	16.5" <sup>10</sup>	16.5" <sup>10</sup>	16.5" <sup>10</sup>
	17" <sup>10</sup>	17" <sup>10</sup>	17" <sup>10</sup>	17" <sup>10</sup>	17" <sup>10</sup>
	17.5" <sup>10</sup>	17.5" <sup>10</sup>	17.5" <sup>10</sup>	17.5" <sup>10</sup>	17.5" <sup>10</sup>
	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"
24" Self Propel <sup>1</sup>	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"
	19.5" <sup>7,8</sup>	19.5"	19.5"	19.5" <sup>9</sup>	19.5" <sup>7,8</sup>
	20" <sup>5</sup>	20" <sup>6</sup>	20" <sup>6</sup>	20"	20"
24" Self Propel with <sup>1</sup> Drum Brake	18.5" <sup>4,11</sup>	18.5" <sup>11</sup>	18.5" <sup>11</sup>	18.5" <sup>11</sup>	18.5" <sup>11</sup>
	19"	19"	19"	19"	19"
	19.5" <sup>11</sup>	19.5" <sup>9,12</sup>	19.5" <sup>9,12</sup>	19.5" <sup>9,12</sup>	19.5" <sup>9,12</sup>
	20" <sup>4</sup>	20" <sup>6</sup>	20"	20"	20"

5 Not available in combination with spoke wheel and Full Poly tire  
 6 Not available in combination with spoke wheel, Full Poly tire and 17-20" Frame Depth  
 7 Not available in combination with Mag wheel and 17-20" Frame Depth  
 8 Not available in combination with spoke wheel, Pneumatic tires and 17-20" Frame Depth  
 9 17-20" Frame Depth requires spoke wheel and Full Poly tire  
 10 Not available with Drum Brake  
 11 Available only with Full Poly tire  
 12 12" Frame Depth requires spoke wheel and Full Poly tire  
 13 Not available with Pneumatic with Airless Insert



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



### Casters

		Caster Size				
		Price	5"	6"	7"	8"
Caster Type	.75" Lighted Rollerblade	\$81				
	1.25" Soft Roll	N/C				
	1" Poly Aluminum	\$123				
	1" Pneumatic	\$96				
	1.5" Soft Roll	\$101				
	1.5" Soft Roll Aluminum	\$198				
	2" Poly	\$101				

### Axles

*No selection needed with Drum Brake*

- Fixed \_\_\_\_\_ N/C
- Not available with Folding Base Frame*
- Quick Release \_\_\_\_\_ N/C

### Axle Plates

- Standard \_\_\_\_\_ N/C
- Self Propel Axle Plate \_\_\_\_\_ N/C

## Wheel Locks

### Wheel Lock

*Rear wheels must be removed to fold chair with Push to Lock, Pull to Lock and Dual wheel lock selections*

- Push to Lock \_\_\_\_\_ N/C
- Pull to Lock \_\_\_\_\_ N/C
- Not available with 12" or 16" Full Poly tires*
- Attendant Foot Lock \_\_\_\_\_ N/C
- Available only with 12" and 16" Mag wheels*
- Drum Brake (E2228) \_\_\_\_\_ \$619
- Available only with 12" Mag, 16" Mag, or MAXX Spoke wheels. Not available with 16" Low Poly tire. Axle is Quick Release. Requires Anti-Tips.*
- Dual - Attendant Foot Lock with Push to Lock \_\_\_\_\_ \$86
- Available only with 12" and 16" Mag wheels*
- Dual - Attendant Foot Lock with Pull to Lock \_\_\_\_\_ \$86
- Available only with 12" and 16" Mag wheels*
- Dual - Drum Brake with Push to Lock (E2228) \_\_\_\_\_ \$731
- Available only with 12" Mag, 16" Mag, or MAXX Spoke wheels. Not available with 16" Low Poly tire. Axle is Quick Release. Requires Anti-Tips.*
- Dual - Drum Brake with Pull to Lock (E2228) \_\_\_\_\_ \$731
- Available only with 12" Mag, 16" Mag, or MAXX Spoke wheels. Not available with 16" Low Poly tire. Axle is Quick Release. Requires Anti-Tips.*
- Omit \_\_\_\_\_ N/C

### Wheels

		Rear Wheel Size				
		12"	16"	20"	22"	24"
Rear Wheel Type	Mag	N/C				
	Colored Mag <sup>1</sup>	\$113				
	MAXX Spoke	\$214				

*1 Mag color selection is on page 7*

### Rear Tire Type

- Low Poly (E2218) \_\_\_\_\_ N/C
- Available only with 16" wheel size*
- Full Poly (E2218) \_\_\_\_\_ N/C
- Pneumatic (E2211,E2212) \_\_\_\_\_ \$96
- Pneumatic with Airless Insert (E2211, E2213) \_\_\_\_\_ \$101

### Handrims

*Handrims not available on 12" or 16" wheel sizes*

- Aluminum Anodized \_\_\_\_\_ N/C
- Plastic Coated \_\_\_\_\_ \$129
- Projection (8 Vertical Projections) \_\_\_\_\_ \$363
- Available only with full profile tires in combination with 22", 24" MAXX Spoke wheels, or 22" and 24" Mag wheels*
- Omit Handrim \_\_\_\_\_ N/C

### Handrim Option

*Not available on Plastic Coated or Projection Handrims*

- 3M Micro Grip \_\_\_\_\_ \$70

### Wheel Lock Option

- Extension Handles (E0961) \_\_\_\_\_ \$86
- Available only with Push to Lock or Pull to Lock Wheel Locks*

## Armrests

Armrest Type	Length	Price
Height Adjustable Pediatric E0973 <sup>1</sup>	Desk(6.5"-9")	\$203
	Desk Tall(8"-10.5")	
Height Adjustable T-Arm Tall (11.5"-15") E0973 <sup>2</sup>	Desk	\$337
	Full	
Height Adjustable T-Arm (8.5"-13.5") E0973 <sup>2</sup>	Desk	\$337
	Full	
Height Adjustable T-Arm Low (6.5"-10.5") E0973 <sup>2,3</sup>	Desk	\$337
	Full	

<sup>1</sup> Best fit with Frame Depths less than 14"

<sup>2</sup> Not available with 12" Frame Depth

<sup>3</sup> If requested configuration is unavailable, Standard Height T-Arm will be used

<sup>4</sup> Not available with Dynamic Back. Available only with 12" and 16" Rear Wheels.

<sup>5</sup> Not available with Dynamic Back, Reclining Back or Arm Pads

Armrest Type	Length	Price
Angle Adjustable Locking Extendable Flip Up K0020 <sup>4</sup>	Desk	\$309
	Full	
Tubular Flip Up <sup>5</sup>		N/C
Omit Armrests		N/C

## Arm Pads

Standard Arm Pad	N/C
Waterfall Arm Pad	\$28
<i>Not available with Angle Adjustable Locking Extendable Flip Up Armrest</i>	
Foam Grip	N/C
<i>Available only with Angle Adjustable Locking Extendable Flip Up Armrest</i>	

## Accessories

### Removable Side Guards

*Not available with Height Adjustable Armrests*

Adult Composite Side Guards (14" - 16") (K0108)	\$155
Pediatric Composite Side Guard (10" - 14") (K0108)	\$155

### Anti-Tips

Rear Anti-Tips (E0971)	\$172
Omit	N/C

### Pelvic Positioning Belts

1" Auto Buckle Pelvic Positioning Belt (E0978)	\$90
1" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$183
1 1/2" Auto Buckle Pelvic Positioning Belt (E0978)	\$91
1 1/2" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$183
2" Airline Buckle Pelvic Positioning Belt (E0978)	\$120
Small Flex-Mount Positioning Belt (16-24") (E0978)	\$208
Medium Flex-Mount Positioning Belt (19-29") (E0978)	\$208
Large Flex-Mount Positioning Belt (25-37") (E0978)	\$208

### Belt Mounting Kit

Belt Mounting Band Clamp	\$37
<i>Clamp includes hardware to mount Belt</i>	
Seat Belt Mounting Hardware	\$32

### Spoke Guards

Clear Spoke Guards (K0065)	\$125
<i>Not available with Mag wheel</i>	
Black Spoke Guards (K0065)	\$128
<i>Available only with 22" and 24" Spoke wheels</i>	
Blue Monster Spoke Guards (K0065) <sup>1</sup>	\$271
Robo Fun Spoke Guards (K0065) <sup>1</sup>	\$271
Tutu Cute Spoke Guards (K0065) <sup>1</sup>	\$271
Tuxedo Time Spoke Guards (K0065) <sup>1</sup>	\$271

<sup>1</sup> Available only with 20" and 22" wheel sizes. Not available with Mag wheel

### Other Accessories

Tool Kit	\$63
Removable Underseat Pouch	\$64
Pediatric Ki Mobility Backpack	\$74
Ki Mobility Backpack	\$83
Black Neoprene Impact Guards	\$99
Tilt Stop Posterior (Pair) (K0108)	\$33
Tilt Stop Anterior (Pair) (K0108)	\$33
<i>Available only with Rigid Base</i>	
IV Pole	\$303
Canopy with Hardware (K0108)	\$308

### Upper Extremity Support Tray

<sup>3/8"</sup> Polycarbonate with ABS ribs. Hardware ordered separately. BC = Body Cut Out. Armrest assembly is required for tray mounting.

Tray	Fits Chair	Retail
Tray, X-Small 21" x 18" (BC 10" x 7") (E0950)	8" - 12"	\$235
Tray, Small 22" x 20" (BC 13" x 8") (E0950)	10" - 14"	\$235
Tray, Medium 24" x 22" (BC 15" x 9") (E0950)	12" - 16"	\$235
Tray, Large 26" x 24" (BC 17" x 10") (E0950)	14" - 18"	\$235
Tray, X-Large 28" x 26" (BC 18" x 10") (E0950)	16" - 20"	\$255
Tray Hardware Assy, Knob Release - Pair		\$222

## Color

Select desired frame and accent colors. Accent color is used on side panel; if no accent color is selected, frame color will be used.

Options	Price	Colors																				
		Gloss													Matte							
		Black	Ruby Red	Radiant Orange	Yellow	Electric Green	Rich Teal	Glacier Blue	Blue Steel	Deep Sapphire Blue	Purple Plum	Opalescent Violet	Dark Pink	Light Pink	Champagne	Shadow Black	Iridescent Navy	Bright Red	Tactical Green	Dry Earth Tan	Platinum	White
Frame	N/C																					
Accent	N/C																					
Mag Color																						

Touch-up paint is available for purchase through the parts manual



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



## Axiom for Kids P - Positioning Cushion

*Weight Capacity 165 lbs / 75 kg. Cushion includes 2 outer covers of the selected color. Cover Color Code Black=BLK, Blue=BLU, Purple=PUR.*

Width	8"	9"	10"	11"	12"	13"	14"	15"	16"	
Retail	\$405									
HCPCS	E2605	E2605	E2605	E2605	E2605	E2605	E2605	E2605	E2605	
Depth	8"	PXP0808BLK	PXP0908BLK	PXP1008BLK						
		PXP0808BLU	PXP0908BLU	PXP1008BLU						
		PXP0808PUR	PXP0908PUR	PXP1008PUR						
	9"	PXP0809BLK	PXP0909BLK	PXP1009BLK	PXP1109BLK					
		PXP0809BLU	PXP0909BLU	PXP1009BLU	PXP1109BLU					
		PXP0809PUR	PXP0909PUR	PXP1009PUR	PXP1109PUR					
	10"	PXP0810BLK	PXP0910BLK	PXP1010BLK	PXP1110BLK	PXP1210BLK	PXP1310BLK			
		PXP0810BLU	PXP0910BLU	PXP1010BLU	PXP1110BLU	PXP1210BLU	PXP1310BLU			
		PXP0810PUR	PXP0910PUR	PXP1010PUR	PXP1110PUR	PXP1210PUR	PXP1310PUR			
	11"			PXP1011BLK	PXP1111BLK	PXP1211BLK	PXP1311BLK	PXP1411BLK		
				PXP1011BLU	PXP1111BLU	PXP1211BLU	PXP1311BLU	PXP1411BLU		
				PXP1011PUR	PXP1111PUR	PXP1211PUR	PXP1311PUR	PXP1411PUR		
	12"			PXP1012BLK	PXP1112BLK	PXP1212BLK	PXP1312BLK	PXP1412BLK		
				PXP1012BLU	PXP1112BLU	PXP1212BLU	PXP1312BLU	PXP1412BLU		
				PXP1012PUR	PXP1112PUR	PXP1212PUR	PXP1312PUR	PXP1412PUR		
	13"			PXP1013BLK	PXP1113BLK	PXP1213BLK	PXP1313BLK	PXP1413BLK	PXP1513BLK	
				PXP1013BLU	PXP1113BLU	PXP1213BLU	PXP1313BLU	PXP1413BLU	PXP1513BLU	
				PXP1013PUR	PXP1113PUR	PXP1213PUR	PXP1313PUR	PXP1413PUR	PXP1513PUR	
	14"			PXP1014BLK	PXP1114BLK	PXP1214BLK	PXP1314BLK	PXP1414BLK	PXP1514BLK	PXP1614BLK
				PXP1014BLU	PXP1114BLU	PXP1214BLU	PXP1314BLU	PXP1414BLU	PXP1514BLU	PXP1614BLU
				PXP1014PUR	PXP1114PUR	PXP1214PUR	PXP1314PUR	PXP1414PUR	PXP1514PUR	PXP1614PUR
	15"					PXP1215BLK	PXP1315BLK	PXP1415BLK	PXP1515BLK	PXP1615BLK
						PXP1215BLU	PXP1315BLU	PXP1415BLU	PXP1515BLU	PXP1615BLU
						PXP1215PUR	PXP1315PUR	PXP1415PUR	PXP1515PUR	PXP1615PUR
16"					PXP1216BLK	PXP1316BLK	PXP1416BLK	PXP1516BLK	PXP1616BLK	
					PXP1216BLU	PXP1316BLU	PXP1416BLU	PXP1516BLU	PXP1616BLU	
					PXP1216PUR	PXP1316PUR	PXP1416PUR	PXP1516PUR	PXP1616PUR	
17"						PXP1317BLK	PXP1417BLK	PXP1517BLK	PXP1617BLK	
						PXP1317BLU	PXP1417BLU	PXP1517BLU	PXP1617BLU	
						PXP1317PUR	PXP1417PUR	PXP1517PUR	PXP1617PUR	
18"						PXP1318BLK	PXP1418BLK	PXP1518BLK	PXP1618BLK	
						PXP1318BLU	PXP1418BLU	PXP1518BLU	PXP1618BLU	
						PXP1318PUR	PXP1418PUR	PXP1518PUR	PXP1618PUR	

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

**DCN1075.1 Pricing and specifications subject to change without notice.**



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



## Axiom for Kids SP - Skin Protection and Positioning Cushion

*Weight Capacity 165 lbs. / 75 kg. Cushion includes 2 outer covers of the selected color. Cover Color Code Black=BLK, Blue=BLU, Purple=PUR.*

Width	8"	9"	10"	11"	12"	13"	14"	15"	16"	
Retail	\$491									
HCPCS	E2607	E2607	E2607	E2607	E2607	E2607	E2607	E2607	E2607	
Depth	8"	PXSP0808BLK	PXSP0908BLK	PXSP1008BLK						
		PXSP0808BLU	PXSP0908BLU	PXSP1008BLU						
		PXSP0808PUR	PXSP0908PUR	PXSP1008PUR						
	9"	PXSP0809BLK	PXSP0909BLK	PXSP1009BLK	PXSP1109BLK					
		PXSP0809BLU	PXSP0909BLU	PXSP1009BLU	PXSP1109BLU					
		PXSP0809PUR	PXSP0909PUR	PXSP1009PUR	PXSP1109PUR					
	10"	PXSP0810BLK	PXSP0910BLK	PXSP1010BLK	PXSP1110BLK	PXSP1210BLK	PXSP1310BLK			
		PXSP0810BLU	PXSP0910BLU	PXSP1010BLU	PXSP1110BLU	PXSP1210BLU	PXSP1310BLU			
		PXSP0810PUR	PXSP0910PUR	PXSP1010PUR	PXSP1110PUR	PXSP1210PUR	PXSP1310PUR			
	11"			PXSP1011BLK	PXSP1111BLK	PXSP1211BLK	PXSP1311BLK	PXSP1411BLK		
				PXSP1011BLU	PXSP1111BLU	PXSP1211BLU	PXSP1311BLU	PXSP1411BLU		
				PXSP1011PUR	PXSP1111PUR	PXSP1211PUR	PXSP1311PUR	PXSP1411PUR		
	12"			PXSP1012BLK	PXSP1112BLK	PXSP1212BLK	PXSP1312BLK	PXSP1412BLK		
				PXSP1012BLU	PXSP1112BLU	PXSP1212BLU	PXSP1312BLU	PXSP1412BLU		
				PXSP1012PUR	PXSP1112PUR	PXSP1212PUR	PXSP1312PUR	PXSP1412PUR		
	13"			PXSP1013BLK	PXSP1113BLK	PXSP1213BLK	PXSP1313BLK	PXSP1413BLK	PXSP1513BLK	
				PXSP1013BLU	PXSP1113BLU	PXSP1213BLU	PXSP1313BLU	PXSP1413BLU	PXSP1513BLU	
				PXSP1013PUR	PXSP1113PUR	PXSP1213PUR	PXSP1313PUR	PXSP1413PUR	PXSP1513PUR	
	14"			PXSP1014BLK	PXSP1114BLK	PXSP1214BLK	PXSP1314BLK	PXSP1414BLK	PXSP1514BLK	PXSP1614BLK
				PXSP1014BLU	PXSP1114BLU	PXSP1214BLU	PXSP1314BLU	PXSP1414BLU	PXSP1514BLU	PXSP1614BLU
				PXSP1014PUR	PXSP1114PUR	PXSP1214PUR	PXSP1314PUR	PXSP1414PUR	PXSP1514PUR	PXSP1614PUR
	15"					PXSP1215BLK	PXSP1315BLK	PXSP1415BLK	PXSP1515BLK	PXSP1615BLK
						PXSP1215BLU	PXSP1315BLU	PXSP1415BLU	PXSP1515BLU	PXSP1615BLU
						PXSP1215PUR	PXSP1315PUR	PXSP1415PUR	PXSP1515PUR	PXSP1615PUR
16"					PXSP1216BLK	PXSP1316BLK	PXSP1416BLK	PXSP1516BLK	PXSP1616BLK	
					PXSP1216BLU	PXSP1316BLU	PXSP1416BLU	PXSP1516BLU	PXSP1616BLU	
					PXSP1216PUR	PXSP1316PUR	PXSP1416PUR	PXSP1516PUR	PXSP1616PUR	
17"						PXSP1317BLK	PXSP1417BLK	PXSP1517BLK	PXSP1617BLK	
						PXSP1317BLU	PXSP1417BLU	PXSP1517BLU	PXSP1617BLU	
						PXSP1317PUR	PXSP1417PUR	PXSP1517PUR	PXSP1617PUR	
18"						PXSP1318BLK	PXSP1418BLK	PXSP1518BLK	PXSP1618BLK	
						PXSP1318BLU	PXSP1418BLU	PXSP1518BLU	PXSP1618BLU	
						PXSP1318PUR	PXSP1418PUR	PXSP1518PUR	PXSP1618PUR	

### Accessories


#### Extra Cover

Extra Cover Black	\$120
Extra Cover Blue	\$120
Extra Cover Purple	\$120

#### Solid Seat Insert

Solid Seat Insert (E0992)	\$73
---------------------------	------


Axiom for Kids Back Width (all models)	Fits back Frame Widths
10" (25.5 cm)	9" - 11" (23 cm - 28 cm)
12" (30.5 cm)	11" - 13" (28 cm - 33 cm)
14" (35.5 cm)	13" - 15" (33 cm - 38 cm)
16" (40.5 cm)	15" - 17" (38 cm - 43 cm)



## Axiom for Kids PL - Posterior Lateral Back

Not available on Ki Mobility Wheelchairs with 2" Inward Offset or 2" Outward Offset back posts. Not available on Ki Mobility Wheelchairs with Tubular Flip-Up Armrests. May have mounting conflicts with low back cane heights and Angle Adjustable Locking Extendable Flip-Up Armrests with Rigidizer Bar.

Back Width	10" (25.5 cm)	12" (30.5 cm)	14" (35.5 cm)	16" (40.5 cm)
<b>Fits Back Frame Widths:</b>	9"-11" (23 cm - 28 cm)	11"-13" (28 cm - 33 cm)	13"-15" (33 cm - 38 cm)	15"-17" (38 cm - 43 cm)
<b>Actual Back Width (F):</b>	10" (25.5 cm)	12" (30.5 cm)	14" (35.5 cm)	16" (40.5 cm)
<b>Fits Trunk Widths (F):</b>	10" (25.5 cm)	12" (30.5 cm)	14" (35.5 cm)	16" (40.5 cm)
<b>Contour Depth (H):</b>	3" (7.5 cm)	3" (7.5 cm)	3" (7.5 cm)	3" (7.5 cm)



### Back

Includes Lumbar Pad. Weight capacity 165 lbs. (75 kg). Includes Axiom DEX for Kids Hardware.

		Wheelchair Frame Width			
		10	12	14	16
Retail		\$755			
HCPCS		E2615	E2615	E2615	E2615
Axiom Back Height	10	PXPL1010			
	12	PXPL1012	PXPL1212	PXPL1412	
	14	PXPL1014	PXPL1214	PXPL1414	PXPL1614
	16		PXPL1216	PXPL1416	PXPL1616
	18				PXPL1618

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.


Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1076.5 Pricing and specifications subject to change without notice.

## Axiom for Kids AL - Posterior Back for Adjustable Laterals

*Not available on Ki Mobility Wheelchairs with 2" Inward Offset or 2" Outward Offset back posts. Not available on Ki Mobility Wheelchairs with Tubular Flip-Up Armrests. May have mounting conflicts with low back cane heights and Angle Adjustable Locking Extendable Flip-Up Armrests with Rigidizer Bar.*

<b>Back Width</b>	10" (25.5 cm)	12" (30.5 cm)	14" (35.5 cm)	16" (40.5 cm)
<b>Fits Back Frame Widths:</b>	9"-11" (23 cm - 28 cm)	11"-13" (28 cm - 33 cm)	13"-15" (33 cm - 38 cm)	15"-17" (38 cm - 43 cm)
<b>Actual Back Width (F):</b>	8" (20cm)	10" (25.5 cm)	12" (30.5 cm)	14" (35.5 cm)
<b>Component Lateral Width Range (J)*:</b> <i>(*-1" with included spacers)</i>	8" - 11" (20 cm - 28 cm)	10" - 13" (25.5 cm - 33 cm)	12" - 15" (30.5 cm - 38 cm)	14" - 17" (35.5 cm - 43 cm)
<b>Contour Depth (H):</b>	1.25" (3 cm)	1.25" (3 cm)	1.25" (3 cm)	1.25" (3 cm)



### Back

*Weight capacity 165 lbs. (75 kg). Includes Axiom DEX for Kids Hardware.*

		Wheelchair Frame Width			
		10"	12"	14"	16"
Axiom Back Height	Retail	\$887			
	HCPCS	E2291	E2291	E2291	E2291
	10"	PXAL1010			
	12"	PXAL1012	PXAL1212	PXAL1412	
	14"	PXAL1014	PXAL1214	PXAL1414	PXAL1614
	16"		PXAL1216	PXAL1416	PXAL1616
	18"			PXAL1618	

### Axiom for Kids - Adjustable Lateral Trunk Supports (for PXAL Back)

*Not available for Posterior Back for Contoured Adjustable Laterals (PXALCL)*

Description	Item Code	HCPCS	Retail
3 x 4 Lateral Trunk Pads, Flat with Fixed Mounting Hardware (Pair) (E0956)	PLF3x4PP	E0956	\$311
4 x 5 Lateral Trunk Pads, Flat with Fixed Mounting Hardware (Pair) (E0956)	PLF4x5PP	E0956	\$311
3 x 4 Lateral Trunk Pads, Flat with Swing Away Mounting Hardware (Pair) (E0956,E1034)	PLTAS3x4PP	E0956, E1034	\$806
4 x 5 Lateral Trunk Pads, Flat with Swing Away Mounting Hardware (Pair) (E0956,E1034)	PLTAS4x5PP	E0956, E1034	\$806
3 x 4 Lateral Trunk Pads, Curved with Swing Away Mounting Hardware (Pair) (E0956,E1034)	PLTASC3x4PP	E0956, E1034	\$806
4 x 5 Lateral Trunk Pads, Curved with Swing Away Mounting Hardware (Pair) (E0956,E1034)	PLTASC4x5PP	E0956, E1034	\$806

### Axiom for Kids - Adjustable Lateral Pelvic Supports (for PXAL Back)


*Not available for Posterior Back for Contoured Adjustable Laterals (PXALCL)*

Description	Item Code	HCPCS	Retail
3 x 4 Lateral Pelvic Pads with Fixed Mounting Hardware (Pair) (E0956)	PLPF3x4PP	E0956	\$311
4 x 5 Lateral Pelvic Pads with Fixed Mounting Hardware (Pair) (E0956)	PLPF4x5PP	E0956	\$311
3 x 4 Lateral Pelvic Pads with Swing Away Mounting Hardware (Pair) (E0956,E1034)	PLASP3x4PP	E0956, E1034	\$806
4 x 5 Lateral Pelvic Pads with Swing Away Mounting Hardware (Pair) (E0956,E1034)	PLASP4x5PP	E0956, E1034	\$806

## Axiom for Kids ALCL - Posterior Back for Contoured Adjustable Laterals

*Not available on Ki Mobility Wheelchairs with 2" Inward Offset or 2" Outward Offset back posts. Not available on Ki Mobility Wheelchairs with Tubular Flip-Up Armrests. May have mounting conflicts with low back cane heights and Angle Adjustable Locking Extendable Flip-Up Armrests with Rigidizer Bar.*

<b>Back Width</b>	10" (25.5 cm)	12" (30.5 cm)	14" (35.5 cm)	16" (40.5 cm)
<b>Fits Back Frame Widths:</b>	9"-11" (23 cm - 28 cm)	11"-13" (28 cm - 33 cm)	13"-15" (33 cm - 38 cm)	15"-17" (38 cm - 43 cm)
<b>Actual Back Width (F):</b>	8" - 11" (20 cm - 28 cm)	10" - 13" (25.5 cm - 33 cm)	12" - 15" (30.5 cm - 38 cm)	14" - 17" (35.5 cm - 43 cm)
<b>Contoured Lateral Width Range (J):</b>	8" - 11" (20 cm - 28 cm)	10" - 13" (25.5 cm - 33 cm)	12" - 15" (30.5 cm - 38 cm)	14" - 17" (35.5 cm - 43 cm)
<b>Contour Depth (K):</b>	5" - 5.5" (12.5 cm - 14 cm)	5" - 5.5" (12.5 cm - 14 cm)	5" - 5.5" (12.5 cm - 14 cm)	5" - 5.5" (12.5 cm - 14 cm)



### Back

*Requires Contoured Adjustable Lateral Trunk Support selection below. Weight capacity 165 lbs. (75 kg). Includes Axiom DEX for Kids Hardware.*

		Wheelchair Frame Width			
		10"	12"	14"	16"
Retail		\$989			
HCPCS		E2293	E2293	E2293	E2293
Axiom Back Height	10"	PXAL1010CL			
	12"	PXAL1012CL	PXAL1212CL	PXAL1412CL	
	14"	PXAL1014CL	PXAL1214CL	PXAL1414CL	PXAL1614CL
	16"		PXAL1216CL	PXAL1416CL	PXAL1616CL
	18"				PXAL1618CL

### Axiom for Kids - Contoured Adjustable Lateral Trunk Supports (for PXALCL Back)

Description	Item Code	HCPCS	Retail
3" Contoured Lateral Wings with Fixed Mounting Hardware (Pair) (E0956)	PCL3	E0956	\$332
4" Contoured Lateral Wings with Fixed Mounting Hardware (Pair) (E0956)	PCL4	E0956	\$332

*4" Contoured Lateral Wings are not available on 10" or 12" Back Heights*

## Accessories

### Axiom for Kids Back Accessories & Modifications

Lumbar Pad (K0108)	\$54
Privacy Flap (K0108)	\$54
<i>Privacy Flap not available for 10" Wide backs</i>	
Axiom DEX for Kids Hardware Extension (K0108)	\$97
<i>Hardware Extension adds an additional 2" of anterior depth adjustment</i>	
Axiom for Kids Adjustable Shoulder Guides (K0108)	\$99

## Headrest Pads

Selection of a Headrest Pad requires a selection from Headrest Mounting Hardware section below



Contour Positioning



Tri-Pad Adjustable



Comfort Soft



TriSoft



Comfort Soft Adjustable

Headrest Pads			
Description	Item Code	HCPCS	Retail
Contour Positioning Headrest Pad, Small, 8W x 6.25H (E0955)	116919	E0955	\$190
Contour Positioning Headrest Pad, Medium, 9W x 6.5H (E0955)	116921	E0955	\$190
Contour Positioning Headrest Pad, Large, 10W x 7.5H (E0955)	116923	E0955	\$190
Tri-Pad Adjustable Headrest, Small (E0955)	116925	E0955	\$460
Tri-Pad Adjustable Headrest, Medium (E0955)	116927	E0955	\$460
Tri-Pad Adjustable Headrest, Large (E0955)	116929	E0955	\$460
Comfort Soft Headrest, 8" (E0955)	116931	E0955	\$204
Comfort Soft Headrest, 10" (E0955)	116933	E0955	\$204
Comfort Soft Headrest, 14" (E0955)	116935	E0955	\$204
Comfort Soft Headrest, 19" (E0955)	116937	E0955	\$235
TriSoft Headrest, 10" (E0955)	116939	E0955	\$204
TriSoft Headrest, 14" (E0955)	116941	E0955	\$204
Comfort Soft Adjustable Headrest, 10" (E0955)	116943	E0955	\$383
Comfort Soft Adjustable Headrest, 14" (E0955)	116945	E0955	\$383
Comfort Soft Adjustable Headrest, 19" (E0955)	116947	E0955	\$477

See Axiom Hardware Sizing Guide for full item sizing

## Headrest Hardware

Selection of Headrest Hardware requires a selection from Headrest Pad section above



Std. Straight Removable Adjustable



Multi-Axis Removable Adjustable



Omnilink Multi-Axis Adjustable



Posalinc Removable Adjustable



Posalinc Flip-Back Adjustable

Headrest Hardware			
Description	Item Code	HCPCS	Retail
Standard Straight Removable Adjustable Hardware (E1033)	120038	E1033	\$198
Multi-Axis Removable Adjustable Hardware (E1033)	120039	E1033	\$295
Omnilink Multi-Axis Adjustable Hardware (E1033)	123603	E1033	\$417
POSALinc Removable Adjustable Hardware (E1033)	120040	E1033	\$359
POSALinc Removable Adjustable Hardware, Tall (E1033)	123528	E1033	\$396
POSALinc Flip-Back Adjustable Hardware (E1033)	123530	E1033	\$430
POSALinc Flip-Back Adjustable Hardware, Tall (E1033)	123531	E1033	\$474

Universal Headrest Mounting Plate  
 Not required with Axiom Headrests

\$54

### Anterior Trunk Supports



2-Piece  
Hook & Loop



1-Piece  
with Buckle



Dynaform Standard,  
Dynamic



Dynaform Slim Cut,  
Dynamic



Airlogic Standard,  
Structure



Airlogic Slim Cut,  
Dynamic

### Anterior Trunk Supports

Anterior Trunk Supports			
Description	Item Code	HCPCS	Retail
Dynaform Chest Support: 2-Piece with Hook&Loop, 3" High Pad, Small	14941	E0960	\$146
Dynaform Chest Support: 2-Piece with Hook&Loop, 4" High Pad, Medium	14942	E0960	\$146
Dynaform Chest Support: 2-Piece with Hook&Loop, 5" High Pad, Large	14943	E0960	\$146
Dynaform Chest Support: 1-Piece with Buckle, 3" Height x 18" Width, Small	14924	E0960	\$146
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 20" Width, Medium	14927	E0960	\$146
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 22" Width, Large	14928	E0960	\$146
Dynaform Anterior Trunk Support: Standard, Dynamic, X-Small	14851	E0960	\$179
Dynaform Anterior Trunk Support: Standard, Dynamic, Small	14852	E0960	\$179
Dynaform Anterior Trunk Support: Standard, Dynamic, Medium	14853	E0960	\$179
Dynaform Anterior Trunk Support: Standard, Dynamic, Large	14854	E0960	\$179
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, X-Small	14856	E0960	\$179
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Small	14857	E0960	\$179
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Medium	14858	E0960	\$179
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Large	14859	E0960	\$179
AirLogic Anterior Trunk Support: Standard, Structure, Small	14282	E0960	\$203
AirLogic Anterior Trunk Support: Standard, Structure, Medium	14283	E0960	\$203
AirLogic Anterior Trunk Support: Standard, Structure, Large	14284	E0960	\$203
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Small	14257	E0960	\$203
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Medium	14258	E0960	\$203

### Gloss Finish Back Shell

If no selection is made, standard Shadow Black will be used. Touch-up paint is available for purchase through the parts manual.

Black	N/C
Ruby Red	N/C
Radiant Orange	N/C
Yellow	N/C
Electric Green	N/C
Rich Teal	N/C
Glacier Blue	N/C
Blue Steel	N/C
Deep Sapphire Blue	N/C
Purple Plum	N/C
Opalescent Violet	N/C
Dark Pink	N/C
Light Pink	N/C
Champagne	N/C

### Matte Finish Back Shell

If no selection is made, standard Shadow Black will be used. Touch-up paint is available for purchase through the parts manual.

Iridescent Navy	N/C
Bright Red	N/C
Tactical Green	N/C
Dry Earth Tan	N/C
Platinum	N/C
White	N/C