

AXIOM

POSITIONING

AXIOM POSITIONING PRODUCTS

SIZING GUIDE


TABLE OF CONTENTS

Axiom Soft Goods	1
Anterior Trunk Supports.....	1
Airlogic 1-Piece Chest Support	1
Dynaform 1-Piece Chest Support	1
Dynaform 2-piece Chest Support.....	1
2-piece Chest Support.....	1
Airlogic Standard	2
Airlogic Slim Cut.....	2
Dynaform Standard.....	3
Dynaform Slim Cut.....	3
Anterior Pelvic Supports.....	4
2-Point, Padded, Center Pull - Push Button Buckle	4
2-Point, Padded, Center Pull - Side Release Buckle.....	4
2-Point, Padded, Rear Pull - Push Button Buckle.....	4
2-Point, Non-Padded, Center Pull - Push Button Buckle.....	5
2-Point, Non-Padded, Center Pull - Side Release Buckle.....	5
4-Point, Padded, Center Pull - Push Button Buckle	6
4-Point, Padded, Center Pull - Side Release Buckle.....	6
4-Point, Padded, Rear Pull - Push Button Buckle.....	6
Accessories - Slip-On Pads	6
Lower Extremity	7
Foot/Ankle Positioning - 3-Point Ankle Positioners, Side Release	7
Foot/Ankle Positioning - 3-Point Foot/Ankle Positioners.....	7
Foot/Ankle Positioning - 2-Point Ankle Straps, Side Release	7
Calf Supports - Airlogic Calf Panel	8
Calf Supports - Dynaform Calf Panel.....	8
Shoe Holders.....	8
Axiom Hardware	9
Seat.....	9
Drop Seat Base	9
Headrest Pads.....	10
Contour Positioning Headrest Pad	10
Tri-Pad Adjustable Headrest Pad	10
Comfort Soft Headrest Pad.....	11
Comfort Soft Adjustable Headrest Pad.....	11
Tri-Soft Headrest Pad.....	11

AXIOM SOFT GOODS

All Dimensions are +/- 1/4".


Anterior Trunk Supports




Airlogic 1-Piece Chest Support

Part Number	Size	Height	Width	Recommended Chest Width	Webbing
14222	XX-Small	3"	14"	6" - 8"	1"
14223	X-Small	3"	16"	8" - 10"	1"
14224	Small	3"	18"	10" - 12"	1"
14227	Medium	4"	20"	12" - 14"	1"
14228	Large	4"	22"	14" - 16"	1"
14229	X-Large	4"	24"	16" - 18"	1"
14230	XX-Large	5"	26"	18" - 20"	1"

Dynaform 1-Piece Chest Support



Part Number	Size	Height	Width	Recommended Chest Width	Webbing
14922	XX-Small	3"	14"	6" - 8"	1"
14923	X-Small	3"	16"	8" - 10"	1"
14924	Small	3"	18"	10" - 12"	1"
14927	Medium	4"	20"	12" - 14"	1"
14928	Large	4"	22"	14" - 16"	1"
14929	X-Large	4"	24"	16" - 18"	1"






Dynaform 2-piece Chest Support

Part Number	Size	Height	Recommended Chest Width	Webbing
14941	Small	3"	6" - 12"	1"
14942	Medium	4"	10" - 16"	1 1/2"
14943	Large	5"	12" - 18"	2"

2-piece Chest Support

Part Number	Size	Width		Webbing
14221	76"		Non-Padded	2"
14220	76"	12"	Padded	2"






Airlogic Standard

Part Number - Dynamic	Part Number - Structured	Size	A	B	C	D	E	F	Webbing
14250		XX-Small	6"	10 1/2"	4 1/2"	5 1/4"	1 1/4"	2"	1"
14251	14281	X-Small	8"	14"	4 1/2"	8 1/2"	4"	2"	1"
14252	14282	Small	9 1/4"	17"	6"	10"	5 1/2"	2"	1"
14253	14283	Medium	10 1/4"	20"	7 1/2"	11 1/4"	6 1/4"	2 1/4"	1"
14254	14284	Large	11 1/4"	22"	8 1/4"	12 1/2"	7"	2 1/2"	1"
14255		X-Large	12 1/4"	24"	9"	13 1/2"	7 1/2"	2 3/4"	1"

Airlogic Slim Cut


Part Number - Dynamic	Part Number - Structured	Size	A	B	C	D	E	F	Webbing
14257	14287	Small	9 1/4"	17"	6"	10"	1 1/4"	2"	1"
14258	14288	Medium	10 1/4"	20"	7 1/2"	11 1/4"	1 1/4"	2 1/4"	1"



Dynaform Standard									
Part Number - Dynamic	Part Number - Structured	Size	A	B	C	D	E	F	Webbing
14851	14881	X-Small	8"	14"	4"	8 1/2"	4"	2"	3/4"
14852	14882	Small	9"	17"	5"	10"	5"	2"	3/4"
14853	14883	Medium	10"	20"	6"	11 1/2"	6"	2"	1"
14854	14884	Large	11"	22"	7"	12"	7"	2 1/4"	1"


Dynaform Slim Cut									
Part Number - Dynamic	Part Number - Structured	Size	A	B	C	D	E	F	Webbing
	14886	X-Small	8"	14"	4"	8 1/2"	1"	2"	3/4"
14857	14887	Small	9"	17"	5"	10"	1"	2"	3/4"
14858	14888	Medium	10"	20"	6"	11 1/2"	1 1/4"	2"	3/4"
	14889	Large	11"	22"	7"	12"	1 1/2"	2 1/2"	1"

Anterior Pelvic Supports



2-Point, Padded, Center Pull		Push Button Buckle		
Part Number	Buckle Type	Size	Pad Width	Webbing
14633	Non-Spring Metal	X-Small	6"	1"
14634	Non-Spring Metal	Small	8"	1"
14643	Spring	Small	8"	1 1/2"
14663	Non-Spring Metal	Small	8"	1 1/2"
14644	Spring	Medium	10"	1 1/2"
14653	Spring	Medium	10"	2"
14654	Spring	Large	12"	2"

2-Point, Padded, Center Pull		Side Release Buckle		
Part Number	Buckle Type	Size	Pad Width	Webbing
14603	Side Release	X-Small	6"	1"
14604	Side Release	Small	8"	1"
14613	Side Release	Small	8"	1 1/2"
14614	Side Release	Medium	10"	1 1/2"
14623	Side Release	Medium	10"	2"
14624	Side Release	Large	12"	2"



2-Point, Padded, Rear Pull		Push Button Buckle		
Part Number	Buckle Type	Size	Pad Width	Webbing
14435	Non-Spring Metal	X-Small	6"	1"
14645	Spring	Small	8"	1 1/2"
14446	Spring	Medium	10"	1 1/2"
14447	Spring	Large	12"	2"



2-Point, Non-Padded, Center Pull		Push Button Buckle		
Part Number	Buckle Type	Size	Webbing Length	Webbing
14630	Non-Spring Metal	Small	60"	1"
14647	Spring	Medium	72"	1 1/2"
14658	Spring	Large	72"	2"
14659	Spring	X-Large	85"	2"



2-Point, Non-Padded, Center Pull		Side Release Buckle		
Part Number	Buckle Type	Size	Webbing Length	Webbing
14606	Side Release	Small	60"	1"
14620	Side Release	Medium	72"	1 1/2"
14621	Side Release	Large	72"	2"
14622	Side Release	X-Large	85"	2"

4-Point, Padded, Center Pull		Push Button Buckle		
Part Number	Buckle Type	Size	Pad Width	Webbing
14749	Non-Spring Metal	Small	5 1/2"	1"
14748	Non-Spring Metal	Medium	7 1/2"	1"
14750	Spring	Medium	7 1/2"	1 1/2"
14751	Spring	Large	9 1/4"	2"
14752	Spring	X-Large	12"	2"



4-Point, Padded, Center Pull		Side Release Buckle		
Part Number	Buckle Type	Size	Pad Width	Webbing
14754	Side Release	Small	5 1/2"	1"
14753	Side Release	Medium	7 1/2"	1"
14755	Side Release	Medium	7 1/2"	1 1/2"
14756	Side Release	Large	9 1/4"	2"
14757	Side Release	X-Large	12"	2"

4-Point, Padded, Rear Pull		Push Button Buckle		
Part Number	Buckle Type	Size	Pad Width	Webbing
14448	Non-Spring Metal	Small	5 1/2"	1"
14449	Non-Spring Metal	Medium	7 1/2"	1"
14450	Spring	Medium	7 1/2"	1 1/2"
14451	Spring	Large	9 1/4"	2"



Accessories		Slip-On Pads		
Part Number	Size	Height	Width	Webbing
14592	Small		8"	1"
14594	Medium		10"	1 1/2"
14596	Large		12"	2"

Lower Extremity

Foot/Ankle Positioning		3-Point Ankle Positioners, Side Release			
Part Number	Size	Ankle Circumference	Foot Base Length	Foot Base Width	Width Between Ankle Slots
13150-2	X-Small	4 1/2" - 7"	7 1/2"	4 1/2"	2 3/4"
13151-2	Small	6 1/2" - 9"	7 1/2"	4 1/2"	2 3/4"
13152-2	Medium	8 1/2" - 11"	9"	5"	3"
13153-2	Large	10 1/2" - 13"	10"	5 1/2"	3 1/2"
13154-2	X-Large	12 1/2" - 15"	10"	5 1/2"	3 1/2"



Foot/Ankle Positioning		3-Point Foot/Ankle Positioners			
Part Number	Size	Ankle Circumference	Foot Base Length	Foot Base Width	Width Between Ankle Slots
13161-2	Small	6 1/2" - 9"	7 1/2"	4 1/2"	2 3/4"
13162-2	Medium	8 1/2" - 11"	9"	5"	3"
13163-2	Large	10 1/2" - 13"	10"	5 1/2"	3 1/2"



Foot/Ankle Positioning		2-Point Ankle Straps, Side Release
Part Number	Size	Ankle Circumference
13130	X-Small	4 1/2" - 7"
13131	Small	6 1/2" - 9"
13132	Medium	8 1/2" - 11"
13133	Large	10 1/2" - 13"
13134	X-Large	12 1/2" - 15"



Calf Supports Airlogic Calf Panel					
Part Number	Size	Height	Width	Frame Size	Webbing
13179	Medium	6 1/2"	23"	16" - 18"	1"



Calf Supports Dynaform Calf Panel					
Part Number	Size	Height	Width	Frame Size	Webbing
13175	Small	5"	19"	12" - 15"	1"
13176	Medium	6 1/2"	23"	16" - 18"	1"
13177	Large	8"	26"	19" - 22"	1"

Shoe Holders					
Part Number	Size	Width at Toe	Width at Ankle	Length	Number of Straps
13011-2	Small	2 3/4"	2 3/4"	6"	1
13012-2	Medium	4"	3"	8"	2
13013-2	Large	4 1/2"	3 1/2"	10"	2




AXIOM HARDWARE

All Dimensions are +/- 1/4".

Seat

Drop Seat Base




Part Number	Wheelchair Frame	Base Depth (A)	Base Width (B)	Crossbar Cut Out (C)
FDSV12/12-K	12"	12"	9"	1 1/2"
FDSV12/14-K	12"	14"	9"	1 1/2"
FDSV14/14-K	14"	14"	11"	2"
FDSV14/16-K	14"	16"	11"	2"
FDSV16/14-K	16"	14"	13"	2 1/2"
FDSV16/16-K	16"	16"	13"	2 1/2"
FDSV16/18-K	16"	18"	13"	2 1/2"
FDSV18/16-K	18"	16"	15"	3"
FDSV18/18-K	18"	18"	15"	3"
FDSV18/20-K	18"	20"	15"	3"
FDSV20/16-K	20"	16"	17"	3 1/2"
FDSV20/18-K	20"	18"	17"	3 1/2"



Headrest Pads

Contour Positioning Headrest Pad

Part Number	Size	HE	HA	Foam Thickness
116919	Small	8"	6 1/4"	1 3/4"
116921	Medium	9"	6 1/2"	1 3/4"
116923	Large	10"	7 1/2"	1 3/4"




Tri-Pad Adjustable Headrest Pad

Part Number	Size	HE	HA	HFL/R	HAL/R	HBL/R	HDL/R	Foam Thickness
116925	Small	5"	3 1/4"	4 1/4"	2 3/4"	2"	1"	1/2"
116927	Regular	5"	4 1/2"	5 1/2"	3 1/2"	2 3/4"	1 1/2"	1/2"
116929	Large	6"	6 1/4"	6 1/4"	4"	3"	1 3/4"	1/2"


Comfort Soft Headrest Pad

Part Number	Size	HE	HA	Foam Thickness
116931	8"	8"	5"	1 3/4"
116933	10"	10"	5 1/2"	1 3/4"
116935	14"	14"	6 1/2"	1 3/4"
116937	19"	19"	6 1/2"	1 3/4"




Comfort Soft Adjustable Headrest Pad

Part Number	Size	HY	HX	HZ	HA	Foam Thickness
116943	10"	5"	3 1/2"	3 1/2"	5 1/2"	1 3/4"
116945	14"	5"	5 1/2"	5 1/2"	6 1/2"	1 3/4"
116947	19"	6 1/2"	7"	7"	6 1/2"	1 3/4"



Tri-Soft Headrest Pad

Part Number	Size	HY	HX	HZ	HA	Foam Thickness
116939	10"	5"	4"	4"	5 1/2"	1 3/4"
116941	14"	6 1/2"	5"	5"	6 1/2"	1 3/4"





Ki Mobility
5201 Woodward Drive
Stevens Point, Wisconsin 54481
715-254-0991
www.kimobility.com