



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com



Company:	Ship To:		
Account Number:	Date:	Address:	
Quote:	PO Number:	City:	State:
ATP:	Zip:	Phone:	
Marked For:	Email/Fax:		
Associated Quote/Order#:	Effective: 1/5/2026		

Spark- Transit Included (E1238)	\$2635
Spark - Transit Not Included (E1238)	\$2294
165 lb weight capacity Not available with upholstery	

Frame

Growth capabilities: Minimum 2" of seat width and minimum 3" of seat depth

Seat Width

Width growth capabilities stated in image to the right

10"	N/C
11"	N/C
12"	N/C
13"	N/C
14"	N/C
15"	N/C
16"	N/C

Seat Depth

Depth growth capabilities stated in image to the right


12"	N/C
13"	N/C
14"	N/C
15"	N/C
16"	N/C
17"	N/C
18"	N/C

Transit Option

Transit System (E1022)	\$342
<i>No charge when selected with Spark - Transit Included</i>	
Omit Transit	N/C


Seat Pan

Removable Aluminum Seat Pan (E2231)	\$312
<i>Pan depth 1/2" shorter than Frame Depth</i>	
Omit	N/C



Frame Width

Spark is set up to allow 2" increase in frame width from initial frame width in 1" increments. Growth is achieved by expanding cross braces. 16" maximum frame width. 10" Frame Width expands only 1".



Frame Depth

Spark is set up to allow 4" increase in frame depth from initial frame depth in 1" increments. 19" maximum frame depth.

Wheelchair Drive Ready

Chair will ship from Ki Mobility with adaption kit installed for selected Alber system. Order system separately from Alber USA. Select Rear Frame Height (full profile selections only) in drive wheel size on Rear Frame Height chart (page 4). Requires Rear Anti-Tips. Not available with Spinerdy wheel.

e-motion Ready _____ N/C



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com



Backrest

Backrest

Back Posts are folding and angle adjustable in 4° increments

Back Post	Price	Height Range	Back Height							
			20"	24"	13"	14"	15"	16"	17"	
Stroller Handle Back Post (Depth Adjustable)	N/C	Short	20"							
		Tall	24"							
Height Adjustable Straight with Push Handle (Depth Adjustable)	\$81	Short ¹	13"	14"	15"	16"	17"			
		Medium	16"	17"	18"	19"	20"			
		Tall	20" ²	21"	22"	23"	24"			

¹ Not available with Dynamic Back

² Not available with Angle Adjustable Locking Extendable Flip Up Armrest

Backrest Options

Removable Stroller Handle Extension	\$84
<i>Not available with Height Adjustable Straight with Push Handle Back Posts</i>	
Dynamic Back (E2398)	\$1451
<i>Available only with T-Arms. Not available with 12" Frame Depths with Rear Frame Heights lower than 14.5".</i>	

Frame Type / Footrest

Frame Type

Swing Away	N/C
Super Low - Swing Away	\$54

Footrest Hanger Type

Hanger	60°	70°	80°	90°		Retail
Extension Tube Mount			1			N/C
Front Mount			1	1		N/C
Pro Elevating Leg Rests E0990 ^{2,3}						\$354
Pediatric Pro ELR ³						\$332
Omit Hangers						(\$20)

¹ Not available with 7" or 8" Casters

² Not available with Seat Widths less than 14"

³ Available only with Classic Hanger Release

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
Hanger Bend	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13"	13.5"	14"	15"	16"	17"
Minimum Front Seat Height with 3" Extension	14.5"	16"	16.5"	18"	18.5"	20"

Footrest / Extension / Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements are based on 70° Hanger.

Footrest	Extension	Seat to Footrest Length											
		11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"
Extension Tube Hangers (60°, 70°, 80°)	Short	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"			
	Medium	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"			
	Long	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"			
Front Mount ¹ Hangers	Short Front Mount Hanger (70°, 80°, 90°)	Standard Selections									Requires 3" Extension		
		4" ⁴	4.5" ⁴	5"	5.5"	6"	6.5"				10"	10.5"	11"
	Medium Front Mount Hanger (70°, 80°, 90°)	7"	7.5"	8"	8.5"	9"	9.5"				11.5"	12"	12.5"
		4" ⁴	4.5" ⁴	5"	5.5"	6"	6.5"	7"	7.5"		12"	12.5"	13"
	Long Front Mount Hanger (70°, 80°)	8"	8.5"	9"	9.5"	10"	10.5"	11"	11.5"		13.5"	14"	14.5"
		5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"	9"	14"	14.5"	15"
PRO ELR Low Pivot Setting	PRO ELR Low Pivot Setting Short (12" - 18")	9.5"	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	15.5"	16"	16.5"
PRO ELR High Pivot Setting	PRO ELR High Pivot Setting Short (10.5" - 16")												
	PRO ELR High Pivot Setting Long (13" - 19.5")												
Peds PRO ELR	Low Pivot Setting	7.5"	8" ³	8.5" ³	9"	9.5"	10"	10.5" ²	11" ²	11.5" ²	12" ²	12.5" ²	13" ²
	High Pivot Setting	7.5" ³	8" ³	8.5" ³	9" ²	9.5" ²	10" ²	10.5" ²	11" ²	11.5" ²			

¹ See Minimum Front Seat Height chart above

² Comes standard with 3" Extension

³ Not available on chairs 14" wide or more

⁴ Available with 90° only



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com



3" Extension

Available only with Front Mount Hangers

3" Extension	\$59
3" Extension (Provided as Parts)	\$59

Footplates

Composite	N/C
<i>Not available with Front Mount Hanger or Frame Widths less than 14"</i>	
Composite Angle Adjustable (K0040)	\$208
<i>Not available with Front Mount Hanger or Frame Widths less than 16"</i>	
Aluminum Angle Adjustable (K0040)	\$214
Aluminum Locking Multi-Angle Adjustable (K0108)	\$321
<i>Available only with Front Mount Hanger</i>	
One Piece Flip Up Angle Adjustable (K0040)	\$208
<i>Not available with Front Mount Hanger, ELR, Peds Pro ELR, Pro ELR or Frame Widths less than 14"</i>	

Hanger Release

If no selection is made Classic Hanger Release will be used

Classic	N/C
4-Way - In-Line	N/C
4-Way - Out	N/C

Heel Loops

Not available with One Piece Flip Up Angle Adjustable Footplate

Heel Loops (E0951)	\$43
<i>Not available with Frame Widths less than 14"</i>	
Adjustable Heel Loops (E0951)	\$54
Omit Heel Loops	N/C

Calf Straps

Velcro Adjustable (K0038)	\$43
Padded Velcro Adjustable (K0038)	\$75
Padded Airmesh Calf Strap (K0038)	\$139

Residual Limb Support

	Limb Support	Price	Omit	Cover		Pad Size				
				Smooth Cover	Air Mesh Cover	5X8 Pad	7X10 Pad	9x10 Pad	9x14 Pad	11x10 Pad
Left	Left Limb Support (E1020)	\$438	Omit Left Hanger	Smooth Cover Left	Air Mesh Cover Left	5X8 Pad Left Side	7X10 Pad Left Side	9x10 Pad Left Side	9x14 Pad Left Side	11x10 Pad Left Side
Right	Right Limb Support (E1020)	\$438	Omit Right Hanger	Smooth Cover Right	Air Mesh Cover Right	5X8 Pad Right Side	7X10 Pad Right Side	9x10 Pad Right Side	9x14 Pad Right Side	11x10 Pad Right Side

Front Seat Height and Casters

Front Seat Height

Select Front Seat Height in row corresponding to desired caster size

Caster Size	Front Seat Height													
	Uses Super Low Front Frame													
3"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"			
4"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"
5"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	
6"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"			
7"	17"	17.5"		18"	18.5"	19"	19.5"	20"	20.5"					
8"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"							

Castors

Caster Type	Caster Size	Price	Caster Size				
			3"	4"	5"	6"	7"
.75" Rollerblade	N/C						
.75" Lighted Rollerblade	\$81						
1.25" Soft Roll	N/C						
1" Poly Aluminum	\$123						
1" Pneumatic	\$96						
1.5" Soft Roll	\$59						
1.5" Soft Roll Aluminum	\$208						
2" Poly	\$101						
2" Pneumatic w/Foam Insert	\$256						

Maximum Seat Dump Reference

(Front Seat Height minus Rear Seat Height)

Seat Depth	12"	13"	14"	15"	16"	17"	18"	19"
Maximum Dump Standard Caster Housing	1.5"	1.5"	1.5"	2"	2"	2"	2.5"	2.5"
Maximum Dump Performance Caster Housing	2.5"	2.5"	3"	3"	3"	3.5"	3.5"	4"



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com



Fork Stem

Best Fit	N/C
<i>Ki Mobility will pick Fork and Stem combination to optimize performance based on Caster Size and Front Seat Height</i>	
Standard Fork Stem	N/C
3/4" Longer Fork Stem	N/C
1 1/2" Longer Fork Stem	N/C
Frog Legs Suspension Fork	\$508
<i>Available only with 4" - 6" Caster Sizes. Not available with Caster Pin Locks, 2" Poly Caster Type or 2" Pneumatic with Foam Insert Caster Type.</i>	

Fork Options

Frog Leg Options

User Weight Under 80 lbs	N/C
User Weight 80-120 lbs	N/C
User Weight 120-160 lbs	N/C

Caster Pin Locks

Caster Pin Locks (K0073)	\$219
<i>Not available with Performance Caster Housing, 1 1/2" Longer Fork Stem or Frog Legs Fork</i>	

Caster Housing

Standard Housing	N/C
Performance Housing	\$54

Wheel Locks

Wheel Lock

Push to Lock	N/C
Pull to Lock	N/C
Omit	N/C

Wheel Lock Options

Extension Handles (E0961)	\$81
Grade Aids (E0974)	\$161
<i>Available only with Push to Lock Wheel Locks and Pneumatic and Pneumatic with Airless Insert tires</i>	

Rear Seat Height and Wheels

Rear Seat Height

Select Rear Seat Height in row corresponding to Wheel Size and Tire Profile. 12"-14" Seat Depths not available with 24" wheel size.

Wheel Size	Tire Profile	Rear Seat Height										
		18"	Low	11" ²	11.5" ¹	12" ²	12.5" ¹	13" ²	13.5"	14" ²	14.5"	15" ²
	Full	11.5" ²	12" ¹	12.5" ²	13" ¹	13.5" ²	14"	14.5" ²	15"	15.5" ²	16" ²	16.5" ²
20"	Low	11.5" ²	12" ¹	12.5" ²	13" ¹	13.5" ²	14"	14.5" ²	15"	15.5" ²	16" ²	16.5" ²
	Full	12" ²	12.5" ¹	13" ²	13.5" ¹	14" ²	14.5"	15" ²	15.5"	16" ²	16.5" ²	17" ²
22"	Low	12.5" ²	13" ¹	13.5" ²	14" ¹	14.5" ²	15"	15.5" ²	16"	16.5" ²	17" ²	17.5" ²
	Full	13" ²	13.5" ¹	14" ²	14.5" ¹	15" ²	15.5"	16" ²	16.5"	17" ²	17.5" ²	18" ²
24"	Low	13.5" ²	14" ¹	14.5" ²	15" ¹	15.5" ²	16"	16.5" ²	17"	17.5" ²	18" ²	18.5" ²
	Full	14" ²	14.5" ¹	15" ²	15.5" ¹	16" ²	16.5"	17" ²	17.5"	18" ²	18.5" ²	19" ²

¹ Amputee Axle Plate not available with lowest two Rear Seat Height positions

² Requires Universal Axle Plate

One Arm Drive

One Arm Drive - Quick Release (E0958)	\$1168
<i>Available only with 20", 22" and 24" spoke wheels and Quick Release Axles. Spoke wheel style is exclusive to OAD. Not available with 4 degrees of camber.</i>	

Axle Plates

No selection needed with One Arm Drive

Universal	N/C
Amputee (E0959)	\$195
<i>Not available with Wheelchair Drive Ready Option</i>	

Drive Handrim Side

Left Side Drive	N/C
Right Side Drive	N/C

Drive Handrim Type

Aluminum Drive Handrim	N/C
Poly Coated Drive Handrim	\$214



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com



Handrims

Standard Handrims

Aluminum Anodized	N/C
Plastic Coated	\$113
Projection (8 Vertical Projections)	\$374
<i>Available only with full profile tires in combination with 22" and 24" MAXX Spoke wheels, or 22" and 24" Mag wheels</i>	
Omit Handrim	N/C

Ergonomic Handrim

Not available with 20" wheels

Lite Regular Grip	\$497
Standard Regular Grip	\$475
Lite Extra Grip	\$550
Standard Extra Grip	\$528
Lite Omit Grip	\$422
Standard Omit Grip	\$369

Wheels

		Rear Wheel Size				
		Price	18"	20"	22"	24"
Rear Wheel Type	MAXX Spoke	\$134				
	Mag	N/C				
	Colored Mag ¹	\$113				
	Spinergy Spox	\$1014				

¹ Color selection on page 7

Spinergy Options

Select spoke color for Spinergy wheels

Black Spokes Red Spokes Blue Spokes White Spokes _____ N/C

Axles

Axles

No selection needed with One Arm Drive

Quick Release	N/C
Quad Quick Release	\$198
<i>Not available with Mag or Spinergy wheels</i>	

Handrim Options

3M Micro Grip	\$75
<i>Not available on Plastic Coated or Projection Handrims</i>	

Handrim Mount

Close Mount	\$33
<i>If no selection is made, Standard Mount will be used. Reduces overall width by .75".</i>	

Rear Tire Type

Full Profile

Street Pneumatic 1 3/8" (E2211, E2212)	N/C
Full Poly 1 3/8" (E2218)	\$59
Pneumatic with Airless Insert 1 3/8" (E2211, E2213)	\$101

Low Profile

Low Poly 1" (E2218)	\$59
Sentinel Solid Tire 1"	\$101
<i>Not available on 18" wheel size or 20" and 22" Mag wheels</i>	
High Pressure 1"	\$101
<i>Not available on Mag wheel or 18" wheel size</i>	

Axle Position

Rear Wheel Camber

0° 2° 4°	N/C
--------------	-----

Center of Gravity

No selection needed with Amputee Axle Plate

Rear Center Forward	N/C
-------------------------	-----

Armrests

Armrest Type	Length	Price
Height Adjustable T-Arm Tall (11.5"-15") E0973 ¹	Desk	\$337
	Full	
Height Adjustable T-Arm (8.5"-13.5") E0973 ¹	Desk	\$337
	Full	
Height Adjustable T-Arm Low (6.5"-10.5") E0973 ²	Desk	\$337
	Full	

Armrest Type	Length	Price
Angle Adjustable Locking Extendable Flip Up K0020	Desk	\$284
	Full	
Pediatric Height Adjustable T-Arm Desk (6.5"-10")		\$241
Tubular Flip Up ³		N/C
Omit Armrests		N/C

¹ Not available with Frame Depths less than 14"

² Not available with Frame Depths less than 14". If requested configuration is unavailable, Standard Height T-Arm will be used.

³ Not available with Arm Pads

Arm Pads

Standard Arm Pad	N/C
Waterfall Arm Pad	\$28
<i>Not available with Angle Adjustable Locking Extendable Flip Up Armrest</i>	
Foam Grip	N/C
<i>Available only with Angle Adjustable Locking Extendable Flip Up Armrest</i>	

Accessories

Removable Side Guards

Not available with Height Adjustable T-Arm Armrests

Adult Composite Side Guards (14" - 16") (K0108)	\$166
Pediatric Composite Side Guard (10" - 14") (K0108)	\$166

Pelvic Positioning Belts

1" Auto Buckle Pelvic Positioning Belt (E0978)	\$90
1" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$183
1 1/2" Auto Buckle Pelvic Positioning Belt (E0978)	\$91
1 1/2" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$183
2" Airline Buckle Pelvic Positioning Belt (E0978)	\$120
Small Flex-Mount Positioning Belt (16-24") (E0978)	\$208
Medium Flex-Mount Positioning Belt (19-29") (E0978)	\$208
Large Flex-Mount Positioning Belt (25-37") (E0978)	\$208
Belt Mounting Band Clamp	\$37
<i>Clamp includes hardware to mount Belt</i>	

Upper Extremity Support Tray

3/8" Polycarbonate with ABS ribs. Hardware ordered separately. BC = Body Cut Out. Armrest assembly is required for tray mounting.

Tray	Fits Chair	Retail
Tray, X-Small 21" x 18" (BC 10" x 7") (E0950)	8" - 12"	\$235
Tray, Small 22" x 20" (BC 13" x 8") (E0950)	10" - 14"	\$235
Tray, Medium 24" x 22" (BC 15" x 9") (E0950)	12" - 16"	\$235
Tray, Large 26" x 24" (BC 17" x 10") (E0950)	14" - 18"	\$235
Tray, X-Large 28" x 26" (BC 18" x 10") (E0950)	16" - 20"	\$255
Tray Hardware Assy, Knob Release - Pair		\$222

Anti-tips


Rear Anti-Tips (E0971)	\$161
Omit	N/C

Spoke Guards

Clear Spoke Guards (K0065)	\$125
<i>Not available with Mag wheel</i>	
Black Spoke Guards (K0065)	\$128
<i>Not available with 18" or 20" wheel sizes or Mag wheel</i>	
Blue Monster Spoke Guards (K0065)	\$271
<i>Available only with 20" and 22" wheel sizes. Not available with Mag wheel.</i>	
Robo Fun Spoke Guards (K0065)	\$271
<i>Available only with 20" and 22" wheel sizes. Not available with Mag wheel.</i>	
Tutu Cute Spoke Guards (K0065)	\$271
<i>Available only with 20" and 22" wheel sizes. Not available with Mag wheel.</i>	
Tuxedo Time Spoke Guards (K0065)	\$271
<i>Available only with 20" and 22" wheel sizes. Not available with Mag wheel.</i>	

Other Accessories

Tool Kit	\$63
Removable Underseat Pouch	\$64
<i>Not available with 10" - 13" Seat Widths</i>	
Pediatric Ki Mobility Backpack	\$74
Ki Mobility Backpack	\$83
Black Neoprene Impact Guards	\$99
Cane and Crutch Holder	\$214
Canopy with Hardware (K0108)	\$308
<i>See image below</i>	



Color

Touch-up paint available for purchase in the parts manual

Options	Price	Color																			
		Gloss													Matte						
		Black	Ruby Red	Radiant Orange	Yellow	Electric Green	Rich Teal	Glacier Blue	Blue Steel	Deep Sapphire Blue	Purple Plum	Opalescent Violet	Dark Pink	Light Pink	Champagne	Shadow Black	Iridescent Navy	Bright Red	Tactical Green	Dry Earth Tan	Platinum
Frame	N/C																				
Mag Color																					
Color Accent Package ¹	N/C																				

¹ Includes axle plates and forks. When One Arm Drive or Amputee Axle Plate is selected, only forks will be painted. When Frog Legs is selected, only axle plate will be painted.

Hydrographics Options

Not available with Frame Color. Hydrographics Base Color required with Hydrographics. Zebra not available with black.

Options	Price	Hydrographics					Hydrographics Base Color						
		AOR1 (Desert Digital Camo)	NEXT Camo	Pink Camo	Carbon Fiber	Muddy Girl	Zebra	Black	Blue	Granny Smith	Pink	Purple	Red
Frame	\$859												
Base Color	N/C												



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com



Spark Front Frame Height Matrix																							
3" Caster				4" Caster				5" Caster				6" Caster				7" Caster				8" Caster			
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork
14.5	13.5	Std	4"	15	13.5	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	16.5	Std	6"	19	17.5	Std	6"
15	14			15.5	14			16.5	15.5	3/4		17.5	16.5	3/4		18.5	17			19.5	18.5	3/4	
15.5	14	15.5 ¹		14 ¹	17.5	16		1 1/2	18.5	17		1 1/2	18.5	17.5		20.5	n/a	1 1/2					
16	14.5	3/4		16	14.5	3/4		16	14.5	Std	17	15.5	Std	19	18	3/4	19	17.5		Std			
16	15	16.5 ¹		15 ¹	16.5	15		17.5	16		19.5	18		20	18.5		19.5	18					
16.5	15.5	1 1/2		17	15.5	1 1/2		16.5	15.5	3/4	18	16.5	3/4	20	18.5	20	18.5						
15.5	14.5	Std		5"	15.5	14		Std	5"	17	16	5"	17.5	16.5	6"	18.5	17	Std		7"	20.5	19	3/4
16	15				16	14.5				17.5	16		1 1/2	18		17	3/4				19	17.5	
16.5	15	3/4			16	15		3/4		18	16.5		1 1/2	18.5		17.5	19.5	18			20	19	
17	15.5	1 1/2			16.5	15.5		3/4		16.5	15	Std	18.5	17	1 1/2	19	17.5	1 1/2			20	18.5	3/4
17	16		17		15.5	17	15.5			19	17.5		19	18		21	n/a						
17.5	16.5	Std	17.5		16	Std	17.5	16		3/4	18.5	17	Std	19.5	18	3/4	20.5	19	1 1/2		20.5	19	
16.5	15.5		17		15.5		18	16.5			18.5	17		20	19		21	n/a					
17	16	3/4	17		16	3/4	18	16.5		1 1/2	18.5	17	1 1/2	19	17.5	1 1/2	20	18.5	1 1/2		21	n/a	
18	16.5		17.5		16.5		18	16.5			18.5	17		20.5	19								
18	17	1 1/2	18		16.5	1 1/2	18.5	17		1 1/2	18.5	17.5	3/4	19	18	3/4	21	n/a	1 1/2		21	n/a	
18.5	17.5		18	16.5	18.5		17	19	17.5		19	18											
				18	16.5	Std	7"	18	16.5	Std	7"	19.5	18.5	1 1/2	6" x 2" caster not available with 5" fork								
				18.5	17			18.5	17			20	19										
				18.5	17.5	3/4		19	17.5	19.5		18											
				19	18	1 1/2		18.5	17.5	1 1/2		20	18.5										
				19.5	18.5			19	18			20.5	19										
				20	19	19.5		18.5	21	n/a													

¹ Not available with 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com



Axiom for Kids P - Positioning Cushion

Weight Capacity 165 lbs / 75 kg. Cushion includes 2 outer covers of the selected color. Cover Color Code Black=BLK, Blue=BLU, Purple=PUR.

Width	8"	9"	10"	11"	12"	13"	14"	15"	16"	
Retail	\$405									
HCPCS	E2605	E2605	E2605	E2605	E2605	E2605	E2605	E2605	E2605	
Depth	8"	PXP0808BLK	PXP0908BLK	PXP1008BLK						
		PXP0808BLU	PXP0908BLU	PXP1008BLU						
		PXP0808PUR	PXP0908PUR	PXP1008PUR						
	9"	PXP0809BLK	PXP0909BLK	PXP1009BLK	PXP1109BLK					
		PXP0809BLU	PXP0909BLU	PXP1009BLU	PXP1109BLU					
		PXP0809PUR	PXP0909PUR	PXP1009PUR	PXP1109PUR					
	10"	PXP0810BLK	PXP0910BLK	PXP1010BLK	PXP1110BLK	PXP1210BLK	PXP1310BLK			
		PXP0810BLU	PXP0910BLU	PXP1010BLU	PXP1110BLU	PXP1210BLU	PXP1310BLU			
		PXP0810PUR	PXP0910PUR	PXP1010PUR	PXP1110PUR	PXP1210PUR	PXP1310PUR			
	11"			PXP1011BLK	PXP1111BLK	PXP1211BLK	PXP1311BLK	PXP1411BLK		
				PXP1011BLU	PXP1111BLU	PXP1211BLU	PXP1311BLU	PXP1411BLU		
				PXP1011PUR	PXP1111PUR	PXP1211PUR	PXP1311PUR	PXP1411PUR		
	12"			PXP1012BLK	PXP1112BLK	PXP1212BLK	PXP1312BLK	PXP1412BLK		
				PXP1012BLU	PXP1112BLU	PXP1212BLU	PXP1312BLU	PXP1412BLU		
				PXP1012PUR	PXP1112PUR	PXP1212PUR	PXP1312PUR	PXP1412PUR		
	13"			PXP1013BLK	PXP1113BLK	PXP1213BLK	PXP1313BLK	PXP1413BLK	PXP1513BLK	
				PXP1013BLU	PXP1113BLU	PXP1213BLU	PXP1313BLU	PXP1413BLU	PXP1513BLU	
				PXP1013PUR	PXP1113PUR	PXP1213PUR	PXP1313PUR	PXP1413PUR	PXP1513PUR	
	14"			PXP1014BLK	PXP1114BLK	PXP1214BLK	PXP1314BLK	PXP1414BLK	PXP1514BLK	PXP1614BLK
				PXP1014BLU	PXP1114BLU	PXP1214BLU	PXP1314BLU	PXP1414BLU	PXP1514BLU	PXP1614BLU
				PXP1014PUR	PXP1114PUR	PXP1214PUR	PXP1314PUR	PXP1414PUR	PXP1514PUR	PXP1614PUR
	15"					PXP1215BLK	PXP1315BLK	PXP1415BLK	PXP1515BLK	PXP1615BLK
						PXP1215BLU	PXP1315BLU	PXP1415BLU	PXP1515BLU	PXP1615BLU
						PXP1215PUR	PXP1315PUR	PXP1415PUR	PXP1515PUR	PXP1615PUR
16"					PXP1216BLK	PXP1316BLK	PXP1416BLK	PXP1516BLK	PXP1616BLK	
					PXP1216BLU	PXP1316BLU	PXP1416BLU	PXP1516BLU	PXP1616BLU	
					PXP1216PUR	PXP1316PUR	PXP1416PUR	PXP1516PUR	PXP1616PUR	
17"						PXP1317BLK	PXP1417BLK	PXP1517BLK	PXP1617BLK	
						PXP1317BLU	PXP1417BLU	PXP1517BLU	PXP1617BLU	
						PXP1317PUR	PXP1417PUR	PXP1517PUR	PXP1617PUR	
18"						PXP1318BLK	PXP1418BLK	PXP1518BLK	PXP1618BLK	
						PXP1318BLU	PXP1418BLU	PXP1518BLU	PXP1618BLU	
						PXP1318PUR	PXP1418PUR	PXP1518PUR	PXP1618PUR	

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1075.2 Pricing and specifications subject to change without notice.



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com



Axiom for Kids SP - Skin Protection and Positioning Cushion

Weight Capacity 165 lbs. / 75 kg. Cushion includes 2 outer covers of the selected color. Cover Color Code Black=BLK, Blue=BLU, Purple=PUR.

Width	8"	9"	10"	11"	12"	13"	14"	15"	16"	
Retail	\$491									
HCPCS	E2607	E2607	E2607	E2607	E2607	E2607	E2607	E2607	E2607	
Depth	8"	PXSP0808BLK	PXSP0908BLK	PXSP1008BLK						
		PXSP0808BLU	PXSP0908BLU	PXSP1008BLU						
		PXSP0808PUR	PXSP0908PUR	PXSP1008PUR						
	9"	PXSP0809BLK	PXSP0909BLK	PXSP1009BLK	PXSP1109BLK					
		PXSP0809BLU	PXSP0909BLU	PXSP1009BLU	PXSP1109BLU					
		PXSP0809PUR	PXSP0909PUR	PXSP1009PUR	PXSP1109PUR					
	10"	PXSP0810BLK	PXSP0910BLK	PXSP1010BLK	PXSP1110BLK	PXSP1210BLK	PXSP1310BLK			
		PXSP0810BLU	PXSP0910BLU	PXSP1010BLU	PXSP1110BLU	PXSP1210BLU	PXSP1310BLU			
		PXSP0810PUR	PXSP0910PUR	PXSP1010PUR	PXSP1110PUR	PXSP1210PUR	PXSP1310PUR			
	11"			PXSP1011BLK	PXSP1111BLK	PXSP1211BLK	PXSP1311BLK	PXSP1411BLK		
				PXSP1011BLU	PXSP1111BLU	PXSP1211BLU	PXSP1311BLU	PXSP1411BLU		
				PXSP1011PUR	PXSP1111PUR	PXSP1211PUR	PXSP1311PUR	PXSP1411PUR		
	12"			PXSP1012BLK	PXSP1112BLK	PXSP1212BLK	PXSP1312BLK	PXSP1412BLK		
				PXSP1012BLU	PXSP1112BLU	PXSP1212BLU	PXSP1312BLU	PXSP1412BLU		
				PXSP1012PUR	PXSP1112PUR	PXSP1212PUR	PXSP1312PUR	PXSP1412PUR		
	13"			PXSP1013BLK	PXSP1113BLK	PXSP1213BLK	PXSP1313BLK	PXSP1413BLK	PXSP1513BLK	
				PXSP1013BLU	PXSP1113BLU	PXSP1213BLU	PXSP1313BLU	PXSP1413BLU	PXSP1513BLU	
				PXSP1013PUR	PXSP1113PUR	PXSP1213PUR	PXSP1313PUR	PXSP1413PUR	PXSP1513PUR	
	14"			PXSP1014BLK	PXSP1114BLK	PXSP1214BLK	PXSP1314BLK	PXSP1414BLK	PXSP1514BLK	PXSP1614BLK
				PXSP1014BLU	PXSP1114BLU	PXSP1214BLU	PXSP1314BLU	PXSP1414BLU	PXSP1514BLU	PXSP1614BLU
				PXSP1014PUR	PXSP1114PUR	PXSP1214PUR	PXSP1314PUR	PXSP1414PUR	PXSP1514PUR	PXSP1614PUR
	15"					PXSP1215BLK	PXSP1315BLK	PXSP1415BLK	PXSP1515BLK	PXSP1615BLK
						PXSP1215BLU	PXSP1315BLU	PXSP1415BLU	PXSP1515BLU	PXSP1615BLU
						PXSP1215PUR	PXSP1315PUR	PXSP1415PUR	PXSP1515PUR	PXSP1615PUR
16"					PXSP1216BLK	PXSP1316BLK	PXSP1416BLK	PXSP1516BLK	PXSP1616BLK	
					PXSP1216BLU	PXSP1316BLU	PXSP1416BLU	PXSP1516BLU	PXSP1616BLU	
					PXSP1216PUR	PXSP1316PUR	PXSP1416PUR	PXSP1516PUR	PXSP1616PUR	
17"						PXSP1317BLK	PXSP1417BLK	PXSP1517BLK	PXSP1617BLK	
						PXSP1317BLU	PXSP1417BLU	PXSP1517BLU	PXSP1617BLU	
						PXSP1317PUR	PXSP1417PUR	PXSP1517PUR	PXSP1617PUR	
18"						PXSP1318BLK	PXSP1418BLK	PXSP1518BLK	PXSP1618BLK	
						PXSP1318BLU	PXSP1418BLU	PXSP1518BLU	PXSP1618BLU	
						PXSP1318PUR	PXSP1418PUR	PXSP1518PUR	PXSP1618PUR	

Accessories


Extra Cover

Extra Cover Black	\$120
Extra Cover Blue	\$120
Extra Cover Purple	\$120

Solid Seat Insert

Solid Seat Insert (E0992)	\$73
---------------------------	------


Axiom for Kids Back Width (all models)	Fits back Frame Widths
10" (25.5 cm)	9" - 11" (23 cm - 28 cm)
12" (30.5 cm)	11" - 13" (28 cm - 33 cm)
14" (35.5 cm)	13" - 15" (33 cm - 38 cm)
16" (40.5 cm)	15" - 17" (38 cm - 43 cm)



Axiom for Kids PL - Posterior Lateral Back

Not available on Ki Mobility Wheelchairs with 2" Inward Offset or 2" Outward Offset back posts. Not available on Ki Mobility Wheelchairs with Tubular Flip-Up Armrests. May have mounting conflicts with low back cane heights and Angle Adjustable Locking Extendable Flip-Up Armrests with Rigidizer Bar.

Back Width	10" (25.5 cm)	12" (30.5 cm)	14" (35.5 cm)	16" (40.5 cm)
Fits Back Frame Widths:	9"-11" (23 cm - 28 cm)	11"-13" (28 cm - 33 cm)	13"-15" (33 cm - 38 cm)	15"-17" (38 cm - 43 cm)
Actual Back Width (F):	10" (25.5 cm)	12" (30.5 cm)	14" (35.5 cm)	16" (40.5 cm)
Fits Trunk Widths (F):	10" (25.5 cm)	12" (30.5 cm)	14" (35.5 cm)	16" (40.5 cm)
Contour Depth (H):	3" (7.5 cm)	3" (7.5 cm)	3" (7.5 cm)	3" (7.5 cm)



Back

Includes Lumbar Pad. Weight capacity 165 lbs. (75 kg). Includes Axiom DEX for Kids Hardware.

		Wheelchair Frame Width			
		10	12	14	16
Retail		\$755			
HCPCS		E2615	E2615	E2615	E2615
Axiom Back Height	10	PXPL1010			
	12	PXPL1012	PXPL1212	PXPL1412	
	14	PXPL1014	PXPL1214	PXPL1414	PXPL1614
	16		PXPL1216	PXPL1416	PXPL1616
	18				PXPL1618

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.


Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1076.5 Pricing and specifications subject to change without notice.

Axiom for Kids AL - Posterior Back for Adjustable Laterals

Not available on Ki Mobility Wheelchairs with 2" Inward Offset or 2" Outward Offset back posts. Not available on Ki Mobility Wheelchairs with Tubular Flip-Up Armrests. May have mounting conflicts with low back cane heights and Angle Adjustable Locking Extendable Flip-Up Armrests with Rigidizer Bar.

Back Width	10" (25.5 cm)	12" (30.5 cm)	14" (35.5 cm)	16" (40.5 cm)
Fits Back Frame Widths:	9"-11" (23 cm - 28 cm)	11"-13" (28 cm - 33 cm)	13"-15" (33 cm - 38 cm)	15"-17" (38 cm - 43 cm)
Actual Back Width (F):	8" (20cm)	10" (25.5 cm)	12" (30.5 cm)	14" (35.5 cm)
Component Lateral Width Range (J)*: <i>(*-1" with included spacers)</i>	8" - 11" (20 cm - 28 cm)	10" - 13" (25.5 cm - 33 cm)	12" - 15" (30.5 cm - 38 cm)	14" - 17" (35.5 cm - 43 cm)
Contour Depth (H):	1.25" (3 cm)	1.25" (3 cm)	1.25" (3 cm)	1.25" (3 cm)



Back

Weight capacity 165 lbs. (75 kg). Includes Axiom DEX for Kids Hardware.

		Wheelchair Frame Width			
		10"	12"	14"	16"
Axiom Back Height	Retail	\$887			
	HCPCS	E2291	E2291	E2291	E2291
	10"	PXAL1010			
	12"	PXAL1012	PXAL1212	PXAL1412	
	14"	PXAL1014	PXAL1214	PXAL1414	PXAL1614
	16"		PXAL1216	PXAL1416	PXAL1616
	18"			PXAL1618	

Axiom for Kids - Adjustable Lateral Trunk Supports (for PXAL Back)

Not available for Posterior Back for Contoured Adjustable Laterals (PXALCL)

Description	Item Code	HCPCS	Retail
3 x 4 Lateral Trunk Pads, Flat with Fixed Mounting Hardware (Pair) (E0956)	PLF3x4PP	E0956	\$311
4 x 5 Lateral Trunk Pads, Flat with Fixed Mounting Hardware (Pair) (E0956)	PLF4x5PP	E0956	\$311
3 x 4 Lateral Trunk Pads, Flat with Swing Away Mounting Hardware (Pair) (E0956,E1034)	PLTAS3x4PP	E0956, E1034	\$806
4 x 5 Lateral Trunk Pads, Flat with Swing Away Mounting Hardware (Pair) (E0956,E1034)	PLTAS4x5PP	E0956, E1034	\$806
3 x 4 Lateral Trunk Pads, Curved with Swing Away Mounting Hardware (Pair) (E0956,E1034)	PLTASC3x4PP	E0956, E1034	\$806
4 x 5 Lateral Trunk Pads, Curved with Swing Away Mounting Hardware (Pair) (E0956,E1034)	PLTASC4x5PP	E0956, E1034	\$806

Axiom for Kids - Adjustable Lateral Pelvic Supports (for PXAL Back)


Not available for Posterior Back for Contoured Adjustable Laterals (PXALCL)

Description	Item Code	HCPCS	Retail
3 x 4 Lateral Pelvic Pads with Fixed Mounting Hardware (Pair) (E0956)	PLPF3x4PP	E0956	\$311
4 x 5 Lateral Pelvic Pads with Fixed Mounting Hardware (Pair) (E0956)	PLPF4x5PP	E0956	\$311
3 x 4 Lateral Pelvic Pads with Swing Away Mounting Hardware (Pair) (E0956,E1034)	PLASP3x4PP	E0956, E1034	\$806
4 x 5 Lateral Pelvic Pads with Swing Away Mounting Hardware (Pair) (E0956,E1034)	PLASP4x5PP	E0956, E1034	\$806

Axiom for Kids ALCL - Posterior Back for Contoured Adjustable Laterals

Not available on Ki Mobility Wheelchairs with 2" Inward Offset or 2" Outward Offset back posts. Not available on Ki Mobility Wheelchairs with Tubular Flip-Up Armrests. May have mounting conflicts with low back cane heights and Angle Adjustable Locking Extendable Flip-Up Armrests with Rigidizer Bar.

Back Width	10" (25.5 cm)	12" (30.5 cm)	14" (35.5 cm)	16" (40.5 cm)
Fits Back Frame Widths:	9"-11" (23 cm - 28 cm)	11"-13" (28 cm - 33 cm)	13"-15" (33 cm - 38 cm)	15"-17" (38 cm - 43 cm)
Actual Back Width (F):	8" - 11" (20 cm - 28 cm)	10" - 13" (25.5 cm - 33 cm)	12" - 15" (30.5 cm - 38 cm)	14" - 17" (35.5 cm - 43 cm)
Contoured Lateral Width Range (J):	8" - 11" (20 cm - 28 cm)	10" - 13" (25.5 cm - 33 cm)	12" - 15" (30.5 cm - 38 cm)	14" - 17" (35.5 cm - 43 cm)
Contour Depth (K):	5" - 5.5" (12.5 cm - 14 cm)	5" - 5.5" (12.5 cm - 14 cm)	5" - 5.5" (12.5 cm - 14 cm)	5" - 5.5" (12.5 cm - 14 cm)



Back

Requires Contoured Adjustable Lateral Trunk Support selection below. Weight capacity 165 lbs. (75 kg). Includes Axiom DEX for Kids Hardware.

		Wheelchair Frame Width			
		10"	12"	14"	16"
Retail		\$989			
HCPCS		E2293	E2293	E2293	E2293
Axiom Back Height	10"	PXAL1010CL			
	12"	PXAL1012CL	PXAL1212CL	PXAL1412CL	
	14"	PXAL1014CL	PXAL1214CL	PXAL1414CL	PXAL1614CL
	16"		PXAL1216CL	PXAL1416CL	PXAL1616CL
	18"				PXAL1618CL

Axiom for Kids - Contoured Adjustable Lateral Trunk Supports (for PXALCL Back)

Description	Item Code	HCPCS	Retail
3" Contoured Lateral Wings with Fixed Mounting Hardware (Pair) (E0956)	PCL3	E0956	\$332
4" Contoured Lateral Wings with Fixed Mounting Hardware (Pair) (E0956)	PCL4	E0956	\$332

4" Contoured Lateral Wings are not available on 10" or 12" Back Heights

Accessories

Axiom for Kids Back Accessories & Modifications

Lumbar Pad (K0108)	\$54
Privacy Flap (K0108)	\$54
<i>Privacy Flap not available for 10" Wide backs</i>	
Axiom DEX for Kids Hardware Extension (K0108)	\$97
<i>Hardware Extension adds an additional 2" of anterior depth adjustment</i>	
Axiom for Kids Adjustable Shoulder Guides (K0108)	\$99

Headrest Pads

Selection of a Headrest Pad requires a selection from Headrest Mounting Hardware section below



Contour Positioning



Tri-Pad Adjustable



Comfort Soft



TriSoft



Comfort Soft Adjustable

Headrest Pads			
Description	Item Code	HCPCS	Retail
Contour Positioning Headrest Pad, Small, 8W x 6.25H (E0955)	116919	E0955	\$190
Contour Positioning Headrest Pad, Medium, 9W x 6.5H (E0955)	116921	E0955	\$190
Contour Positioning Headrest Pad, Large, 10W x 7.5H (E0955)	116923	E0955	\$190
Tri-Pad Adjustable Headrest, Small (E0955)	116925	E0955	\$460
Tri-Pad Adjustable Headrest, Medium (E0955)	116927	E0955	\$460
Tri-Pad Adjustable Headrest, Large (E0955)	116929	E0955	\$460
Comfort Soft Headrest, 8" (E0955)	116931	E0955	\$204
Comfort Soft Headrest, 10" (E0955)	116933	E0955	\$204
Comfort Soft Headrest, 14" (E0955)	116935	E0955	\$204
Comfort Soft Headrest, 19" (E0955)	116937	E0955	\$235
TriSoft Headrest, 10" (E0955)	116939	E0955	\$204
TriSoft Headrest, 14" (E0955)	116941	E0955	\$204
Comfort Soft Adjustable Headrest, 10" (E0955)	116943	E0955	\$383
Comfort Soft Adjustable Headrest, 14" (E0955)	116945	E0955	\$383
Comfort Soft Adjustable Headrest, 19" (E0955)	116947	E0955	\$477

See Axiom Hardware Sizing Guide for full item sizing

Headrest Hardware

Selection of Headrest Hardware requires a selection from Headrest Pad section above



Std. Straight Removable Adjustable



Multi-Axis Removable Adjustable



Omnilink Multi-Axis Adjustable



Posalinc Removable Adjustable



Posalinc Flip-Back Adjustable

Headrest Hardware			
Description	Item Code	HCPCS	Retail
Standard Straight Removable Adjustable Hardware (E1033)	120038	E1033	\$198
Multi-Axis Removable Adjustable Hardware (E1033)	120039	E1033	\$295
Omnilink Multi-Axis Adjustable Hardware (E1033)	123603	E1033	\$417
POSALinc Removable Adjustable Hardware (E1033)	120040	E1033	\$359
POSALinc Removable Adjustable Hardware, Tall (E1033)	123528	E1033	\$396
POSALinc Flip-Back Adjustable Hardware (E1033)	123530	E1033	\$430
POSALinc Flip-Back Adjustable Hardware, Tall (E1033)	123531	E1033	\$474

Universal Headrest Mounting Plate
 Not required with Axiom Headrests

\$54

Anterior Trunk Supports



2-Piece
Hook & Loop



1-Piece
with Buckle



Dynaform Standard,
Dynamic



Dynaform Slim Cut,
Dynamic



Airlogic Standard,
Structure



Airlogic Slim Cut,
Dynamic

Anterior Trunk Supports

Anterior Trunk Supports			
Description	Item Code	HCPCS	Retail
Dynaform Chest Support: 2-Piece with Hook&Loop, 3" High Pad, Small	14941	E0960	\$146
Dynaform Chest Support: 2-Piece with Hook&Loop, 4" High Pad, Medium	14942	E0960	\$146
Dynaform Chest Support: 2-Piece with Hook&Loop, 5" High Pad, Large	14943	E0960	\$146
Dynaform Chest Support: 1-Piece with Buckle, 3" Height x 18" Width, Small	14924	E0960	\$146
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 20" Width, Medium	14927	E0960	\$146
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 22" Width, Large	14928	E0960	\$146
Dynaform Anterior Trunk Support: Standard, Dynamic, X-Small	14851	E0960	\$179
Dynaform Anterior Trunk Support: Standard, Dynamic, Small	14852	E0960	\$179
Dynaform Anterior Trunk Support: Standard, Dynamic, Medium	14853	E0960	\$179
Dynaform Anterior Trunk Support: Standard, Dynamic, Large	14854	E0960	\$179
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, X-Small	14856	E0960	\$179
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Small	14857	E0960	\$179
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Medium	14858	E0960	\$179
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Large	14859	E0960	\$179
AirLogic Anterior Trunk Support: Standard, Structure, Small	14282	E0960	\$203
AirLogic Anterior Trunk Support: Standard, Structure, Medium	14283	E0960	\$203
AirLogic Anterior Trunk Support: Standard, Structure, Large	14284	E0960	\$203
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Small	14257	E0960	\$203
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Medium	14258	E0960	\$203

Gloss Finish Back Shell

If no selection is made, standard Shadow Black will be used. Touch-up paint is available for purchase through the parts manual.

Black	N/C
Ruby Red	N/C
Radiant Orange	N/C
Yellow	N/C
Electric Green	N/C
Rich Teal	N/C
Glacier Blue	N/C
Blue Steel	N/C
Deep Sapphire Blue	N/C
Purple Plum	N/C
Opalescent Violet	N/C
Dark Pink	N/C
Light Pink	N/C
Champagne	N/C

Matte Finish Back Shell

If no selection is made, standard Shadow Black will be used. Touch-up paint is available for purchase through the parts manual.

Iridescent Navy	N/C
Bright Red	N/C
Tactical Green	N/C
Dry Earth Tan	N/C
Platinum	N/C
White	N/C