

12-Port Drying Chamber



View online: <https://www.labconco.com/product/12-port-drying-chamber/3127>

Overview



Drying Chambers provide a means to attach flasks to your freeze dryer or contain samples in bulk.

The 12-Port Drying Chamber includes 12 neoprene valves with plastic knobs that accommodate both 1/2" and 3/4" adapters for connection of flasks, a 1/2" thick acrylic lid with neoprene gasket, and one latching pin. Made of type 304 stainless steel. It is compatible with FreeZone 1, 2.5, 4.5, 6, 8, 12 or 18 Liter Freeze Dryers.

Catalog Number: 7522800

Specifications

- **Weight:** 11.0 lbs
- **Weight metric:** 5.0 kg
- **Dimensions:** 8.0" diameter x 9.75" h
- **Applicable Products:** Freeze Dryers

For Product Assistance



[Andy Miller](#)
Senior Sales Representative



[Vivian Goldmon](#)
Customer Service Representative



Labconco Corporation
8811 Prospect Avenue
Kansas City, MO 64132-2696
Phone: (816) 333-8811

