


REPLACEMENT PARTS

FR587WD DELUXE CLINICAL CARE WITH
DROP ARM - WIDE



21



22



23

Heat and Massage Feature —
Control Module and Transformer Unit
Showing All Connections Correctly Made



24



25

Item #	UOM	Description
5792A001G	1 ea	Table, Side
5792A002	4 ea / set	Hardware, Clevis Pin w Washer, Side Table
5792A003	10 ea / set	Hardware, Clevis Pin w Cotter Clip
566DG002W	1 ea	Drop Arm Assembly, Left
566DG003W	1 ea	Drop Arm Assembly, Right
566DG05100	1 ea	Hardware, Pull Pin, Drop Arm (not shown)
58737J(color code)FR	1 ea	Ottoman, Upholstered
58735J(color code)FR	1 ea	Armrest, Upholstered
58733J(color code)FR	1 ea	Footrest, Upholstered
58734J(color code)FR	1 ea	Legrest, Upholstered
58731J(color code)FR	1 ea	Seat, Upholstered
58731JH(color code)FR	1 ea	Seat, Upholstered, Heat & Massage / Heat Only
56605J(color code)FR	1 ea	Side Panels, Upholstered
58732J(color code)FR	1 ea	Back, Upholstered
58732JH(color code)FR	1 ea	Back, Upholstered, Heat & Massage

Item #	UOM	Description
58732JH0(color code)FR	1 ea	Back, Upholstered, Heat Only
AP020004-8	8 ea / set	Hardware, Push Pin & Spring, Side Panels
CSTRTRGN01	1 ea	Caster, Directional Green Tab Tente
CSTRBLK01	1 ea	Caster, Swivel, Locking Black Tab Tente
5873C400A	1 ea	Push Handle w Grip
9205	1 ea	Grip Only, Push Handle
5610F200	1 ea	Mechanism, Recliner
5770A073	1 ea	Hardware, Gas Spring Bracket
EI000001	1 ea	Hardware, Gas Spring
LRHM0114	1 ea	Pendant Holster, Heat & Massage
HMO0G3CORD	1 ea	Powercord Supply, Heat & Massage Gen 3
ECMCGF2	1 ea	Control Box, Heat & Massage Gen 3
HMO0G3PEND	1 ea	Pendant, Heat & Massage Gen 3
HMO0G3PENDHO	1 ea	Pendant, Heat Only Gen 3
5772A001	1 ea	Hardware, Cable Assemble (not shown)
LRMP0001	1 ea	Paint, Touch-Up, Gray (not shown)



GF Health Products, Inc. is an
ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2023, GF Health Products, Inc. All Rights Reserved.
Lumex is a trademark of GF Health Products, Inc.
GF Health Products, Inc. is an ISO 13485 Certified Company.

