

Features:

- Sturdy 1" aluminum tubing provides maximum strength while remaining lightweight
- Dual-release folding mechanism enables user to fold walker sides independently
- Accommodates all leg attachments
- Designed to meet the requirements of HCPCS code: E0135
- One (1) year limited warranty

Specifications:

| | <u>Adult</u> | <u>Junior</u> |
|---|--------------|---------------|
| Width Inside Handgrips | 16 ½" | 16 ½" |
| Width Inside Rear Legs | 19 ½" | 19" |
| Overall Width* | 19 ½" | 19 ½" |
| Depth Front to Rear* | 15" | 15" |
| Depth Folded | 4" | 4" |
| Minimum Height | 32 ¼" | 25 ½" |
| Maximum Height | 39 ¼" | 34 ¼" |
| Weight Capacity - EVENLY DISTRIBUTED | 300 lb | 300 lb |



| | | |
|------------------|--------|---------|
| 716070A-2 | Adult | 2 ea/cs |
| 716070A-4 | Adult | 4 ea/cs |
| 716170A-4 | Junior | 4 ea/cs |