

Disposable Ear Curettes


McKesson Disposable Ear Curettes are made of smooth, bendable polypropylene, and are designed as a single-use alternative to traditional, reusable ear curettes made of stainless steel. Available in five different tip styles to accommodate a variety of patients, including adult, child and infant sizes. Disposable Ear Curettes are individually packaged, designed for single-use, and color-coded for easy product identification.

Features and Benefits

- Designed to flex and bend, increasing patient comfort
- Disposable to reduce the risk of cross-contamination and the need to autoclave
- Available in a variety of sizes
- Color-coded for easy identification
- Not made with natural rubber latex



Infant Scoop Tip



Child Tip



Spoon Tip



Round Tip



Oval Tip

ORDERING INFORMATION

DESCRIPTION	MFR #	QTY
Ear Curette, Gray, Infant Scoop Tip, 2mm	532	50ea/bx, 10bx/cs
Ear Curette, Pink, Child Tip, 3mm	529	50ea/bx, 10bx/cs
Ear Curette, Blue, Spoon Tip, 4mm	520	50ea/bx, 10bx/cs
Ear Curette, Green, Round Tip, 4mm	530	50ea/bx, 10bx/cs
Ear Curette, White, All-purpose, Oval Tip, 4mm	531	50ea/bx, 10bx/cs

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them – making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.