

ONE-CARE®

NON CONTACT INFRARED
BODY THERMOMETER



For detecting body temperature of patients of all ages, in-home and clinical environment...



REF
4020
HTD8813C



REF
4021
HTD8818A

NON - CONTACT

Measure body temperature without touching the ear or the forehead¹

CONVENIENT

For patients of all ages, in-home and clinical environment¹

COMFORTABLE

High definition LCD display with easy to read measurements outcome (even in the dark)¹

ACCURATE READING

Equipped with advanced infrared technology and high-precision sensor¹

www.one-care.com

ONE-CARE®
NON CONTACT INFRARED
BODY THERMOMETER




Product features:

- **Non-contact:** allows to measure body temperature without touching the forehead (or the ear), preventing cross-infection
- **Accurate Reading:** equipped with advanced infrared technology and high precision sensor, provides instant reading, which only takes 1 second to read the temperature. The accuracy of temperature measurement is within 0.2°C.
- **Intelligent:** enables storing up to 50 measurements for model HTD8813, and up to 10 measurements for model HT8818
- **Silent mode:** allows to take the temperature without disturbing the baby's sleep.
- **Adaptable:** the device has two units - Fahrenheit °F and Celsius °C
- **Color backlight representing the measurement results:**
 - Green** - normal body temperature;
 - Yellow** - mild fever;
 - Red** - high fever (accompanied by rapid alarm sound).

Instruction For Use

- Turn on the thermometer by pressing the On/Scan button.
- Align the sensor with the middle of your forehead, and then press the On/Scan button to start the measurement.



Ordering Information

| REF | PRODUCT | MODEL | PACKAGING |
|------|--|----------|----------------------------|
| 4020 | ONE-CARE® Non-Contact Infrared Thermometer | HTD8813C | 1 pc / box 40 bx / case |
| 4021 | ONE-CARE® Non-Contact Infrared Thermometer | HTD8818A | 1 pc / box 40 bx / case |

Sources:
1. Data on file

www.one-care.com

20241006_V1