

# ONE-CARE<sup>®</sup>

BLOOD COLLECTION TUBE HOLDERS



## Collection Holders

Standard Tube Holder



Standard holder for regular tube

Wide Holder



Holder for wider tube

Access Holder with luer adapter



Male luer adapter pre-attached

Blood Transfer Holder



Female luer adapter pre-attached

Wide Neck Transfer Holder



Male Luer adapter pre-attached  
Fit for regular and wider tube

Wide Neck Access Holder



Female luer adapter pre-attached  
Fit for regular and wider tube

Safety Holder



Safety Holder

All holder are designed to resistant to isopropyl alcohol, no breaking in connection.

## Features and Benefits

- Compatible with standard blood collection needles.
- All Materials are latex &DEHP free.

# ONE-CARE®

BLOOD COLLECTION TUBE HOLDERS

## Collection Holders

### Package Configuration

Caina Code	Description of goods	Individual Pack	Packing Config	Shipper Case
			ea/bx/cs	
5015	Standard Tube Holder	Bulk Pack	250/1000	450*420*210mm
5015-W	Wide Holder	Bulk Pack	250	400*285*290mm
5016	Safety Holder	Bulk Pack	500/500	435*365*180mm
5017	Blood Transfer Holder	Soft Pack	1/200/800	470*430*320mm
5017-W	Wide Neck Transfer Holder	Soft Pack	1/60/240	470*430*320mm
5018	Access Holder with luer adapter	Soft Pack	1/200/800	470*430*320mm
5018-W	Wide Neck Access Holder	Soft Pack	1/60/240	470*430*320mm

### Refer Standards

ISO 23908: 2011 sharps injury protection-requirements and test methods-sharps protection features for single-use hypodermic needles, introducers for catheters and needles used for blood sampling.

Biocompatibility per ISO 10993-5,7,10,11.

### Regulatory Compliance

US FDA 510K: K200027.

CE certificate MDD Annex V.

### Shelf Life

5 years



# www.one-care.com

**MediVena®**  
inspired by innovation

[www.medivena.com](http://www.medivena.com)

AMERICAS:  
MediVena  
1 Glenlake Parkway, Ste 700, Atlanta, GA 30328, USA  
+1 404 514 2586 [info@medivena.com](mailto:info@medivena.com)  
**CUSTOMER SERVICE TOLL FREE: 888 270 6457**

EUROPE:  
MediVena Sp. z o.o.  
Al. Jana Pawła II 27, 00-867 Warsaw, POLAND  
+48 695 995700