

MediVena[®]

INSPIRED BY INNOVATION



PRODUCT CATALOG

2026

MediVena®

Who are we?



MediVena is committed to bringing innovative technologies to the market, that improve patient care and the daily work of healthcare professionals. Supported by a team of experienced professionals with deep expertise in the medical device industry, our mission is to develop and deliver innovative medical solutions for blood collection, rapid diagnostics, fluid injections, and personal protection that are safe, effective, affordable, and comfortable for both healthcare providers and patients worldwide.

As the medical device market evolves, MediVena remains at the forefront by continuously expanding our portfolio and seeking out new, innovative solutions to meet emerging healthcare needs.

Sales and Marketing

We offer a comprehensive product portfolio and a unique value proposition tailored to various market segments, including Acute Care, Long-Term Care, Physician Offices, Blood Banks, Home Health, and Animal Health across the US and international markets.

Design and Innovation

MediVena places a strong emphasis on customer and market insights, driving innovation from ideation to commercialization. Our capabilities include first concept development, modeling, fast prototyping, design verification and validation, clinical trials, and seamless implementation into production.

Manufacturing & Supply Chain

Our best-in-class manufacturing facilities, supply chain, and skilled workforce span the US, Europe, Latin America, and Asia, ensuring product leadership, operational excellence, and adherence to the highest regulatory standards

MediVena has established itself as a trusted and reliable partner in the healthcare industry, with a strong presence in the Americas and Europe. Our products are distributed through partnerships with primary care providers, hospitals, long-term care facilities, blood banks, reference labs, physician offices, and the home care market. We are committed to quality and customer service, offering reliable products and prompt delivery, with orders shipped within 48 hours from our Atlanta warehouse and complimentary ground shipping for our partners.

Our Brands

One-Care[®]

SAFETY MEDICAL DEVICES

ONE-CARE[®] is our flagship brand, covering a wide range of products across several key medical fields, including blood micro-sampling, insulin delivery, blood collection, glucose monitoring, syringes and needles, and personal protective equipment. Our goal is to provide both healthcare professionals and lay users with effective, affordable, safe, and comfortable solutions.

Our products are widely distributed to various sectors of the healthcare industry, including primary care, hospitals, long-term care facilities, physician offices, and others, as well as the home care market.

We are systematically enlarging our portfolio in response to consumer needs and expectations. ONE-CARE[®] brand offers a variety of choices tailored to meet the market's specific needs. We are dedicated to delivering solutions for every requirement, ensuring that our users have access to the best possible care and support.

Gluco-Care[®]

DIABETES CARE SOLUTIONS

The GLUCO-CARE[®] brand was created with the mission of providing comprehensive support and innovative solutions to people living with diabetes. Our aim is to enhance their comfort, safety, and overall quality of life.

We focus on several key areas to achieve this:

- Supporting Accurate Diagnosis: Ensuring that individuals receive precise and timely diagnoses, which is critical for effective diabetes management.
- Assisting in Achieving Treatment Goals: Offering products and services that help individuals manage their condition and reach their treatment objectives.
- Preventing Diabetes Complications: Providing tools and resources designed to minimize the risk of complications, promoting long-term health and well-being.

GLUCO-CARE[®] is dedicated to empowering people with diabetes to lead healthier and more fulfilling lives through our commitment to quality, innovation, and patient-centric care.

Contents

Blood Microsampling

- Safety Lancets..... 6

Blood Collection

- Blood Collection Needles..... 13
- Safety Blood Collection Winged Sets..... 15
- Tube Holders..... 18

Diagnostics

- Glucose Monitors..... 20
- Thermometers..... 22

Drug Delivery

- Syringes and Safety Needles..... 24
- Luer-Lock Syringes..... 26
- Safety Insulin Syringes..... 27
- Insulin Syringes..... 29
- Safety Pen Needles..... 35
- Pen Needles..... 39

Feed and Saline Flush Syringes

- Oral Syringes..... 42
- Enteral Syringes..... 43
- Saline Flush Syringes..... 44

Skin Preparation Solutions

- Alcohol Prep Pads..... 46

USERS:



For Professional use only



For Home and Personal use



Blood Microsampling



Blood Collection



Diagnostics



Drug Delivery



Feed and Saline Flush Syringes

SAFETY

Blood Microsampling



Our Blood Microsampling line is designed to offer safe, reliable, and efficient capillary blood collection. Featuring single-use safety lancets -both contact-activated and push-button activated - Medivena's solutions help protect healthcare workers while ensuring patient comfort and sampling precision.

Gluco-Care®

SAFETY LANCETS

BLOOD MICROSAMPLING Safety Lancets

USERS:



PRODUCT FEATURES

Safe

- Single-use mechanism prevents sharps injuries and reuse of the device.
- Fully enclosed housing ensures that needle is hidden before and after use.
- All lancets are gamma-sterilized and protected with a sterility tab.

Comfortable

- Silicone-coated, ultra-thin needles ideally positioned during skin penetration with double-spring mechanism to reduce vibrations.
- Color-coded product versions to meet specific blood sample requirements and address a variety of patient skin types.
- Ergonomic design for comfortable grip and secure blood collection.

Simple

- Easy to use, intuitive operation in 3 steps only, no pre-loading required.
- Wide safety cap designed for easy removal - simply twist off the cap and pull it out.
- Optimal skin contact area for precise positioning on the desired test site and adequate blood flow.



ORDERING INFORMATION

REF	Diameter (G) x Depth (mm/in)	Color Coding	Packaging
1111	30G x 1.8mm		100 pcs/box 50 bx (5,000 pcs)/case
1112	28G x 1.8mm		100 pcs/box 50 bx (5,000 pcs)/case
1113	26G x 1.8mm		100 pcs/box 40 bx (4,000 pcs)/case
1114	21G x 1.8mm		100 pcs/box 50 bx (5,000 pcs)/case
1115	21G x 2.4mm		100 pcs/box 50 bx (5,000 pcs)/case
1116	21G x 2.8mm		100 pcs/box 50 bx (5,000 pcs)/case
1117	23G x 1.8mm		100 pcs/box 50 bx (5,000 pcs)/case
1118	23G x 2.4mm		100 pcs/box 50 bx (5,000 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

USERS:



T-Shaped body
for a secure grip.



PRODUCT FEATURES

Safe:

- All lancets are gamma-sterilized and protected with a sterility tab.
- Fully enclosed housing ensures that needle is hidden before and after use.
- Single-use mechanism prevents sharps injuries and reuse of the device.

Comfortable:

- Silicone-coated, ultra-thin needles ideally positioned during skin penetration with double-spring mechanism to reduce vibrations.
- Color-coded product versions to meet specific blood sample requirements and address a variety of patient skin types.
- Ergonomic design for comfortable grip and secure blood collection.

Simple:

- Easy to use, intuitive operation in 3 steps only, no pre-loading required.
- Long safety tab designed for easy removal - simply twist off the tab and pull it out.
- Optimal skin contact area for precise positioning on the desired test site and adequate blood flow.

ORDERING INFORMATION | BOX 100 PCS

REF	Diameter (G) x Depth (mm/in)	Color Coding	Packaging
1041	30G x 0.7mm		100 pcs/box 20 bx (2,000 pcs)/case
1042	28G x 1.0mm		100 pcs/box 20 bx (2,000 pcs)/case
1043	25G x 1.1mm		100 pcs/box 20 bx (2,000 pcs)/case
1044	21G x 1.8mm		100 pcs/box 20 bx (2,000 pcs)/case
1045	25G x 2.0mm		100 pcs/box 20 bx (2,000 pcs)/case
1046	21G x 2.4mm		100 pcs/box 20 bx (2,000 pcs)/case
1047	21G x 2.8mm		100 pcs/box 20 bx (2,000 pcs)/case
1048	23G x 1.8mm		100 pcs/box 20 bx (2,000 pcs)/case
1049	23G x 2.4mm		100 pcs/box 20 bx (2,000 pcs)/case

ORDERING INFORMATION | BOX 200 PCS

REF	Diameter (G) x Depth (mm/in)	Color Coding	Packaging
1031	30G x 1.5mm		200 pcs/box 18 bx (3,600 pcs)/case
1032	28G x 1.6mm		200 pcs/box 18 bx (3,600 pcs)/case
1033	25G x 1.6mm		200 pcs/box 18 bx (3,600 pcs)/case
1034	23G x 1.8mm		200 pcs/box 18 bx (3,600 pcs)/case
1035	21G x 1.8mm		200 pcs/box 18 bx (3,600 pcs)/case
1036	21G x 2.4mm		200 pcs/box 18 bx (3,600 pcs)/case
1037	21G x 2.8mm.		200 pcs/box 18 bx (3,600 pcs)/case
1038	17G BLADE x 1.0mm		200 pcs/box 18 bx (3,600 pcs)/case
1039	17G BLADE x 1.5mm		200 pcs/box 18 bx (3,600 pcs)/case
1040	17G BLADE x 2.0mm		200 pcs/box 18 bx (3,600 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

USERS:



T-Shaped body
for a secure grip.



PRODUCT FEATURES

Enhanced Safety:

- Single-use mechanism prevents reuse and minimizes the risk of sharps injuries.
- Fully enclosed needle before and after activation ensures maximum user and patient protection.
- Gamma-sterilized lancets with a protective sterility tab support safe and hygienic use..

Precision and comfort:

- Ultra-thin, silicone-coated needles combined with a controlled double-spring mechanism reduce vibration and help minimize pain during skin penetration.
- A dedicated activation button ensures reliable, consistent activation with controlled force, supporting precise puncture depth and dependable blood sampling.
- Ergonomic, T-shaped design provides a secure grip, stable positioning, and confident handling.
- Color-coded options allow easy selection based on blood sampling requirements and patient skin type.

Simple and intuitive use:

- Intuitive, three-step operation with no pre-loading required.
- Wide safety cap designed for easy removal with a simple twist-off motion.
- Optimized skin contact area supports accurate positioning and consistent blood flow.

ORDERING INFORMATION | BOX 100 PCS

REF	Diameter (G) x Depth (mm/in)	Color Coding	Packaging
1011	30G x 1.5mm		100 pcs/box 20 bx (2,000 pcs)/case
1012	28G x 1.5mm		100 pcs/box 20 bx (2,000 pcs)/case
1013	26G x 1.6mm		100 pcs/box 20 bx (2,000 pcs)/case
1014	21G x 1.8mm		100 pcs/box 20 bx (2,000 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

USERS:



PRODUCT FEATURES

Designed for Precision Detection

- Stable puncture depth ensures accurate and consistent blood sampling.
- Reduces skin trauma and minimizes tissue fluid contamination.

Safe & Sterile

- Fully enclosed needle for maximum safety – automatically retracts after use.
- Prevents reuse and accidental needle sticks.
- Gamma sterilized and chemical-residue free. Meets OSHA's safety standards.

Comfort & Convenience

- Side push-button activation for quick, easy operation.
- Wide, non-slip button allows for precise and easy activation.
- Ergonomic, narrowing design for precise blood collection.
- Two-step intuitive operation.
- Available in multiple gauge and depth options for different testing needs.

ORDERING INFORMATION

REF	Diameter (G) x Depth (mm/in)	Color Coding	Packaging
1061	30G x 1.8mm		100 pcs/box 20bx (2,000 pcs)/case
1062	28G x 1.8mm		100 pcs/box 20bx (2,000 pcs)/case
1063	26G x 1.8mm		100 pcs/box 20bx (2,000 pcs)/case
1064	21G x 1.8mm		100 pcs/box 20bx (2,000 pcs)/case
1065	21G x 2.2mm		100 pcs/box 20bx (2,000 pcs)/case
1066	21G x 2.8mm		100 pcs/box 20bx (2,000 pcs)/case
1067	17G x 2.0mm		100 pcs/box 20bx (2,000 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

USERS:



PRODUCT FEATURES

SAFE

- Single-use mechanism prevents sharps injuries and device reuse.
- Fully enclosed housing keeps the needle safely covered before and after use.
- Lancets are sterilized by irradiation and protected with a sterility tab

CONVENIENT

- Color-coded versions support different blood sample requirements and accommodate various skin types.
- Non-slip body ensures a comfortable grip and secure blood collection.
- One-hand activation allows for quick and easy testing.

SIMPLE

- Intuitive, easy-to-use design with no pre-loading required.
- Wide safety cap can be easily removed
- Optimized skin contact area ensures precise positioning and an adequate blood sample.



ORDERING INFORMATION

REF	Diameter (G) x Depth (mm/in)	Color Coding	Packaging
1071	30G x 1.5mm		100 pcs/box 20bx (2,000 pcs)/case
1072	28G x 1.8mm		100 pcs/box 20bx (2,000 pcs)/case
1073	26G x 1.8mm		100 pcs/box 20bx (2,000 pcs)/case
1074	21G x 1.8mm		100 pcs/box 20bx (2,000 pcs)/case
1075	21G x 2.4mm		100 pcs/box 20bx (2,000 pcs)/case
1076	21G x 2.8mm		100 pcs/box 20bx (2,000 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact Medivena directly.

SAFETY

Blood Collection



Our Blood Collection portfolio is designed to meet the highest standards of safety and efficiency in clinical settings. It includes blood collection needles engineered for smooth venipuncture, and safety winged sets that combine patient comfort with needlestick injury prevention. These products support accurate diagnostic results while promoting healthcare worker and patient safety.




One-Care[®]

SAFETY BLOOD COLLECTION NEEDLES

BLOOD COLLECTION

Safety Blood Collection Needles

USERS:



- ✓ One-hand operation - thumb or surface activation.
- ✓ Audible “click” to indicate that safety shield is activated and the needle is permanently locked.

PRODUCT FEATURES

Widely Accepted Safety Solution

- The side-shield flipped over the needle - minimal training required.

Universal Luer

- Allows for universal and tight fit to tube holders and multi-samplig.

Ultra-Sharp, lubricated needle

- For patient comfort.

Ergonomic Design

- Concave shield to fit thumb shape and textured surface to improve grip.

High Fill Rate

- Thin-Wall Technology for optimal flow and quick blood collection.



Transparent Tube Holder

- Improved visibility of collection tubes filling up with blood.
- Provided separately or pre-attached.

ORDERING INFORMATION

REF	Description	Color Coding	Packaging
5020	Safety Blood Collection Needles - without holder 22G x 1 1/4"		48 pcs/box 50 bx (2,400 pcs)/case
5021	Safety Blood Collection Needles - without holder 21G x 1 1/4"		48 pcs/box 50 bx (2,400 pcs)/case
5022	Safety Blood Collection Needles - without holder 20G x 1 1/4"		48 pcs/box 50 bx (2,400 pcs)/case
5020-H	Safety Blood Collection Needles - with holder 22G x 1 1/4"		48 pcs/box 50 bx (2,400 pcs)/case
5021-H	Safety Blood Collection Needles - with holder 21G x 1 1/4"		48 pcs/box 50 bx (2,400 pcs)/case
5022-H	Safety Blood Collection Needles - with holder 20G x 1 1/4"		48 pcs/box 50 bx (2,400 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact Medivena directly.

One-Care[®]

SAFETY BLOOD COLLECTION WINGED SET **PUSH BUTTON**

BLOOD COLLECTION

Safety Blood Collection Winged Sets

USERS:



PUSH BUTTON AUTO-RETRACTABLE

Before use

After use

PRODUCT FEATURES

Thin-Wall Technology

- For optimal flow and higher tube fill rates.

Universal Luer

- Allows for universal and tight fit to tube holders and multi-samplig.

Tri-Bevel

- Electropolished, Ultra-sharp and lubricated needle.

DEPF-Free

- Low memory and flexible tubing assists handling and preserves blood composition.

PUSH BUTTON AUTO-RETRACTABLE

ONE-HAND OPERATION

Provides the ultimate protection for healthcare professionals with push button and automated needle retraction technology. It allows for one-hand activation of safety shield in-vein or immediately after the blood draw.

ORDERING INFORMATION



REF	Description	Color Coding	Packaging
5070	PUSH BUTTON Auto-Retractable Safety Blood Collection Winged Sets 25G x 3/4", Luer-Adapter, 12" Tube	Orange	50 pcs/box, 4bx/inner pack 16 bx (800 pcs)/case
5071	PUSH BUTTON Auto-Retractable Safety Blood Collection Winged Sets 23G x 3/4", Luer-Adapter, 12" Tube	Blue	50 pcs/box, 4bx/inner pack 16 bx (800 pcs)/case
5072	PUSH BUTTON Auto-Retractable Safety Blood Collection Winged Sets 21G x 3/4", Luer-Adapter, 12" Tube	Green	50 pcs/box, 4bx/inner pack 16 bx (800 pcs)/case
5070-H	PUSH BUTTON Auto-Retractable Safety Blood Collection Winged Sets 25G x 3/4", Pre -attached tube holder, 12" Tube	Orange	20 pcs/box, 5bx/inner pack 20 bx (400 pcs)/case
5071-H	PUSH BUTTON Auto-Retractable Safety Blood Collection Winged Sets 23G x 3/4", Pre -attached tube holder, 12" Tube	Blue	20 pcs/box, 5bx/inner pack 20 bx (400 pcs)/case
5072-H	PUSH BUTTON Auto-Retractable Safety Blood Collection Winged Sets 21G x 3/4", Pre -attached tube holder, 12" Tube	Green	20 pcs/box, 5bx/inner pack 20 bx (400 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

One-Care®


SAFETY BLOOD COLLECTION WINGED SET

GLIDE & LOCK

BLOOD COLLECTION

Safety Blood Collection Winged Sets

USERS:



GLIDE & LOCK SAFETY SHIELD



PRODUCT FEATURES

Thin-Wall Technology

- For optimal flow and higher tube fill rates.

Universal Luer

- Allows for universal and tight fit to tube holders and multi-samplig.

Tri-Bevel

- Electropolished, Ultra-sharp and lubricated needle.

DEPF-Free

- Low memory and flexible tubing assists handling and preserves blood composition.

ORDERING INFORMATION

REF	Description	Color Coding	Packaging
5040	GLIDE & LOCK Safety Blood Collection Winged Sets 25G x 3/4", Luer-Adapter, 12" Tube		50 pcs/box 20 bx (1,000 pcs)/case
5041	GLIDE & LOCK Safety Blood Collection Winged Sets 24G x 3/4", Luer-Adapter, 12" Tube		50 pcs/box 20 bx (1,000 pcs)/case
5042	GLIDE & LOCK Safety Blood Collection Winged Sets 23G x 3/4", Luer-Adapter, 12" Tube		50 pcs/box 20 bx (1,000 pcs)/case
5043	GLIDE & LOCK Safety Blood Collection Winged Sets 22G x 3/4", Luer-Adapter, 12" Tube		50 pcs/box 20 bx (1,000 pcs)/case
5044	GLIDE & LOCK Safety Blood Collection Winged Sets 21G x 3/4", Luer-Adapter, 12" Tube		50 pcs/box 20 bx (1,000 pcs)/case
5045	GLIDE & LOCK Safety Blood Collection Winged Sets 20G x 3/4", Luer-Adapter, 12" Tube		50 pcs/box 20 bx (1,000 pcs)/case
5046	GLIDE & LOCK Safety Blood Collection Winged Sets 19G x 3/4", Luer-Adapter, 12" Tube		50 pcs/box 20 bx (1,000 pcs)/case
5040-H	GLIDE & LOCK Safety Blood Collection Winged Sets 25G x 3/4", Pre-Attached Holder, 12" Tube		25 pcs/box 20 bx (500 pcs)/case
5041-H	GLIDE & LOCK Safety Blood Collection Winged Sets 24G x 3/4", Pre-Attached Holder, 12" Tube		25 pcs/box 20 bx (500 pcs)/case
5042-H	GLIDE & LOCK Safety Blood Collection Winged Sets 23G x 3/4", Pre-Attached Holder, 12" Tube		25 pcs/box 20 bx (500 pcs)/case
5043-H	GLIDE & LOCK Safety Blood Collection Winged Sets 22G x 3/4", Pre-Attached Holder, 12" Tube		25 pcs/box 20 bx (500 pcs)/case
5044-H	GLIDE & LOCK Safety Blood Collection Winged Sets 21G x 3/4", Pre-Attached Holder, 12" Tube		25 pcs/box 20 bx (500 pcs)/case
5045-H	GLIDE & LOCK Safety Blood Collection Winged Sets 20G x 3/4", Pre-Attached Holder, 12" Tube		25 pcs/box 20 bx (500 pcs)/case
5046-H	GLIDE & LOCK Safety Blood Collection Winged Sets 19G x 3/4", Pre-Attached Holder, 12" Tube		25 pcs/box 20 bx (500 pcs)/case



*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

One-Care®

SAFETY BLOOD COLLECTION WINGED SET **PRESS & LOCK**

BLOOD COLLECTION

Safety Blood Collection Winged Sets

USERS:



PRESS & LOCK SAFETY SHIELD



PRESS & LOCK SAFETY SHIELD

Allows for activation of safety shield in-vein or immediately after the blood draw

- Press the side buttons and pull back to retract the needle into the shield.
- Audible click confirms that needle is securely locked.

PRODUCT FEATURES

Thin-Wall Technology

- For optimal flow and higher tube fill rates.

Universal Luer

- Allows for universal and tight fit to tube holders and multi-samplig.

Tri-Bevel

- Electropolished, Ultra-sharp and lubricated needle.

DEPF-Free

- Low memory and flexible tubing assists handling and preserves blood composition.

ORDERING INFORMATION

REF	Description	Color Coding	Packaging
5060	PRESS & LOCK Safety Blood Collection Winged Sets 25G x 3/4", Luer-Adapter, 12" Tube		50 pcs/box 20 bx (1,000 pcs)/case
5061	PRESS & LOCK Safety Blood Collection Winged Sets 23G x 3/4", Luer-Adapter, 12" Tube		50 pcs/box 20 bx (1,000 pcs)/case
5062	PRESS & LOCK Safety Blood Collection Winged Sets 21G x 3/4", Luer-Adapter, 12" Tube		50 pcs/box 20 bx (1,000 pcs)/case
5060-H	PRESS & LOCK Safety Blood Collection Winged Sets 25G x 3/4", Pre-Attached Holder, 12" Tube		25 pcs/box 20 bx (500 pcs)/case
5061-H	PRESS & LOCK Safety Blood Collection Winged Sets 23G x 3/4", Pre-Attached Holder, 12" Tube		25 pcs/box 20 bx (500 pcs)/case
5062-H	PRESS & LOCK Safety Blood Collection Winged Sets 21G x 3/4", Pre-Attached Holder, 12" Tube		25 pcs/box 20 bx (500 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

One-Care®

TUBE HOLDERS

BLOOD COLLECTION Tube Holders

USERS:



Access Holder



Blood Transfer Device



Standard Tube Holder



Safety Holder



Blood Collection Tube Holder
with Male Luer Adapter



Blood Transfer Device
with Female Luer Adapter

PRODUCT FEATURES

- Size: 0.9 x 19mm.
- Transparent material for better visibility of collection tubes filling up with blood.
- Universal design for compatibility with blood collection needles and winged sets.
- Tight fit connectors to prevent leakage.
- DEHP-free and Latex-free materials.

Regulatory Compliance

- US FDA, CE, ISO 13485.

Materials

- Medical grade PP,
- SUS304,
- ABS,

ORDERING INFORMATION

REF	Description	Color Coding	Packaging
5015-250	ONE-CARE® Tube Holder	Transparent	250 pcs/box 4 bx (1,000 pcs)/case
5016	ONE-CARE® Safety Tube Holder	Transparent	50/box 20 bx (1,000 pcs)/case
5017	ONE-CARE® Blood Collection Tube Holder with Female Luer Adapter	Transparent/Pink	100 pcs/box 12 bx (1,200 pcs)/case
5018	ONE-CARE® Blood Collection Tube Holder with Male Luer Adapter	Transparent/Blue	100 pcs/box 12 bx (1,200 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

Diagnostics



Our Diagnostics range features essential tools for everyday clinical and home health monitoring. Glucose monitors deliver fast, reliable blood glucose readings to support diabetes management, while our digital thermometers provide accurate body temperature measurements for a variety of patient populations. All devices are user-friendly and built for consistent performance.



One-Care®

GLUCOSE MONITOR & STRIPS
FOR PROFESSIONAL USE
& SELF TESTING KIT

DIAGNOSTICS

Blood Glucose Monitors

USERS:



4001 ONE-CARE® PRO Glucose Monitor for Professional Use

PRODUCT FEATURES

- **Fast results:** Get accurate readings in just 5 seconds.
- **Large display:** The clear, easy-to-read screen is great for quick interpretation.
- **Easy operation:** A simple two-button interface makes testing intuitive for everyone.
- **Hygienic strip ejector:** Reduces contact with blood, helping to maintain hygiene.
- **Memory function:** Stores up to 500 test results to help you track trends.
- **Alternative site testing (AST):** Allows for blood sampling from various body parts, which can improve comfort.
- **Auto-coding:** Simply insert the strip and start testing – no manual coding required.
- **Compatible strips & lancing device:** Works with Gluco-Care test strips and a user-friendly lancing device.
- **Low blood sample required:** Needs just 0.5 µL of blood, minimizing discomfort.

TEST KIT: For Self Testing

- + 1 Blood Glucose Meter with batteries
- + 10 Blood Glucose Test Strips
- + 1 Adjustable Depth Lancing Device
- + 10 Sterilized Lancets
- + 1 Protective Carrying Case
- + User Manual, Package Inserts

Note: The contents of a Complete Kit may vary and are subject to change according to the actual purchase.



4006 ONE-CARE® Glucose Monitoring Kit for Self-Testing

ORDERING INFORMATION

REF	Description	PCS per BOX	BOXES per CASE
4001	ONE-CARE® PRO Glucose Monitor for Professional Use	1	94
4002	ONE-CARE® PRO Glucose Strips for Professional Use	50	560
4010	ONE-CARE® PRO Control Solution Level 1 & 3 (Low & High)	2	900
4006	ONE-CARE® Glucose Monitoring Kit for Self-Testing	1	80
4007	ONE-CARE® Glucose Strips for Self-Testing	50	560

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

Gluco-Care®

GLUCOSE MONITORING SYSTEM

DIAGNOSTICS

Blood Glucose Monitors

USERS:



PRODUCT FEATURES

- **Fast Results** – Delivers accurate readings in just 5 seconds.
- **Small Sample Size** – Requires only 0.5 µL of blood.
- **Auto Coding** – No need for manual calibration.
- **Large LCD Display** – Easy to read, even for users with poor eyesight.
- **Compact & Portable** – Lightweight design ideal for travel and everyday use.
- **500 Test Memory** – Stores past results for easy tracking and analysis.
- **Simple 2-Button Operation** – Hassle-free and intuitive to use.
- **Ejector Function** – Strip ejector promotes hygienic disposal.
- **Reliable Accuracy** – Meets ISO 15197:2013 standards.

Complete Test Kit

- + 1 Blood Glucose Meter
- + 10 Blood Glucose Test Strips
- + 1 Control Solution
- + 1 Adjustable Depth Lancing Device
- + 10 Sterilized Lancets
- + User Manual
- + Protective carrying case

Note: The contents of a Complete Kit may vary and are subject to change according to the actual purchase.



ORDERING INFORMATION

REF	Description	Packaging
4011	GLUCO-CARE® Blood Glucose Monitoring Kit (meter, 10 lancets, 1 lancing device, case, IFUs)	1 pcs/box 72 bx (72 pcs)/case
4012	GLUCO-CARE® Blood Glucose Test Strips	50 pcs/box 10 bx (500 pcs)/case
4013	GLUCO-CARE® Blood Glucose Monitoring Kit (meter, 10 lancets, 1 lancing device, case, 10 Glucose Test Strips, IFUs)	1 pcs/box 72 bx (72 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

One-Care®

NON-CONTACT INFRARED BODY THERMOMETERS

DIAGNOSTICS

Thermometers

USERS:



model:
HTD8813



model:
HT8818

PRODUCT FEATURES

- **Non-contact measurement:** Measures body temperature without touching the forehead or ear, which prevents cross-infection.
- **Accurate reading:** Provides an instant reading in just 1 second, with an accuracy within $\pm 0.2^{\circ}\text{C}$.
- **Intelligent memory:** Stores up to 50 measurements for model HTD8813 and up to 10 for model HT8818.
- **Silent mode:** Allows for temperature taking without disturbing a sleeping baby.
- **Adaptable units:** Can display temperature in both Fahrenheit ($^{\circ}\text{F}$) and Celsius ($^{\circ}\text{C}$).
- **Color-coded backlight:** The screen changes color to indicate temperature levels: green for normal, yellow for mild fever, and red for high fever (accompanied by a rapid alarm sound).
- **High-definition LCD display:** The screen is easy to read, even in the dark.
- **Convenient:** Suitable for patients of all ages, both at home and in a clinical environment.

COLOR BACKLIGHT

Representing results:

- Green** normal body temperature
- Yellow** mild fever
- Red** high fever (accompanied by rapid sound alarm)

ORDERING INFORMATION

REF	Description	Model	Packaging
4020	ONE-CARE® Non-Contact Infrared Thermometer	HTD8813	1 PC / Box 40 BX / Case
4021	ONE-CARE® Non-Contact Infrared Thermometer	HTD8818	1 PC / Box 40 BX / Case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

Drug Delivery



This category offers a wide array of solutions tailored for the safe and precise administration of injectable medications. Our products include standard and safety syringes, insulin and tuberculin syringes, as well as pen needles and safety pen needles designed to reduce the risk of accidental needlestick injuries. The range supports both professional and self-administration use, emphasizing safety, comfort, and accuracy.




One-Care®


LUER-LOCK SYRINGES
WITH SAFETY NEEDLES

DRUG DELIVERY
Syringes and
Safety needles


USERS:



THUMB ACTIVATION



FINGER ACTIVATION



SURFACE ACTIVATION



PRODUCT FEATURES

Silicone lubrication

- Ensures smooth plunger movement.

High-transparency barrel

- Makes it easy to see syringe content.

Safe plunger backstop

- Prevents accidental withdrawal of plunger.

Large and stable plunger grip

- For Stability and comfort.

Sterilized

- With Ethylene-Oxide.

One-Care® Safety Shield

- ✓ Easy one-hand activation.
- ✓ Audible click indicates safety.
- ✓ Locking mechanism.

ORDERING INFORMATION

REF	Description	Color Coding	Packaging
2110	ONE-CARE® Luer-Lock Syringes 1ml with Safety Needles 23G x 1" (25mm)		100 pcs/box 12 bx (1200 pcs)/case
2112	ONE-CARE® Luer-Lock Syringes 1ml with Safety Needles 25G x 1" (25mm)		100 pcs/box 12 bx (1200 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

One-Care[®]

SAFETY NEEDLES

DRUG DELIVERY Syringes and Safety needles

USERS:



23G x 1"(25mm)



25G x 1"(25mm)



21G x 1.5"(38mm)



18G x 1.5"(38mm)

PRODUCT FEATURES

- Reliable safety mechanism integrated with the needle.
- Convenient, one-hand use.
- Compatible with all Luer-Lock and Luer-Slip Syringes as per ISO 80369.
- Bevel oriented needle for easy injection.
- Compliant with performance requirements of ISO 7864.
- Color coding according to ISO 6009.



ORDERING INFORMATION

REF	Description	Color Coding	Packaging
2111	ONE-CARE [®] Safety Needles 23G x 1" (25mm)		100 pcs/box 20 bx (2,000 pcs)/case
2113	ONE-CARE [®] Safety Needles 25G x 1" (25mm)		100 pcs/box 20 bx (2,000 pcs)/case
2114	ONE-CARE [®] Safety Needles 21G x 1.5" (38mm)		100 pcs/box 20 bx (2,000 pcs)/case
2115	ONE-CARE [®] Safety Needles 18G x 1.5" (38mm)		100 pcs/box 20 bx (2,000 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

One-Care®

LUER-LOCK DISPOSABLE SYRINGES

DRUG DELIVERY Syringes and Safety needles

USERS:



PRODUCT FEATURES

Silicone lubrication

- Ensures smooth plunger movement.

High transparency barrel

- Easy to see syringe content.

Safe plunger backstop

- Prevents accidental withdrawal of plunger.

Large and stable plunger grip

- Ensures stability and comfort during injection.




Low Dead Space
To minimize medication loss

ORDERING INFORMATION

REF	Description	Size	Packaging
2002	ONE-CARE® Luer-Lock Syringes, 3-Part	1ml	100 pcs/box 32 bx (3,200 pcs)/case
2004	ONE-CARE® Luer-Lock Syringes, 3-Part	3ml	100 pcs/box 24 bx (2,400 pcs)/case
2005	ONE-CARE® Luer-Lock Syringes, 3-Part	5ml	100 pcs/box 16 bx (1,600 pcs)/case
2006	ONE-CARE® Luer-Lock Syringes, 3-Part	10ml	100 pcs/box 12 bx (1,200 pcs)/case
2007	ONE-CARE® Luer-Lock Syringes, 3-Part	20ml	100 pcs/box 9 bx (900 pcs)/case
2009	ONE-CARE® Luer-Lock Syringes, 3-Part	60ml	60 pcs/box 4 bx (240 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

USERS:



PRODUCT FEATURES

SAFE

- Reliable shielding mechanism for both transportation position and permanent lock.
- ETO-Sterilized and protected with color-coded cover and blister pack.
- Covering of needle can be done with only one hand by pushing the shield forward towards needle tip to reduce risk of sharps injuries.

SIMPLE

- Easy to use, intuitive operation as per industry standard.
- Safety mechanism engagement indicator via “Click” sound once the shield is pushed forward and turned either direction.
- Easy to read marking for precise dosing.

GENTLE

- Precisely ground, electro-polished, silicone-coated, ultra-thin needles on the end of patients comfort.
- Smooth and easy movement of plunger for gentle injection.

CONVENIENT

- Minimal dead space to reduce medication loss.
- Smart design with temporary needle shielding for transportation and permanent lock after drug delivery to prevent sharps injury and reuse.
- Easy to recognize color-coded product versions address a variety of dose requirements.

ORDERING INFORMATION

REF	Description	Size	Color Coding	Packaging
2101	ONE-CARE® Insulin Safety Syringes 29G x 1/2" (0.3mm x 12.7mm)	1ml		100 pcs/box 10 bx (1,000 pcs)/case
2102	ONE-CARE® Insulin Safety Syringes 29G x 1/2" (0.3mm x 12.7mm)	0.5ml		100 pcs/box 10 bx (1,000 pcs)/case

REF	Description	Size	Color Coding	Packaging
2107	ONE-CARE® Tuberculin Safety Syringes 28G x 1/2" (0.36mm x 12.7mm)	1ml		100 pcs/box 10 bx (1,000 pcs)/case
2108	ONE-CARE® Tuberculin Safety Syringes 27G x 1/2" (0.40mm x 12.7mm)	1ml		100 pcs/box 10 bx (1,000 pcs)/case
2109	ONE-CARE® Tuberculin Safety Syringes 25G x 5/8" (0.50mm x 16.0mm)	1ml		100 pcs/box 10 bx (1,000 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

One-Care®

INSULIN SAFETY SYRINGES OPTI

DRUG DELIVERY Syringes and Safety needles

USERS:



PRODUCT FEATURES

Safe:

- Durable shielding system for a permanent lock for secure disposal after use.
- Quick thumb-activated shield covers the needle instantly, minimizing sharp injury risks.
- ETO-sterilized and sealed in a blister pack with a color-coded protective cap.

Simple:

- User-friendly, intuitive handling consistent with Long Term Care facility standards.
- Audible “click” confirms safety mechanism engagement when the shield is pushed forward and turned in either direction.
- Clear, easy-to-read scale markings for accurate dosing.

Gentle:


- Ultra-thin, precisely ground, electro-polished, silicone-coated needles designed for maximum patient comfort.
- Multiple color-coded versions to suit different skin types.
- Smooth, controlled plunger action for a comfortable injection experience.

Convenient:

- Reduced dead space to limit medication waste.
- Innovative design offers temporary needle shielding during transport and a permanent lock after injection to prevent sharps injuries and reuse.
- Color-coded product versions for quick identification and correct selection.



One-step
and one-hand
activation



ORDERING INFORMATION

REF	Description	Size	Color Coding	Packaging
2141	ONE-CARE® OPTI Insulin Safety Syringes 29G x 1/2" (0.33mm x 12.7mm)	1ml		100 pcs/box 8 bx (800 pcs)/case
2142	ONE-CARE® OPTI Insulin Safety Syringes 29G x 1/2" (0.33mm x 12.7mm)	0.5ml		100 pcs/box 8 bx (800 pcs)/case
2145	ONE-CARE® OPTI Insulin Safety Syringes 29G x 1/2" (0.33mm x 12.7mm)	0.3ml		100 pcs/box 8 bx (800 pcs)/case
2146	ONE-CARE® OPTI Insulin Safety Syringes 30G x 5/16" (0.33mm x 8.0mm)	1ml		100 pcs/box 8 bx (800 pcs)/case
2147	ONE-CARE® OPTI Insulin Safety Syringes 30G x 5/16" (0.33mm x 8.0mm)	0.5ml		100 pcs/box 8 bx (800 pcs)/case

REF	Description	Size	Color Coding	Packaging
2143	ONE-CARE® OPTI Tuberculin Safety Syringes 28G x 1/2" (0.36mm x 12.7mm)	1ml		100 pcs/box 8 bx (800 pcs)/case
2144	ONE-CARE® OPTI Tuberculin Safety Syringes 27G x 1/2" (0.36mm x 12.7mm)	1ml		100 pcs/box 8 bx (800 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

One-Care[®]

INSULIN SYRINGES

DRUG DELIVERY

Syringes and Safety needles

USERS:



PRODUCT FEATURES

MediVena ONE-CARE[®] Insulin Syringes are available in multiple sizes and configurations to match leading brands, supporting comfortable, accurate, and compliant insulin delivery.

CONVENIENT

- Wide flaps for a secure, stable grip.
- Smooth plunger movement for easy injection.

GENTLE


- Tri-beveled, electro-polished needle with medical-grade lubricant.
- Thin-wall design with small-diameter options for delicate injections.

PRECISE

- Clear graduation markings for accurate dosing.
- Reduced dead space to minimize medication loss.

VERSATILE

- Available in 0.3 mL, 0.5 mL, and 1 mL sizes.
- Multiple needle options: 29G to 31G.



Low dead space

Precise needle orientation

Flat base

ORDERING INFORMATION

REF	Description	Size	Packaging
2271	ONE-CARE [®] Insulin Syringes 31G x 5/16" (8mm)	0.3ml/cc	10 pcs/bag (100pcs/ box) 40 bx (4,000 pcs)/case
2272	ONE-CARE [®] Insulin Syringes 30G x 5/16" (8mm)	0.3ml/cc	10 pcs/bag (100pcs/ box) 40 bx (4,000 pcs)/case
2273	ONE-CARE [®] Insulin Syringes 30G x 5/16" (8mm)	0.5ml/cc	10 pcs/bag (100pcs/ box) 40 bx (4,000 pcs)/case
2274	ONE-CARE [®] Insulin Syringes 31G x 5/16" (8mm)	1.0ml/cc	10 pcs/bag (100pcs/ box) 40 bx (4,000 pcs)/case
2275	ONE-CARE [®] Insulin Syringes 30G x 5/16" (8mm)	1.0ml/cc	10 pcs/bag (100pcs/ box) 40 bx (4,000 pcs)/case
2276	ONE-CARE [®] Insulin Syringes 30G x 1/2" (12.7mm)	1.0ml/cc	10 pcs/bag (100pcs/ box) 40 bx (4,000 pcs)/case
2277	ONE-CARE [®] Insulin Syringes 29G x 1/2" (12.7mm)	1.0ml/cc	10 pcs/bag (100pcs/ box) 40 bx (4,000 pcs)/case
2278	ONE-CARE [®] Insulin Syringes 27G x 1/2" (12.7mm)	1.0ml/cc	10 pcs/bag (100pcs/ box) 40 bx (4,000 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

One-Care[®]

INSULIN SYRINGES

DRUG DELIVERY

Syringes and Safety needles

USERS:



PRODUCT FEATURES

MediVena ONE-CARE[®] Insulin Syringes come in various sizes and configurations to align with top brands, ensuring comfortable, precise, and compliant insulin delivery. They are individually packed in blisters.

CONVENIENT

- Wide flaps for a secure, stable grip.
- Smooth plunger movement for easy injection.
- Individually packed in sterile blisters.

GENTLE

- Tri-beveled, electro-polished needle with medical-grade lubricant.
- Thin-wall design with small-diameter options for delicate injections.

PRECISE

- Clear graduation markings for accurate dosing.
- Reduced dead space to minimize medication loss.

VERSATILE

- Available in 0.3 mL, 0.5 mL, and 1 mL sizes.
- Multiple needle options: 29G to 31G.

INDIVIDUALLY PACKED



1pcs Sterile Blisters



Low dead space



Precise needle orientation



Flat base

ORDERING INFORMATION

REF	Description	Size	Packaging
2288	ONE-CARE [®] Insulin Syringes 31G x 6mm	0.3ml/cc	10 pcs/bag (100pcs/ box) 30 bx (3,000 pcs)/case
2289	ONE-CARE [®] Insulin Syringes 31G x 6mm	0.5ml/cc	10 pcs/bag (100pcs/ box) 30 bx (3,000 pcs)/case
2290	ONE-CARE [®] Insulin Syringes 31G x 6mm	1.0ml/cc	10 pcs/bag (100pcs/ box) 30 bx (3,000 pcs)/case
2291	ONE-CARE [®] Insulin Syringes 31G x 8mm	0.3ml/cc	10 pcs/bag (100pcs/ box) 30 bx (3,000 pcs)/case
2292	ONE-CARE [®] Insulin Syringes 31G x 8mm	0.5ml/cc	10 pcs/bag (100pcs/ box) 30 bx (3,000 pcs)/case
2293	ONE-CARE [®] Insulin Syringes 31G x 8mm	1.0ml/cc	10 pcs/bag (100pcs/ box) 30 bx (3,000 pcs)/case
2294	ONE-CARE [®] Insulin Syringes 30G x 12.7mm	0.3ml/cc	10 pcs/bag (100pcs/ box) 30 bx (3,000 pcs)/case
2295	ONE-CARE [®] Insulin Syringes 30G x 12.7mm	0.5ml/cc	10 pcs/bag (100pcs/ box) 30 bx (3,000 pcs)/case
2296	ONE-CARE [®] Insulin Syringes 30G x 12.7mm	1.0ml/cc	10 pcs/bag (100pcs/ box) 30 bx (3,000 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

Gluco-Care®

INSULIN SYRINGES

DRUG DELIVERY

Syringes and Safety needles

USERS:



PRODUCT FEATURES



100 Pcs/box

Engineered for precision, comfort, and ease of use, GLUCO-CARE® Insulin Syringes deliver reliable performance that supports patient compliance and confidence in daily diabetes care.

Ergonomic Control

- Wide, textured finger flanges for a steady grip.
- Smooth-glide plunger enables controlled, low-resistance injections.

Comfort-Focused Design

- Ultra-thin, tri-bevel needles are electro-polished and silicone-coated for a gentler injection.
- Available in small gauge sizes and thin-wall designs to reduce pressure at the injection site.

Dose Accuracy

- Clear, bold graduation marks for accurate dose visualization.
- Low dead space design helps reduce insulin waste.

Versatile Options

- Offered in 0.3 mL, 0.5 mL, and 1 mL volumes.
- Available needle sizes range from 29G to 31G to suit various patient needs.

ORDERING INFORMATION

REF	Description	Size	Packaging
2241	GLUCO-CARE® Insulin Syringes 31G x 5/16" (8mm), U-100	0.3ml/cc	10 pcs/bag (100pcs/ box) 40 bx (4,000 pcs)/case
2242	GLUCO-CARE® Insulin Syringes 30G x 5/16" (8mm), U-100	0.3ml/cc	10 pcs/bag (100pcs/ box) 40 bx (4,000 pcs)/case
2247	GLUCO-CARE® Insulin Syringes 31G x 5/16" (8mm), U-100	0.5ml/cc	10 pcs/bag (100pcs/ box) 40 bx (4,000 pcs)/case
2248	GLUCO-CARE® Insulin Syringes 30G x 5/16" (8mm), U-100	0.5ml/cc	10 pcs/bag (100pcs/ box) 40 bx (4,000 pcs)/case
2250	GLUCO-CARE® Insulin Syringes 29G x 1/2" (12.7mm), U-100	0.5ml/cc	10 pcs/bag (100pcs/ box) 40 bx (4,000 pcs)/case
2254	GLUCO-CARE® Insulin Syringes 31G x 5/16" (8mm), U-100	1.0ml/cc	10 pcs/bag (100pcs/ box) 40 bx (4,000 pcs)/case
2255	GLUCO-CARE® Insulin Syringes 30G x 5/16" (8mm), U-100	1.0ml/cc	10 pcs/bag (100pcs/ box) 40 bx (4,000 pcs)/case
2256	GLUCO-CARE® Insulin Syringes 30G x 1/2" (12.7mm), U-100	1.0ml/cc	10 pcs/bag (100pcs/ box) 40 bx (4,000 pcs)/case
2257	GLUCO-CARE® Insulin Syringes 29G x 1/2" (12.7mm), U-100	1.0ml/cc	10 pcs/bag (100pcs/ box) 40 bx (4,000 pcs)/case
2259	GLUCO-CARE® Insulin Syringes 27G x 1/2" (12.7mm), U-100	1.0ml/cc	10 pcs/bag (100pcs/ box) 40 bx (4,000 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

Gluco-Care®

INSULIN SYRINGES

DRUG DELIVERY

Syringes and Safety needles

USERS:



PRODUCT FEATURES



50 Pcs/box

Engineered for precision, comfort, and ease of use, GLUCO-CARE® Insulin Syringes deliver reliable performance that supports patient compliance and confidence in daily diabetes care.

Ergonomic Control

- Wide, textured finger flanges for a steady grip.
- Smooth-glide plunger enables controlled, low-resistance injections.

Comfort-Focused Design

- Ultra-thin, tri-bevel needles are electro-polished and silicone-coated for a gentler injection.
- Available in small gauge sizes and thin-wall designs to reduce pressure at the injection site.

Dose Accuracy

- Clear, bold graduation marks for accurate dose visualization.
- Low dead space design helps reduce insulin waste.

Versatile Options

- Offered in 0.3 mL, 0.5 mL, and 1 mL volumes.
- Available needle sizes range from 29G to 31G to suit various patient needs.

ORDERING INFORMATION

REF	Description	Size	Packaging
2241-50	GLUCO-CARE® Insulin Syringes 31G x 5/16" (8mm), U-100	0.3ml/cc	10 pcs/bag (50pcs/ box) 60 bx (3,000 pcs)/case
2242-50	GLUCO-CARE® Insulin Syringes 30G x 5/16" (8mm), U-100	0.3ml/cc	10 pcs/bag (50pcs/ box) 60 bx (3,000 pcs)/case
2247-50	GLUCO-CARE® Insulin Syringes 31G x 5/16" (8mm), U-100	0.5ml/cc	10 pcs/bag (50pcs/ box) 60 bx (3,000 pcs)/case
2248-50	GLUCO-CARE® Insulin Syringes 30G x 5/16" (8mm), U-100	0.5ml/cc	10 pcs/bag (50pcs/ box) 60 bx (3,000 pcs)/case
2250-50	GLUCO-CARE® Insulin Syringes 29G x 1/2" (12.7mm), U-100	0.5ml/cc	10 pcs/bag (50pcs/ box) 60 bx (3,000 pcs)/case
2254-50	GLUCO-CARE® Insulin Syringes 31G x 5/16" (8mm), U-100	1.0ml/cc	10 pcs/bag (50pcs/ box) 60 bx (3,000 pcs)/case
2255-50	GLUCO-CARE® Insulin Syringes 30G x 5/16" (8mm), U-100	1.0ml/cc	10 pcs/bag (50pcs/ box) 60 bx (3,000 pcs)/case
2256-50	GLUCO-CARE® Insulin Syringes 30G x 1/2" (12.7mm), U-100	1.0ml/cc	10 pcs/bag (50pcs/ box) 60 bx (3,000 pcs)/case
2257-50	GLUCO-CARE® Insulin Syringes 29G x 1/2" (12.7mm), U-100	1.0ml/cc	10 pcs/bag (50pcs/ box) 60 bx (3,000 pcs)/case
2259-50	GLUCO-CARE® Insulin Syringes 27G x 1/2" (12.7mm), U-100	1.0ml/cc	10 pcs/bag (50pcs/ box) 60 bx (3,000 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.



USERS:



PRODUCT FEATURES

**U-40
&
U-100
AVAILABLE**



Accurate Care

- Half-Unit Markings for accurate dosing.
- Clear, bold graduation marks for precise dose visualization.
- Low dead space design helps reduce insulin waste.


Versatile Options

- U-40 and U-100 options available.
- Offered in 0.3 mL, 0.5 mL, and 1 mL volumes.
- Available needle sizes range from 29G to 31G to suit various patient needs.

Comfortable Design

- Ultra-thin, tri-bevel needles are electro-polished and silicone-coated for a gentler injection.
- Available in small gauge sizes and thin-wall designs to reduce pressure at the injection site.

Half-Unit Markings
for accurate dosing



100 PCS boxes

50 PCS boxes

ORDERING INFORMATION

REF	Description	Size	Packaging
2501 2501-50	Insulin Syringes 29G x 1/2" (12.7mm), U-40, 12 Units Scale	0.3ml/cc	10pcs/bag, 100pcs/box, 40bx/case (4,000pcs) 10pcs/bag, 50pcs/box, 80bx/case (4,000pcs)
2502 2502-50	Insulin Syringes 29G x 1/2" (12.7mm), U-40, 20 Units Scale	0.5ml/cc	10pcs/bag, 100pcs/box, 40bx/case (4,000pcs) 10pcs/bag, 50pcs/box, 80bx/case (4,000pcs)
2503 2503-50	Insulin Syringes 29G x 1/2" (12.7mm), U-40, 40 Units Scale	1.0ml/cc	10pcs/bag, 100pcs/box, 40bx/case (4,000pcs) 10pcs/bag, 50pcs/box, 80bx/case (4,000pcs)
2505 2505-50	Insulin Syringes 31G x 5/16" (8mm), U-100, 30 Units Scale	0.3ml/cc	10pcs/bag, 100pcs/box, 40bx/case (4,000pcs) 10pcs/bag, 50pcs/box, 80bx/case (4,000pcs)
2506 2506-50	Insulin Syringes 29G x 1/2" (12.7mm), U-100, 50 Units Scale	0.3ml/cc	10pcs/bag, 100pcs/box, 40bx/case (4,000pcs) 10pcs/bag, 50pcs/box, 80bx/case (4,000pcs)
2507 2507-50	Insulin Syringes 29G x 1/2" (12.7mm), U-100, 100 Units Scale	0.5ml/cc	10pcs/bag, 100pcs/box, 40bx/case (4,000pcs) 10pcs/bag, 50pcs/box, 80bx/case (4,000pcs)

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.




INSULIN SYRINGES


USERS:



**U-40
&
U-100
AVAILABLE**



Half-Unit Markings
for accurate dosing



PRODUCT FEATURES

Accurate Care

- Half-Unit Markings for accurate dosing.
- Clear, bold graduation marks for precise dose visualization.
- Low dead space design helps reduce insulin waste.

Versatile Options

- U-40 and U-100 options available.
- Offered in 0.3 mL, 0.5 mL, and 1 mL volumes.
- Available needle sizes range from 29G to 31G to suit various patient needs.

Comfortable Design

- Ultra-thin, tri-bevel needles are electro-polished and silicone-coated for a gentler injection.
- Available in small gauge sizes and thin-wall designs to reduce pressure at the injection site.



Blister Packed



ORDERING INFORMATION

REF	Description	Size	Packaging
2511	Insulin Syringes 29G x 1/2" (12.7mm), U-40, 12 Units Scale	0.3ml/cc	1pcs/blister (100pcs/ box) 30bx (3,000 pcs)/case
2512	Insulin Syringes 29G x 1/2" (12.7mm), U-40, 20 Units Scale	0.5ml/cc	1pcs/blister (100pcs/ box) 30bx (3,000 pcs)/case
2513	Insulin Syringes 29G x 1/2" (12.7mm), U-40, 40 Units Scale	1.0ml/cc	1pcs/blister (100pcs/ box) 30bx (3,000 pcs)/case
2515	Insulin Syringes 31G x 5/16" (8mm), U-100, 30 Units Scale	0.3ml/cc	1pcs/blister (100pcs/ box) 30bx (3,000 pcs)/case
2516	Insulin Syringes 29G x 1/2" (12.7mm), U-100, 30 Units Scale	0.3ml/cc	1pcs/blister (100pcs/ box) 30bx (3,000 pcs)/case
2517	Insulin Syringes 29G x 1/2" (12.7mm), U-100, 50 Units Scale	0.5ml/cc	1pcs/blister (100pcs/ box) 30bx (3,000 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

One-Care[®]

SAFETY PEN NEEDLES

DRUG DELIVERY

Safety Pen Needles

USERS:



PRODUCT FEATURES

Safe

- Automated shielding mechanism ensures the needle is covered before and after use.
- Concealed needle design helps to relieve the patient's anxiety.
- Single-use design prevents accidental needle stick injuries and reuse of the device for the safety of patients and healthcare professionals.
- All pen needles are ETO-sterilized and protected with an outer cover and sterility peel-back tab.

Simple

- Easy to use, intuitive handling with no pre-loading required.
- Clear shield with optimized skin contact area allows for improved visibility of insulin flow during priming and precise positioning on the injection site.
- Safety mechanism engagement can be easily recognized by the "Click" sound when the needle is fully inserted into the skin.
- Color-coded inner shield automatically retracts, allowing the user to determine that the needle was used and is locked permanently.

Ergonomic

- Electro-polished, silicone-coated, ultra-thin needles guarantee enhanced user comfort.
- Easy flow cannula minimizes the time of insulin delivery.
- Smooth and quick insert/retract mechanism makes the injection a barely noticeable experience.
- Ergonomic design with ribbed housing ensures a comfortable grip, secure attachment to the pen injector, as well as easy removal.

Universal

- Universal design for compatibility with most injection pens available on the market.
- Tight fit prevents leakages and medication loss.
- Rigorously tested per international standards for fit, functionality, and reliable performance.



BEFORE USE



AFTER USE

ORDERING INFORMATION

REF	Description	Color Coding	Packaging
3101	ONE-CARE [®] Safety Pen Needles, Micro, 30G x 5mm (3/16")		100 pcs/box 24 bx (2,400 pcs)/case
3102	ONE-CARE [®] Safety Pen Needles, Normal, 30G x 8mm (5/16")		100 pcs/box 24 bx (2,400 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

Gluco-Care®

SAFETY PEN NEEDLES

DRUG DELIVERY

Safety Pen Needles

USERS:



PRODUCT FEATURES



Safe

- Automated shielding mechanism covers the needle before and after use.
- Concealed needle design helps relieve patient anxiety.
- Single-use design prevents sharps injuries and reuse of the device.
- Gamma-sterilized needles, protected with a cover and peel-off tab.

Simple

- Intuitive handling—no pre-loading required.
- Clear shield improves visibility of insulin flow and helps precise injection placement.
- Audible “click” confirms needle is fully inserted and safety mechanism engaged.
- Color-coded inner shield retracts automatically, showing the needle has been used and is permanently locked.

Comfortable

- Ergonomic ribbed housing for a comfortable grip.
- Ultra-thin, polished, silicone-coated needles for smooth injections.
- Easy-flow cannula minimizes insulin delivery time.
- Compatible with most injection pens available on the market.



BEFORE USE



DURING USE



AFTER USE

ORDERING INFORMATION

REF	Description	Color Coding	Packaging
3111	GLUCO-CARE® Safety Pen Needles, 31G x 5mm (3/16")		100 pcs/box 40 bx (4000 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

USERS:



PRODUCT FEATURES

Safe:

- Automated shielding mechanism ensures the needle is covered both before and after use.
- Single-use design prevents sharps injuries and eliminates the risk of device reuse.
- All products are ETO-sterilized and protected with a cover and sterile safety tab.

Simple:

- Safety mechanism engagement is indicated by a “Click” sound and the disappearance of the color-coded inner shield once the needle is fully inserted into the skin.
- Clear, optimally designed shield allows precise placement on the injection site, enhancing visibility.

Gentle:

- Concealed needle design helps reduce patient anxiety.
- Polished, silicone-coated ultra-thin needles are carefully positioned for smooth skin penetration, supported by precise tracks in the safety mechanism for greater comfort.
- Ribbed housing offers a secure and comfortable grip, simplifying attachment and removal from the pen.

Universal:

- Compatible with most injection pens available on the market.
- Tight fit prevents insulin leakage.
- Color-coded versions make it easy to identify the right product for different patient skin types and injection needs.



BEFORE USE

AFTER USE

ORDERING INFORMATION

REF	Description	Color Coding	Packaging
3301	ONE-CARE® OPTI Safety Pen Needles, 30G x 5mm (3/16")		100 pcs/box 20 bx (2,000 pcs)/case
3302	ONE-CARE® OPTI Safety Pen Needles, 30G x 8mm (5/16")		100 pcs/box 20 bx (2,000 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.


One-Care[®]

SAFETY PEN NEEDLES DUO

DRUG DELIVERY

Safety Pen Needles

USERS:



Universal

- Compatible with the most Pen Injectors available on the market.
- Available in multiple needle configurations.

PRODUCT FEATURES

Safe:

- Needle front is pre-covered before and after injection to reduce the risk of needle-stick injuries.
- Back-end needle automatically covered after removal from the pen.
- Red visual indicator for easy distinction between used and unused devices.

Convenient:

- Thin-wall technology for optimal insulin flow.
- Reduced resistance during safety mechanism activation and needle penetration.
- Ribbed main housing for a comfortable grip and secure attachment to the pen injector.



BEFORE USE



AFTER USE

ORDERING INFORMATION

REF	Description	Color Coding	Packaging
3350	ONE-CARE [®] DUO Safety Pen Needles, 30G x 5mm (3/16")		100 pcs/box 8 bx (800 pcs)/case
3351	ONE-CARE [®] DUO Safety Pen Needles, 30G x 8mm (5/16")		100 pcs/box 8 bx (800 pcs)/case
3360	ONE-CARE [®] DUO Safety Pen Needles, 31G x 5mm (3/16")		100 pcs/box 16 bx (1,600 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

One-Care[®]

PEN NEEDLES

DRUG DELIVERY

Pen Needles

USERS:



PRODUCT FEATURES

Gentle:

- Electro-polished, silicone-coated, ultra-thin needles maximize patient comfort.
- Thin-wall and extra-thin wall technology allow smooth insulin flow, reducing injection time.
- Special lubrication ensures effortless penetration and minimal sensation during injections.

Comfortable:

- Ergonomic, ribbed housing provides a secure grip and firm connection to the pen injector.
- Color-coded indicators make it easy to choose the correct needle size.
- ETO-sterilized and protected with an outer cover and easy-peel seal for safety.

Universal

- Fits most pen injectors available on the market.
- Leak-proof design prevents medication loss; 100% automated inspection guarantees quality.
- Rigorously tested to meet international standards for safety, functionality, and reliability.



ORDERING INFORMATION

REF	Description	Color Coding	Packaging
3011	ONE-CARE [®] Pen Needles, 33G x 4mm		100 pcs/box 100 bx (10,000 pcs)/case
3012	ONE-CARE [®] Pen Needles, 32G x 4mm		100 pcs/box 100 bx (10,000 pcs)/case
3013	ONE-CARE [®] Pen Needles, 32G x 5mm		100 pcs/box 100 bx (10,000 pcs)/case
3014	ONE-CARE [®] Pen Needles, 32G x 6mm		100 pcs/box 100 bx (10,000 pcs)/case
3015	ONE-CARE [®] Pen Needles, 32G x 8mm		100 pcs/box 100 bx (10,000 pcs)/case
3016	ONE-CARE [®] Pen Needles, 31G x 4mm		100 pcs/box 100 bx (10,000 pcs)/case
3017	ONE-CARE [®] Pen Needles, 31G x 5mm		100 pcs/box 100 bx (10,000 pcs)/case
3018	ONE-CARE [®] Pen Needles, 31G x 6mm		100 pcs/box 100 bx (10,000 pcs)/case
3019	ONE-CARE [®] Pen Needles, 31G x 8mm		100 pcs/box 100 bx (10,000 pcs)/case
3020	ONE-CARE [®] Pen Needles, 30G x 5mm		100 pcs/box 100 bx (10,000 pcs)/case
3021	ONE-CARE [®] Pen Needles, 29G x 12.7mm		100 pcs/box 100 bx (10,000 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

Gluco-Care[®]

PEN NEEDLES PLUS

DRUG DELIVERY

Pen Needles

USERS:



Universal

- Universally compatible with nearly all pen injectors on the market.
- Precision fit prevents leakage and ensures proper medication dosing.
- Comprehensively tested to meet global standards for safety, performance, and reliability.

PRODUCT FEATURES


Smooth

- Ultra-thin, silicone-coated, electro-polished needles deliver maximum injection comfort.
- Advanced thin-wall technology supports quick insulin delivery and less discomfort.
- Enhanced lubrication makes each insertion smooth, consistent, and barely noticeable.

Gentle:

- Flat needle base enlarges the skin contact surface, providing a more comfortable injection experience.
- Ribbed, ergonomic design ensures stable handling and secure pen connection.
- Clear color coding helps patients identify the right size at a glance.
- ETO-sterilized with protective outer cover and easy-peel tab for hygiene and convenience.

PEN NEEDLES PLUS OTHER PRODUCTS



ORDERING INFORMATION

REF	Description	Color Coding	Packaging
3022	GLUCO-CARE [®] PLUS Pen Needles, 32G x 4mm (5/32")		100 pcs/box 100 bx (10,000 pcs)/case
3023	GLUCO-CARE [®] PLUS Pen Needles, 32G x 5mm (3/16")		100 pcs/box 100 bx (10,000 pcs)/case
3024	GLUCO-CARE [®] PLUS Pen Needles, 32G x 6mm (1/4")		100 pcs/box 100 bx (10,000 pcs)/case
3025	GLUCO-CARE [®] PLUS Pen Needles, 33G x 4mm (5/32")		100 pcs/box 100 bx (10,000 pcs)/case
3026	GLUCO-CARE [®] PLUS Pen Needles, 31G x 5mm (3/16")		100 pcs/box 100 bx (10,000 pcs)/case
3027	GLUCO-CARE [®] PLUS Pen Needles, 31G x 6mm (1/4")		100 pcs/box 100 bx (10,000 pcs)/case
3028	GLUCO-CARE [®] PLUS Pen Needles, 31G x 8mm (5/16")		100 pcs/box 100 bx (10,000 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MedVena directly.

Feed and Saline Flush Syringes



Our feeding and flushing syringes provide reliable and hygienic solutions for both feeding and flushing applications. The portfolio includes oral syringes for accurate medication dosing, ENFit-compatible enteral syringes for safe tube feeding, and pre-filled saline flush syringes to maintain IV line patency and reduce the risk of infection. Each product is developed with patient safety and clinical efficiency in mind.



One-Care[®]

ORAL SYRINGES

FEED AND SALINE
FLUSH SYRINGES

Oral Syringes

USERS:



PRODUCT FEATURES

- Clear, easy-to-read scale ensures the accuracy of dosing.
- Wide surface of the plunger for comfortable grip and ease of use.
- Special slip tip secured with a cap.
- Sterilized with Ethylene Oxide, non-toxic, non-pyrogenic.
- Made of latex-free, medical-grade plastic.

ORDERING INFORMATION

REF	Description	Size	Color Coding	Packaging
2011	ONE-CARE [®] Oral Syringes, Mounted Cap	1ml		100 pcs/box 18 bx (1,800 pcs)/case
2012	ONE-CARE [®] Oral Syringes, Mounted Cap	3ml		100 pcs/box 16 bx (1,600 pcs)/case
2013	ONE-CARE [®] Oral Syringes, Mounted Cap	5ml		100 pcs/box 12 bx (1,200 pcs)/case
2014	ONE-CARE [®] Oral Syringes, Mounted Cap	10ml		100 pcs/box 9 bx (900 pcs)/case
2015	ONE-CARE [®] Oral Syringes, Mounted Cap	20ml		100 pcs/box 6 bx (600 pcs)/case
2017	ONE-CARE [®] Oral Syringes, Mounted Cap	60ml		25 pcs/box 8 bx (200 pcs)/case
2018	ONE-CARE [®] Oral Syringes, Transparent Cap, Black Scale	10ml		100 pcs/box 12 bx (1200 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

One-Care®

ENTEREAL SYRINGES

FEED AND SALINE FLUSH SYRINGES

ENFit Syringes

USERS:



PRODUCT FEATURES

- Mounted Cap.
- ISO80369-3 Entereal Connector Standard.
- Easy to read Scale.
- Comfortable and stable grip.
- ENFit tip ISO80369-3 with low dose version.

The single-use, Entereal Syringes are FDA cleared for direct oral administration and enteral delivery.

Enfit™ system is the new connection standard, developed by an international group of clinicians, manufacturers and regulators, as part of the Stay Connected initiative to prevent tubing misconnections and to improve patients safety.



ORDERING INFORMATION

REF	Description	Size	Color Coding	Packaging
2021	ONE-CARE® ENFit Syringes, with Cap, Blister-packed, Sterile, Low Dose	1ml		30 pcs/box 64 bx (1920 pcs)/case
2022	ONE-CARE® ENFit Syringes, with Cap, Blister-packed, Sterile, Low Dose	2.5ml		30 pcs/box 96 bx (2880 pcs)/case
2023	ONE-CARE® ENFit Syringes, with Cap, Blister-packed, Sterile	5ml		30 pcs/box 64 bx (1920 pcs)/case
2024	ONE-CARE® ENFit Syringes, with Cap, Blister-packed, Sterile	10ml		30 pcs/box 48 bx (1440 pcs)/case
2025	ONE-CARE® ENFit Syringes, with Cap, Blister-packed, Sterile	20ml		30 pcs/box 24 bx (720 pcs)/case
2027	ONE-CARE® ENFit Syringes, with Cap, Blister-packed, Sterile	60ml		30 pcs/box 10 bx (300 pcs)/case
2028	ONE-CARE® ENFit Syringes, with Cap, Blister-packed, Sterile	100ml		10 pcs/box 14 bx (140 pcs)/case

REF	Description	Size	Product Photo	Packaging
2029	ONE-CARE® ENFit Transition Connector	N/A		30 pcs/box 72 bx (720 pcs)/case
2040	ONE-CARE® ENFit Syringes with Transition Connector	60ml		25 pcs/box 8 bx (200 pcs)/case
2041 2041-NS	ONE-CARE® Thumb Control Ring Syringe, ENFit Tip, Protector Cap	60ml		Sterile 50 pcs/case Non-Sterile 30 pcs/box, 6bx/case
2042 2042-NS	ONE-CARE® Thumb Control Ring Syringe, Syringe Tip Adapter, Protector Cap	60ml		Sterile 50 pcs/ case Non-Sterile 30 pcs/box, 6bx/ case
2044 2044-NS	ONE-CARE® Flat Top Piston Syringes, with cap & tip adapter	60ml		Sterile 25pcs/box, 10bx/case Non-Sterile 30 pcs/box, 12bx/ case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

One-Care[®]

SALINE FLUSH SYRINGES

FEED AND SALINE
FLUSH SYRINGES

Saline Flush Syringes

USERS:



PRODUCT FEATURES

- **Compatible with needleless systems available on the market**
- **Prefilled with 0.9% sodium chloride injection**
(Conforming with USP standards)
- **Sterile Field or Sterile Fluid Path***
*see ORDERING INFORMATION

Flushing is vital for catheter care - ensuring functionality, maintaining patency, and preventing drug interactions or blockages. Saline flushing helps reduce the risk of thrombus formation and keeps the line clear.

ORDERING INFORMATION


REF	Description	Size	Color Coding	Packaging
2301	ONE-CARE [®] Saline Flush Syringes, Sterile Field	3ml		30 pcs/box 20 bx (600 pcs)/case
2302	ONE-CARE [®] Saline Flush Syringes, Sterile Field	5ml		30 pcs/box 20 bx (600 pcs)/case
2304	ONE-CARE [®] Saline Flush Syringes, Sterile Fluid Path	10ml		30 pcs/box 16 bx (480 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

Skin Preparation SOLUTIONS



Our personal protective equipment provides essential hygiene and safety protection for healthcare professionals and patients. Designed to minimize contamination and ensure a clean, safe environment.



One-Care®

ALCOHOL PREP PADS

PERSONAL PROTECTIVE EQUIPMENT

Alcohol Prep Pads

USERS:



PRODUCT FEATURES



- Soft, absorbent, nonwoven pads saturated with 70% v/v isopropyl alcohol.
- Exceeds CDC's recommended 60% alcohol content for optimal antiseptic action, helping to reduce bacteria.
- Suitable for individual, professional, and hospital use.
- Reliable 4-layer packaging.
- Compact box available in 100 or 200 pads.

ONE-CARE® Alcohol Prep Pads are saturated with 70% Isopropyl Alcohol for optimal antiseptic action (anti-bacteria, anti-virus).

They are intended for topical use in preparation of skin prior to an injection or blood collection, cleaning of minor cuts and scrapes, or cleaning of medical equipment and supplies. Reliable 4-layer packaging provides air-tight seal and prevents leakage or dry-out.



ORDERING INFORMATION

REF	Description	Size	Packaging
1022	ONE-CARE® Alcohol Prep Pads, 70% Isopropyl Alcohol, 200/box	M	200 pcs/box 50 bx (10,000 pcs)/case
1023	ONE-CARE® Alcohol Prep Pads, 70% Isopropyl Alcohol, 200/box	L	200 pcs/box 50 bx (10,000 pcs)/case

*Product packaging and configurations may differ from those shown in this catalogue. For confirmation of current shipment configurations and pricing, please contact MediVena directly.

MediVena®

Comprehensive Support Beyond the Product

At MediVena, we believe that the value of our medical products extends far beyond their use. That's why we are committed to providing ongoing educational and practical support to healthcare professionals, caregivers, and distributors alike.

To ensure the best possible user experience, we offer a wide range of educational materials, including detailed product brochures, step-by-step instructions, and clinical application guides. These resources are designed to support safe and effective product use in real-world medical environments.



Our YouTube channel features a library of instructional videos that demonstrate proper techniques, product features, and safety tips - ideal for training sessions or quick refreshers on the go. Whether you're learning how to use a safety lancet or a new insulin delivery device, our videos provide clear, practical guidance.

In addition, we provide dedicated product web pages with downloadable documents, technical specifications, FAQs, and access to cross-reference charts. These tools help users easily compare MediVena products with equivalent offerings from other leading brands, ensuring informed purchasing decisions and seamless product substitution.

For distributors and partners, we also supply marketing assets and custom support, helping you promote MediVena solutions effectively within your own markets.

Stay Connected

Follow **MediVena**[®]

Stay up to date with our latest product launches, clinical tips, and company news by following us on social media. Join the MediVena community on:

 **YouTube** – for product tutorials and demos

 **Facebook** – for updates, events, and news

 **LinkedIn** – for professional insights and industry engagement

Coming Soon:

**Advanced
Wound
Care**



**Retractable
Insulin
Syringes**



**Neonatal
Safety
Lancets**



**TB &
Allergy
Syringes**



VISIT the **One-Care STORE**



GO TO:

www.one-care.com

One-Care[®]

SAFETY MEDICAL DEVICES

Collect Blood Sample

- Blood Collection Needles & Winged Sets
- Safety Lancets
- Consumer Lancets

Protect Healthcare Professionals and Patients

- Alcohol Prep Pads

Perform Test

- Glucose Meters
- Contactless Thermometers

Deliver Medication, Nutrition; Maintain Device Patency

- Oral and Enteral Feed Syringes
- Saline Flush Syringes

Inject Medication / Vaccination

- Pen Needles
- Insulin & TB Syringes
- Syringes & Safety Needles

www.one-care.com

MediVena[®]
INSPIRED BY INNOVATION

AMERICAS:

MediVena LLC,

1 Glenlake Parkway, Ste 700, Atlanta, GA 30328, USA

 info@medivena.com

EUROPE:

MediVena Sp. z o.o.

Al. Jana Pawla II 27, 00-867 Warsaw, POLAND