

merits

VISION SPORT

K0822/K0823

FULL SIZE



FEATURES

- True mid-wheel drive: 20.6" turning radius
- Robust independent suspension can absorb up to 250 lbs on the footplate without tipping forward
- Premium, semi-reclining ergonomic seat; up to 135°
- 10" knobby drive tires for superior traction
- Dynamic LiNX controller for convenient Bluetooth programming and diagnostics
- Standard seat sizes: 16"W x 16-18"D x 18"H, 18"W x 18-20"D x 20"H, 20"W x 18-20"D x 20"H.
Optional: 22"W x 20"D x 22"H.

SPECIFICATIONS

Weight Capacity	300 lbs
Maximum Speed	5 mph
Per-Charge Range	18 miles
Batteries	U1 - 34Ah
Ground Clearance	2"
Turning Radius	20.6"
Overall Length	34"
Overall Width	24"
Drive Wheels	10"
Front Casters	6"
Rear Casters	6"
Suspension	TruStep
Seat-to-Floor Height	21"-23" (captain), 17.5"-19.5" (pan)
Seat Width	16" , 18" or 20"
Seat Depth	16" , 17" , 18" , 19" or 20"
Backrest Height	18" or 20"
Power Seat Lift	P326D : 7.5" single-post
Controller	P326A - Dynamic LiNX P326D - PG VR2
Warranty	18 Months : Motor and Electronics 5 Year/ Limited Lifetime Frame

MODELS

Vision Sport - Captain	P326-2S C
Vision Sport - Pan	P326-2S S
Vision Sport w/ PSL - Captain	P326D
Vision Sport w/ PSL - Pan	P326D

OPTIONS

Unoccupied transit tie down kit
22"W x 20"-22"D x 22"H seat
7.5" power seat lift
Swing-away joystick
USB Charger
Elevating legrests
Oxygen tank holder
Foldable rear basket
Walker holder
Cane & crutch holder

COLOR OPTIONS



Red

Blue



Pink



Green



Champagne



Black



Merits Health Products



Meritshealthusa



Merits USA

Merits Health Products Inc.

4245 Evans Avenue, Fort Myers FL, 33901

239-772-0579 | sales@meritsusa.com

merits