

merits

ROADSTER S3

S731RS3

Three Wheel Scooter



Enhancing Life Through Mobility.

Contact your Territory Sales Manager or email us at sales@meritsusa.com

Features

- Tiller angle adjustable for ergonomic driving position and easy transferring on and off the unit
- Loop tiller with throttle lever to maneuver by palm or finger
- Transport and store easily with auto latching lockup design
- Detachable battery pack with built in charging port
- Easily disassembles into 5 pieces
- Extra charging port under dashboard for easy plug in
- No tools needed for set-up
- Adjustable seat height and armrest width
- Headlight
- Rear anti-tip wheels prevent a tip-over backwards

Options

- 40" seat belt
- Multi-accessory receiver adapter
- Oxygen tank holder
- Oxygen tank bag holder
- Crutch and cane holder
- Walker holder
- Medium or Large foldable rear basket
- Cup holder
- Safety flag
- Cover



Specifications

Product#	S731RS3
Weight Capacity	260 lbs
Drivetrain	RWD
Maximum Speed	4 mph
Per-Charge Range	up to 12 mi.
Battery Size (2 required)	12Ah, ES12 AGM
Ground Clearance	2.2"
Turning Radius	47"
Overall Length	40.5"
Overall Width	25.5"
Overall Height	36.6"
Seat-to-Ground Height	19" - 21"
Seat-to-Deck Height	16" - 18"
Seat Width	18"
Seat Depth	18"
Backrest Height from seat	14.5"
Basket Size	12.5"W x 10"H x 6"D
Weight of Heaviest Piece	29 lbs.
Total Weight	103 lbs.
Motor	DC 24V 180W
Controller	PG S-DRIVE 45A
Charger	1.8A off-board
Gradient	6°
Front Wheels	7"
Drive Wheels	8"
HCPCS Code	n/a
Color Options	Red Blue

Features MK batteries



merits



merits Enhancing life through mobility.

©Merits Health Products, Inc.
4245 Evans Avenue Fort Myers, FL 33901
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0525