

merits

ROADSTER S4 S741

Four Wheel Travel-Friendly Scooter



Enhancing Life Through Mobility.

Contact your Territory Sales Manager or email us at sales@meritsusa.com

Features

- 4-wheel design provides increased stability
- Easily disassembles into 5 pieces for transport and storage.
- 18" wide deluxe stadium seat
- Detachable battery pack
- Tiller angle adjustable for ergonomic driving position and easy transferring on and off the unit
- Loop tiller with throttle lever to maneuver by palm or finger
- Transport and store easily with auto latching lockup design
- Adjustable seat height and armrest width
- No tools needed for set-up

Options

- 40" seat belt
- Multi-accessory receiver adapter
- Oxygen tank holder
- Oxygen tank bag holder
- Crutch and cane holder
- Walker holder
- Medium or Large foldable rear basket
- Cup holder
- Safety flag
- Cover



Specifications

Product#	S741
Weight Capacity	300 lbs.
Drivetrain	RWD
Maximum Speed	4 mph
Per-Charge Range	up to 6 mi.
Battery Size (2 required)	12Ah, ES12 AGM
Ground Clearance	2.5"
Turning Radius	47"
Overall Length	40.5"
Overall Width	22"
Overall Height	36.6"
Seat-to-Ground Height	19" - 21"
Seat-to-Deck Height	16" - 18"
Seat Width	18"
Seat Depth	18"
Backrest Height from Seat	14.5"
Total Weight	107 lbs.
Weight of Heaviest Piece	36 lbs.
Motor	DC 24V, 100W
Controller	PG S-DRIVE 45A
Charger	1.8A off-board
Gradient	6°
Front Wheels	8"
Drive Wheels	8"
HCPCS Code	n/a
Color Options	 blue

Features MK batteries



merits



merits Enhancing life through mobility.

©Merits Health Products, Inc.
4245 Evans Avenue Fort Myers, FL 33901
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0525