

MICRO-TOUCH[®] Micro-Thin[™] Nitrile

E s s e n t i a l S e r i e s

Nitrile, powder-free examination glove with a thin film formulation for unrestricted tactile sensitivity

This non-sterile glove is made with a thinner formulation to provide excellent tactile sensitivity and dexterity along with exceptional comfort, fit and feel. High quality manufacturing ensures this glove will hold up for a variety of clinical procedures.

KEY FEATURES AND BENEFITS:

- Excellent tactile sensitivity is ideal for palpating veins and starting IVs
- The textured surface of the glove promotes sure handling of materials and instruments
- 300 gloves per box reduces restocking and reorders

Recommended For

- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Protection from Type I latex allergy in HCW's or patients
- Brief, low-risk examination procedures
- Environmentally friendly facilities

FEATURED TECHNOLOGIES



MICRO-TOUCH[®] Everyday, Everywhere Protection

MICRO-TOUCH[®] premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH[®] Micro-Thin[™] Nitrile

Essential Series

PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	Blue	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.06 / 2.4
Finger Thickness (mm/mil)	0.09 / 3.5	Allergy Prevention	Latex (Type I)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Sterile	No	Antistatic	Not Tested
Glove Length (mm/inches)	230 / 9.1	Product Segmentation	Essential Series

ORDERING INFORMATION

Size	L (8.5 - 9)	M (7.5 - 8)	S (6.5 - 7)	XL (9.5 - 10)	XS (5.5 - 6)
Ansell Product Code	6034313	6034312	6034311	6034314	6034310

PACKAGING AND STORAGE

Packaging	300 gloves per dispenser 10 dispensers per case 3000 gloves per case
Shelf Life	5 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

ASTM D6319, Category III, EN 420:2003 + A1:2009, EN 455 1-4, EN ISO 374-1:2016, EN ISO 374-5:2016, FDA21 CFR 177-2600, ISO 10993-10:2002, ISO 9001, ISO 13485

Sensitisers

Zinc Dibutylthiocarbamate (ZDBC)

Zinc Diethylthiocarbamate (ZDEC)

Only a very small number of users may be sensitized to this ingredient in the glove and hence may develop irritant and/or allergic contact reactions.

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email:
insidesalesus@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email:
luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email:
apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email:
protection@ap.ansell.com