

## Evolution Series • Model EVOx-U1F-ADA-BB

### Dimensions & Specifications

|                              |                         |
|------------------------------|-------------------------|
| Exterior Size WxDxH (Inches) | 24" x 28" x 32"         |
| Interior Size WxDxH (Inches) | 19.5" x 21.5" x 17.375" |
| Hinged Solid Door            | 1                       |
| Interior Capacity (Cu/Ft)    | 4.2                     |
| Stainless Steel Drawers      | (2) Full Extension      |
| Shipping Weight              | 175 lbs.                |

### Refrigeration System

|                      |                             |
|----------------------|-----------------------------|
| Refrigerant          | R-404a                      |
| Refrigerating Cycle  | -30° C ~ -15° C             |
| Refrigeration System | Expansion Valve, Forced Air |
| Defrost Method       | Automatic                   |

### Electrical System

|              |                     |
|--------------|---------------------|
| Horse Power  | 1/3+                |
| Power Supply | 115/60/1            |
| Amps         | FLA (8.0) RLA (6.1) |
| BTU          | 1100 (REV 3.1)      |
| NEMA Plug    | 5-15P               |

### Controller Specifications

|                        |               |
|------------------------|---------------|
| Digital Microprocessor | LCD Screen    |
| Touch Key Pad          | Multi-program |
| High/Low Visual Alarm  | Yes           |
| High/Low Audible Alarm | Yes           |
| Migali® IntelVU®       | Optional*     |

### Construction

|                       |                                 |
|-----------------------|---------------------------------|
| Exterior              | 304 Series Stainless Steel      |
| Interior              | 304 Series Stainless Steel      |
| (1) Hinged Solid Door | 304 Series Stainless Steel Door |
| Heavy-Duty Hinges     | Chrome, Side-mount              |
| Positive Latch Handle | Chrome w/ lock and key          |

### Other Features

|                              |                             |
|------------------------------|-----------------------------|
| Access Port                  | Right Side, 2.125" Diameter |
| Product Simulator            | Glycol & Probe Included     |
| Migali® Paper Chart Recorder | 4", 7 Day                   |
| Monitoring Alarm             | Hi/Lo Alarm w/ Dry Contacts |
| Leveling Feet                | (4) 1/2-13 rubber feet      |

### Warranty Protection Plan

|            |         |
|------------|---------|
| Compressor | 7 years |
| Parts      | 2 years |
| Labor      | 1 year  |

### Certifications

|                         |        |
|-------------------------|--------|
| ETL Electrical-Intertek | cETLus |
|-------------------------|--------|

### Migali® IntelVU®

\*Migali® IntelVU®: Touch Screen Microprocessor, 24/7 temperature monitoring, downloadable via USB port, full control over refrigerator performance with many additional features. See Migali® IntelVU® spec sheet.



Made in the USA. Trusted Since 1955


MIGALI® SCIENTIFIC

Toll Free: 855-464-4254 Fax: 856-963-3604

Email: [contact@migaliscientific.com](mailto:contact@migaliscientific.com) Web: [www.migaliscientific.com](http://www.migaliscientific.com)



Evolution Series • Model EVOx-U1F-BB-ADA



Made in the USA. Trusted Since 1955

MIGALI® SCIENTIFIC

Toll Free: 855-464-4254 Fax: 856-963-3604

Email: [contact@migaliscientific.com](mailto:contact@migaliscientific.com) Web: [www.migaliscientific.com](http://www.migaliscientific.com)

