

Genesis Series • Model G-2RG-HC

Dimensions & Specifications

Exterior Size WxDxH (Inches)	51.73" x 31.69" x 81.3"
Interior Size WxDxH (Inches)	47" x 25" x 61"
Hinged Solid Door	2
Interior Capacity (Cu/Ft)	42.2
Powder Coated Shelves	10
Shipping Weight	454 lbs.

Refrigeration System

Refrigerant	R290
Refrigerating Cycle	+ 2°C ~ + 8° C
Refrigeration System	Capillary Tube, Forced Air
Defrost Method	Automatic

Electrical System

Horse Power	1/5
Power Supply	115/60/1
Amps	3.2
kWh	4.66
NEMA Plug	5-15P

Controller Specifications

Digital Microprocessor	LCD Screen
Touch Key Pad	Multi-program
High/Low Visual Alarm	Yes
High/Low Audible Alarm	Yes

Construction

Exterior	304 Series Stainless Steel
Interior	304 Series Stainless Steel
(2) Hinged Glass Door	ABS Extruded Frame, 2 Pane Glass
Hinges	Torque-rod, Flush Mount
Door Handle	Recessed Handle

Other Features

Access Port	Top Side, 0.875" Diameter
Product Simulator	Glycol & Probe Included
Casters	(2) Locking, (2) Non-Locking

Warranty Protection Plan

Compressor	5 years
Parts	1 years
Labor	1 year
Optional Extended Warranty	Total 2 Years Parts, 2 Years Labor

Certifications

ETL Electrical-Intertek	cETLus
ISO	9000, 14000



Trusted Since 1955.

MIGALI® SCIENTIFIC




Toll Free: 855-464-4254 Fax: 856-963-3604

Email: contact@migaliscientific.com

Web: www.migaliscientific.com



Genesis Series • Model G-2RG-HC



Trusted Since 1955.

MIGALI® SCIENTIFIC

Toll Free: 855-464-4254 Fax: 856-963-3604

Email: contact@migaliscientific.com

Web: www.migaliscientific.com

