




Mobile Cart with Integrated Compressor



The portable oil-free compressor is suitable for any dental, veterinary, and hygiene practice.



Features:

- Fully Integrated & Portable Mobile Cart Delivery System
- Clean Water Bottle for Self Contained System
- 3 HP Tubings for High and Low Speed Handpieces
- 120V Plug in for Compressor
- Ultra-Quiet Oilless Compressor – Low Maintenance
- 1 HP Compressor
- Great for Hygienist, Mobile Clinics, Vet Clinics
- NO Additional Plumbing Required

COMPRESSOR

SPECIFICATIONS

Noise Level	60db
Tank Size	2.0 Gallon Aluminum Tank (Rust Free)
Motor	1.0HP (Rated/Running -2.0HP(PEAK))
Maximum Pressure	120 psi
CFM	3.10 CFM @ 40 PSI 2.30 CFM @ 90 PSI



**QUIET & DURABLE
OIL-FREE
COMPRESSOR**

