

Product Name:	Cell Culture Dish	
Product Application:	NEST cell culture dishes are ideal for all types of culturing where microscopic observation is required. Numeric marks on the bottom of each dish allow users to identify the location of cells.	
Product Materials:	Polystyrene (PS) Meets USP, Class VI standards	
product Specifications:	Temperature range: stored at room temperature Shelf life: 3 Years after date of production (ensure package is in good)	
Sterilization:	Yes. Sterilized by E-beam, Sterility Assurance Level: SAL=10 ⁻⁶ . The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137	
Pyrogens:	Non-Pyrogenic	
RNase/DNase:	DNase/RNase free	
BSE/TSE:	These products are deemed animal free.	
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: qualified Flatness test: qualified Hydrophilic test: qualified Cell culture test: qualified Packaging inspection: qualified	
Features	1. High clarity polystyrene 2. Flat transparent surface for distortion-free observation 3. TC treated: Vacuum Plasma tissue culture treatment, excellent cell adherence 4. Stackable for easy storage and handling	

Product Range

Cat. No.	Spec		Cell Growth Area (cm ²)	Recommended Medium Volume (mL)	/Pack	/Case
	Type	Height (mm)				
706001	35	12	8.5	1.8-2.7	20	500
706201	35(With Gripping Ring)	12	8.5	1.8-2.7	20	500
705001	60	15	22.9	4.2-6.3	20	500
705201	60(With Gripping Ring)	15	22.9	4.2-6.3	20	500

Head Office No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China
Tel: +86+ 510-6800 6788 Email: info@nest-wuxi.com
Online: www.cell-nest.com

Oversea NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)