

Product Name:	Cell Culture Flask	
Product Application:	NEST cell culture flasks are an ideal solution for enhanced cell attachment and growth of a variety of primary cells and transformed cells in serum free or serum-containing cultures. These flasks are available as tissue culture treated or non-treated as well as with a vent cap or plug seal cap to meet your requirements.	
Product Materials:	PS (body) + PP (cap) Meets USP, Class VI standards	
product Specifications:	Temperature range: stored at room temperature Shelf life: 3 Years after date of production (ensure package is in good).	
Sterilization:	Yes. Sterilized by E-beam, Sterility Assurance Level: SAL=10⁻⁶. The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137	
Pyrogens:	Non-Pyrogenic	
RNase/DNase:	DNase/RNase free	
BSE/TSE:	These products are deemed animal free.	
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: qualified Particle Test: qualified Hydrophilic test: qualified Cell culture test: qualified Packaging inspection: qualified	
Features	1. High clarity medical grade polystyrene material. 2. Two types of cap style: vent caps & plug seal cap. 3. Vent caps with 0.22um hydrophobic filters for gas exchange without contamination. 4. TC treated: Vacuum Plasma tissue culture treatment, excellent cell adherence. 6. Multiple size options: cell growth areas ranging from 25cm² to 225cm² are available. 7. A brand new specification, T150 U-Shaped Canted Neck Cell Culture Flask: • Without influencing the regular usage, the bottleneck is enlarged, allowing the pipettes or cell scrapers to reach every corner of the flask freely and flexibly. • The bottleneck is enlarged to the size of the large mouth of NEST BioFactory, facilitating accessories' matching and offering better tubes' connection and closed transfer. • The unique U-shaped body is easier to hold and enhances the utilization rate of the bottom.	