

Product Range

3-layer Cell Culture Flask


Cat. No.	Growth area	Cap Style	Size (mm)			TC Treated	/Pack	/Case
			Length	Width	Height			
731301	520 cm ²	Plug Seal Cap	203	120.5	60.1	Yes	1	12
731302		Vent Cap						

5-layer Cell Culture Flask

Cat. No.	Growth area	Cap Style	Size (mm)			TC Treated	/Pack	/Case
			Length	Width	Height			
731001	870 cm ²	Plug Seal Cap	203.6	120.5	84.3	Yes	1	8
731002		Vent Cap						

Technical Drawing

3-layer Cell Culture Flask



Head Office

No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi,
Jiangsu, China
Tel: +86+ 510-6800 6788 Email: info@nest-wuxi.com
Online: www.cell-nest.com

Overseas

NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab
Emirates)