

## Nitro® 3-Wheel Rollator

HCPCS: E0143,E0143



FIGURE A



FIGURE B



FIGURE C



FIGURE D




FIGURE E



FIGURE F



FIGURE G



FIGURE H

test The Nitro 3-Wheel Rollator is the newest addition to our premium line of Nitro rollators. Featuring a hidden inner brake cable, the functional and modern Nitro 3 is designed with safety in mind while indoors or outdoors. The Nitro 3 easily navigates narrow spaces and folds flat while in the standing position so it can roll, park and store anywhere. Also included is a large, removable pouch to keep personal items close. The Nitro 3 comes pre-assembled, so you can unbox, unfold and get your customers going!

- Stylish and modern design features lightweight aluminum frame with built-in brake cable for extra safety and protection
- Ergonomic hand grips provide comfort and control, push down on handbrakes to engage parking brake feature
- Quickly and easily adjust the height of the handles with the push-button height adjustment
- Large, removable storage pouch is perfect for your personal belongings, or transporting items to and from
- Easy side to side fold allows the rollator to stay in the standing position, making it easy to roll, park and store anywhere
- Dual front fork design enhances turning radius while providing added safety
- Large 9" front wheel allows for optimal steering and conveniently moves over most kinds of terrain, whether indoors or out
- Out of box assembled
- Patent RE49,916

| PRODUCTS               | Description                       |
|------------------------|-----------------------------------|
| RTL10266TWHL           | Nitro 3-Wheel Rollator, Red, 1/EA |
| SPECIFICATIONS         | Description                       |
| Warranty               | Limited Lifetime                  |
| Product Weight         | 13.5 lb                           |
| Handle Height          | 32" - 37"                         |
| Weight Capacity        | 300 lb                            |
| Carton Shipping Weight | 17.9 lb                           |
| Dimensions             | 25" (W) x 29" (D)                 |
| Wheel Size             | 9.8" Front   7.5" Back            |
| Width Between Arms     | 16"                               |