

Nitro Aluminum Rollator, Tall Height, 10" Casters

HCPCS: E0143 E0156-1,E0143 E0156-1



- Tall height ideal for users 6'0" and above
- Attractive, Euro-style design
- Brake cable inside frame for added safety
- Handle height easily adjusts with unique push button
- Back support height easily adjusts with tool-free thumb screw
- Caster fork design enhances turning radius
- Large 10" front casters allow optimal steering and rolling comfort
- Easily folds with one hand to ultra-compact size for storage
- Plastic clip holds the Nitro closed when frame is folded
- Cross-brace design allows for side-to-side folding and added stability
- Seat is durable and comfortable
- Handsome and removable zippered storage bag with unique attachment to keep bag securely in place when Nitro is open or folded
- Lightweight, aluminum frame

PRODUCTS

Description

RTL10266BK-T	Nitro Rollator,Tall, Black frame,1EA
RTL10266-T	Rollator NITRO size L Color: red

SPECIFICATIONS

Description

Warranty	Frame: Limited Lifetime Brake and Handle Housing: 5 Year Limited
Seat Dimensions	18" (W) x 8.5" (D)
Handle Height	36" - 41"
Weight Capacity	300 lb
Seat To Floor Height	23.5"
Carton Shipping Weight	30 lb
Weight	18 lb

ACCESSORIES

Description

10266-CH	Nitro Cup Holder, 1/ea, EA
----------	----------------------------