

Pediatric Viper Plus Reclining Wheelchair

HCPCS: E1236 E1226-1,E1236 E1226-1



FIGURE A




FIGURE B

- State-of-the art hydraulic reclining mechanism
- Backrest reclines from 90° to 180°
- Double-coated chrome carbon steel frame for maximum durability
- Durable flame retardant nylon upholstery
- Padded flip-back detachable desk armrests
- Composite lightweight 20" Mag-style rear wheels with urethane tires and 6" front rubber tires
- Offset rear wheel position provides additional stability and evenly distributes user's weight in all reclining positions
- Precision sealed front and rear wheel bearings
- Front caster forks are adjustable in three positions
- Standard dual axle provides easy transition of seat to hemi-level
- Push-to-lock wheel locks
- Spreader bar functions as a back stabilizer and a push handle for caregiver
- Standard with height adjustable swing-away elevating leg rests
- Durable and lightweight composite footplates
- Flip-up anti-tippers with wheels

| PRODUCTS | Description |
|------------------------|--|
| PL412RBDDA | ViperPlusRecl,12"DetDskArm,1/cs |
| PL414RBDDA | Viper,14"Recl.FlpBkDsk1/cs |
| SPECIFICATIONS | Description |
| Warranty | Limited Lifetime |
| Product Weight | 42 lb 43 lb |
| Product Dimensions | 38" (L) w/ riggings x 20" 22.5" (W) x 40.5" (H) |
| Seat Dimensions | 12" 14" (W) x 12" 14" (D) |
| Weight Capacity | 250 lbs |
| Seat To Floor Height | 15", 16" |
| Carton Shipping Weight | 48 lb 49 lb |
| Folded Dimensions | 38" (L) w/ riggings x 13" 13.5" (W) x 40.5" (H) |
| Back Height | 26" |
| CA Prop 65 | [WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov] |