

Poly-Fly High Strength, Lightweight Wheelchair/Flyweight Transport Chair Combo

HCPCS:



FIGURE A



FIGURE B



FIGURE C




FIGURE D

- Hybrid wheelchair that can be used as both self-propelled wheelchair or a transport chair standard
- Quick-release 24" wheels can be removed to transition self-propelled chair to a transport chair
- All-aluminum frame
- Comes with two sets of installed wheel locks: One for use with manual wheelchair configuration, one for use with transport chair configuration
- Back folds down with deluxe back release
- Nylon upholstery
- Comes standard with swing-away footrests with composite foot plates and heel loops
- 8" casters
- Anti-tippers included
- Compatible with optional Swing-away elevated leg rests (Material# STDELR-TF - Swing-away Elevated Leg rest, Tool-free)
- Weight Capacity: 250 lb

PRODUCTS	Description
TR18	PolyFlyW/C/TrsPrt,w/SwgFt,18"AntiTp,1/cs
TR20	PolyFlyWC/TPrt,20w/SgwFt.AntiTp,1/cs
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Product Weight	34 lb 36 lb w/ 24" wheels
Product Dimensions	42" (L) w/ riggings x 24.5" 25.5" (W) x 37" (H) w/ 24" wheels
Seat Dimensions	18" 20" (W) x 16" (D)
Weight Capacity	250 lb
Seat To Floor Height	19.5"
Carton Shipping Weight	39 lb 41 lb
Folded Dimensions	42" (L) w/ riggings x 12" (W) x 37" (H) w/ 24" wheels
Back Height	16"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]