

# Protekt® Aire 3500

Low Air Loss/Alternating Pressure Mattress System  
with 3" Densified Fiber Support Base

HCPSC  
E0277



## FEATURES AND BENEFITS:

### PROTEKT® AIRE 3000 PUMP

- Low air loss and alternating pressure therapies in one system.
- Air flow dial offers custom patient comfort and pressure levels.
- Static function for safe and easy patient transfer.
- Preset 10 minute alternating cycle.
- Visual low pressure alarm indicates air flow may be compromised.

|         |                                      |
|---------|--------------------------------------|
| 83500   | Protekt® Aire 3500 36"x80"x8"        |
| 83500RR | Protekt® Aire 3500 with Raised Rails |

### MATTRESS & COVER

- 5" mattress with 3" densified fiber base provides a safety support base in case of power failure.
- 18 laser holes keeps the patient cool and dry.
- Nylon top cover is fluid resistant and quilted for added comfort.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Optional raised side rails.
- Weight Capacity: 350 lbs.
- Warranty: 18 Month Non-Prorated.



270 Washington Street, Mount Vernon, NY 10553

**Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005**

[www.proactivemedical.com](http://www.proactivemedical.com) • [sales@proactivemedical.com](mailto:sales@proactivemedical.com)