

Protekt® Aire 4600DXAB

Low Air Loss/Alternating Pressure Mattress System with Digital Pump, Raised Air Bolsters and Cell-On-Cell Support Base



FEATURES AND BENEFITS:

PROTEKT® AIRE 4000DX PUMP

- Extended Warranty Digital Pump with Auto Compressor Sleep Technology (ACST) shuts off compressor at desired pressure level extending pump life cycle.
- Low air loss and alternating pressure therapies in one system.
- Adjustable weight settings allows pressure to be customized according to each patient's therapeutic requirements.
- Preset 10 minute alternating cycle.
- Static function for safe and easy patient transfer.
- Audio/visual low pressure alarm indicates air flow may be compromised.

84600DXAB Protekt® Aire 4600DXAB 36"x80"x8"

MATTRESS & COVER (Available in 84" length)

- Raised side perimeter air bolsters mitigates accidental patient roll-outs, and reduces the risk of entrapment.
- Center cut-out in air bolsters enables safe ingress and egress from the mattress.
- Unique cell-on-cell design ensures that the mattress will hold air in the lower cells in the event of a power failure.
- 18 laser holes keeps the patient cool and dry.
- Stretch nylon top cover is fluid resistant and vapor permeable.
- Static head section prevents unnecessary head/neck movement.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Weight Capacity: 400 lbs.
- Warranty: 2 Year Non-Prorated.



270 Washington Street, Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com