

# Protekt® Aire 6000AB

Low Air Loss/Alternating Pressure Mattress  
System with Deluxe Digital Pump, Raised Air Bolsters and  
Cell-On-Cell Support Base



## FEATURES AND BENEFITS:

### PROTEKT® AIRE 6000 PUMP

- Dual therapies; low air loss and alternating pressure.
- 4 alternating cycles 10, 15, 20 & 25 minutes allows clinicians to customize therapy.
- Seat-inflate feature helps protect the patient when they are in the fowler position.
- Max firm mode for quick inflation.
- Static mode provides firm support for egress or ingress from mattress or side edge seating.
- 3 audible and visual alarms; low pressure, power failure and alternating failure.

### MATTRESS & COVER (Available in 84" length)

- Raised side perimeter air bolsters mitigates accidental patient roll-outs, and reduces the risk of entrapment.
- Center cut-out in air bolsters enables safe ingress and egress from the mattress.
- Unique cell-on-cell design ensures that the mattress will hold air in the lower cells in the event of a power failure.
- 18 laser holes keeps the patient cool and dry.
- Stretch nylon top cover is fluid resistant and vapor permeable.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Weight Capacity: 450 lbs.
- Warranty: 2 Year Non-Prorated.

80060AB	Protekt® Aire 6000AB 36"x80"x8"
---------	---------------------------------



590 Franklin Ave. | Mount Vernon, NY | 10550

**Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005**

[www.proactivemedical.com](http://www.proactivemedical.com) • [sales@proactivemedical.com](mailto:sales@proactivemedical.com)