

Protekt® Aire 6400

Low Air Loss/Alternating Pressure Mattress System w/ Deluxe Digital Pump



FEATURES AND BENEFITS:

PROTEKT® AIRE 6400 PUMP

- Dual therapies; low air loss and alternating pressure.
- 4 alternating times of 10, 15, 20 & 25 minutes allows clinicians to customize therapy.
- Seat-inflate feature helps protect the patient when they are in the fowler position.
- Max firm mode for quick inflation.
- Static mode provides firm support for egress or ingress from mattress or side edge seating.
- 3 audible and visual alarms; low pressure, power failure and alternating failure.

MATTRESS & COVER

- 18 laser holes keeps the patient cool and dry.
- Nylon top cover is fluid resistant and quilted for added comfort.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Available in 84" length.
- Optional raised side rails.
- Weight Capacity: 450 lbs.
- Warranty: 2 Year Non-Prorated.

86400	Protekt® Aire 6400 36"x80"x8"
86400RR	Protekt® Aire 6400 with Raised Rails



590 Franklin Ave. | Mount Vernon, NY | 10550

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com