

Protekt® Aire 8000

Low Air Loss/Alternating Pressure Bariatric Mattress System

HCPCS
E0277



FEATURES AND BENEFITS:

- Alternating pressure cycles of 10, 15, 20 & 25 minutes allows clinicians to customize pressure redistribution and skin surface interface.
- Cell-on-cell mattress design prevents patient from bottoming out in case of power failure.
- 10" air cells are constructed of Nylon/PU for added comfort, support and durability.
- Low air loss function controls and manages moisture to prevent skin maceration.
- Audio/visual low pressure alarms alerts clinicians when preset pressure levels decrease.
- Static mode allows for easy transfer in and out of bed.
- Seat inflate increases air flow by 20% to middle bladders to prevent bottoming out.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Quilted top cover for added patient comfort.

- Lock out key completely locks the control panel to eliminate tampering of the patient settings, including power switch.
- CPR pull tab quickly releases air from mattress.
- Optional raised side rails available on all sizes.

SPECIFICATIONS

- Air Flow: 18 liters per minute.
- Mattress Dimensions: 42"x80"x10".
- Available in 48", 54" and 60" widths.
- Weight Capacity: 750 lbs.
- Warranty: 2 Year Non-Prorated.

80080	Protekt® Aire 8000BA-42
80080RR	Protekt® Aire 8000BA-42 with Raised Rails
80085	Protekt® Aire 8000BA-48
80080-54	Protekt® Aire 8000BA-54
80080-60	Protekt® Aire 8000BA-60



590 Franklin Ave., Mount Vernon, NY 10550

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com