

# Protekt® Aire 8600AB

## Low Air Loss/Alternating Pressure Bariatric Mattress System with Raised Air Bolsters



### FEATURES AND BENEFITS:

- Raised side perimeter air bolsters mitigates accidental patient roll-outs, and reduces the risk of entrapment.
- Alternating pressure cycles of 10, 15, 20 & 25 minutes allows clinicians to customize pressure redistribution and skin surface interface.
- Center cut-out in the air bolsters enables safe ingress and egress from the mattress.
- Cell-on-cell mattress design prevents patient from bottoming out in case of power failure.
- Air cells are constructed of Nylon/PU for added comfort, support and durability.
- 18 laser air holes to keep the patient cool and dry.
- Audio/visual low pressure alarms alerts clinicians when preset pressure levels decrease.
- Static mode allows for easy transfer in and out of bed.
- Seat inflate increases air flow by 20% to middle bladders to prevent bottoming out.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- Lock out key completely locks the control panel to eliminate tampering of the patient settings.
- Stretch nylon top cover is fluid resistant and vapor permeable.
- CPR pull tab quickly releases air from mattress.

### SPECIFICATIONS

- Air Flow: 18 liters per minute.
- Mattress Dimensions: 42"x80"x8" & 48"x80"x8".
- Weight Capacity: 750 lbs.
- Warranty: 2 Year Non-Prorated.

86080AB-42	Protekt® Aire 8600AB w/Raised Air Bolsters 42"
86080AB-48	Protekt® Aire 8600AB w/Raised Air Bolsters 48"



590 Franklin Ave., Mount Vernon, NY 10550

**Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005**

[www.proactivemedical.com](http://www.proactivemedical.com) • [sales@proactivemedical.com](mailto:sales@proactivemedical.com)