

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 7/1/2026

Caribou

mm (inch)	kg (lb)				
					
 Overall width		510 mm (20")	570 mm (22½")	650 mm (25½")	750 mm (29½")
 Length		750 mm (29½")	920 mm (36¼")	1110 mm (43¾")	1110 mm (43¾")
 Overall height, vertical		460–570 mm (18–22½")	530–640 mm (20¾–25¼")	680–790 mm (26¾–31")	880–990 mm (34¾–39")
 Overall height, horizontal		700 mm (27½")	700 mm (27½")	800 mm (31½")	800 mm (31½")
 Angling		0°–87°	0°–87°	0°–87°	0°–87°
 Foot plate angle		(-20°)–(20°)	(-20°)–(20°)	(-20°)–(20°)	(-20°)–(20°)
 Total mass		20.5 kg (45.1 lb)	24 kg (52.9 lb)	30.5 kg (67.2 lb)	35 kg (77.1 lb)
 Maximum rated load		62.5 kg (137.7 lb)	86 kg (189.6 lb)	112.5 kg (248 lb)	137 kg (302 lb)
 Total mass, with power		22.5 kg (49.6 lb)	25.5 kg (56.2 lb)	32.5 kg (71.7 lb)	37 kg (81.6 lb)
 Maximum rated load, with power		64.5 kg (142.1 lb)	87.5 kg (192.9 lb)	114.5 kg (252.4 lb)	139 kg (306.4 lb)
 Maximum user mass		40 kg (88.1 lb)	60 kg (132.2 lb)	80 kg (176.3 lb)	100 kg (220.4 lb)
 User height		850–1120 mm (33½–44")	980–1290 mm (38½–50¾")	1250–1560 mm (49¾–61½")	1530–1840 mm (60¾–72½")
 Static stability downhill		max. 5°	max. 5°	max. 5°	max. 5°
 Static stability uphill		max. 5°	max. 5°	max. 5°	max. 5°
 Sidewards stability sideways		max. 5°	max. 5°	max. 5°	max. 5°

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1231.4 Pricing and specifications subject to change without notice.

Caribou Base

Required selections. Choose from Standard or Power

Caribou Base

Caribou Base, Size 1 (860101-FP)	\$1487
Caribou Base, Size 2 (860102-FP)	\$1534
Caribou Base, Size 3 (860113-FP)	\$1709
Caribou Base, Size 4 (860114-FP)	\$2016

Caribou Power Base

Power Actuator recommended for Size 3 and 4

Caribou Base for Power Actuator Options, Size 3 (860103-ELFP)	\$3620
Caribou Base for Power Actuator Options, Size 4 (860104-ELFP)	\$3953



Caribou Mast with Torso - Lower Body Pad

Required Selection

Lower Body Pad, Size 1 (X860101-P)	\$1276
Lower Body Pad, Size 2 (X860102-P)	\$1315
Lower Body Pad, Size 3 (X860113-P)	\$1467
Lower Body Pad, Size 4 (X860114-P)	\$1729


Multi Adjustable Footrests

Required Selection

Multi Adjustable Footrest, Size 1 (X860101-F)	\$443
Multi Adjustable Footrest, Size 2 (X860102-F)	\$477
Multi Adjustable Footrest, Size 3 (X860113-F)	\$505
Multi Adjustable Footrest, Size 4 (X860114-F)	\$534

Top Pad

Chest Pad is a required selection




Chest Pad, Size 1, A: 6", B: 7", C: 7.75" (861101-1)	\$450
Chest Pad, Size 2, A: 7.5", B: 9", C: 9.75" (861101-2)	\$475
Chest Pad, Size 3 A: 8.75", B: 10.5", C: 12.5" (861101-3)	\$500
Chest Pad, Size 4, A: 9.75", B: 12.25", C: 15.75" (861101-4)	\$525

Extra Cushion


Choose the same size as the Caribou Base, then state desired quantity in the box below

Extra Cushion



Extra Cushion, Size 1, A: 9.75", B: 12.25", C: 15.75" (861101-4) (Each)	\$295
Extra Cushion, Size 2: W: 11.75" x H: 4" (Each) (861106-2)	\$295
Extra Cushion, Size 3, W: 13.75" x H: 4" (Each) (861106-3)	\$325
Extra Cushion, Size 4, W: 15.75" x H: 4" (Each) (861106-4)	\$325

Supports and Wedges



Wedges

Wedges for Side Slots, Size 1 (Pair) (861107-1)	\$22
Wedges for Side Slots, Size 2 (Pair) (861107-2)	\$22
Wedges for Side Slots, Size 3 (Pair) (861107-3)	\$22



Back Supports


Back Supports, Size 1 and 2 (Each) (861115-2)	\$350
Back Supports, Size 3 and 4 (Each) (861115-3)	\$375



Upper Arm Supports

Can not be used with shoulder supports. Can not be used with Lateral Side Supports.

Upper Arm Support, Small (861109-10)	\$450
Upper Arm Support, Large (861109-20)	\$485



Shoulder Support


Must match Chest Pad size

Shoulder Support, Size 1 (Pair) (861126-1)	\$300
Shoulder Support, Size 2 (Pair) (861126-2)	\$315
Shoulder Support, Size 3 (Pair) (861126-3)	\$330
Shoulder Support, Size 4 (Pair) (861126-4)	\$350

Vest Options

Measuring Shoulder to Pelvis Diagonally. See pictural under Belts and Straps. Vests require fixlocks for mounting.

Cross Vest




Cross Vest, XS, S-P: 7.75" - 12.5" (Each) (85720-XS)	\$220
Cross Vest, Small, S-P: 10.75" - 12.5" (Each) (85720-S)	\$220
Cross Vest, Medium, S-P: 12.5" - 16" (Each) (85720-M)	\$220
Cross Vest, Large, S-P: 16" - 18.5" (Each) (85720-L)	\$220
Cross Vest, XL, S-P: 20.5" - 23.5" (Each) (85720-XL)	\$220

Torso Vest



Torso Vest, XS, S-P: 7.75" - 10.75" (Each) (85705-XS)	\$232
Torso Vest, Small, S-P: 10.75" - 12.5" (Each) (85705-S)	\$232
Torso Vest, Medium, S-P: 12.5" - 16" (85705-M)	\$232
Torso Vest, Large, S-P: 16" - 18.5" (Each) (85705-L)	\$232
Torso Vest X-Large - S-P: 18.5" - 20.5" (E0960) (85705-XL)	\$232
Torso Vest, XXL, S-P: 20.5" - 23.5" (Each) (85705-XXL)	\$232

Vest Mounting Hardware




Fix Locks

Required selection for all vests

Fix-Locks, Black, Bag of 6 Pcs (862109-35)	\$46
--	------

Side Supports



Side Supports, Fixed

Supports, Fixed, Size 1 and 2, A: 4.75" x B: 4.25" (Pair) (861102-10)	\$273
Side Supports, Fixed, Type 5, Size 3 and 4, A: 6" x 6.25" (Pair) (861102-20)	\$328



Removable Lateral Supports

Select desired size(s) state desired quantity in the corresponding box below

Removable Lateral Supports, Size 1, Bar L: 5.75" Pad H: 5.5" Pad W: 3.75" (Pair) (861102-1)	\$425
--	-------

Removable Lateral Supports, Size 2, Bar L: 6.75" Pad H: 7" Pad W: 5" (Pair) (861102-2)	\$425
---	-------

Removable Lateral Supports, Size 3, Bar L: 7.75" Pad H: 7" Pad W: 5" (Pair) (861102-3)	\$425
---	-------

Belts and Straps



Cross Straps

Cross Strap, Size 1, 15" - 19.25" (Pair) (862102-1)	\$305
Cross Strap, Size 2, 17" - 21.5" (Pair) (862102-2)	\$305
Cross Strap, Size 3, 20.75" - 27.5" (Pair) (862102-3)	\$322
Cross Strap, Size 4, 24.75" - 33.75" (Pair) (862102-4)	\$322

Chest Strap

If multiple straps are required, please add to notes section below or contact customer service



Chest Strap and Pelvic Strap, Small, W: 16.5" - 24.75" (Each) (862101-2)	\$135
Chest Strap and Pelvic Strap, Medium, W: 19" - 30.25" (Each) (862101-3)	\$135
Chest Strap and Pelvic Strap, Large, W: 22.75" - 37.75" (Each) (862101-4)	\$135

Adjustable Plate for Fixed Side Supports

Adjustable plate used with fixed side support to increase width

Adjustable Plate for Fixed Supports (Pair) (861102-PL)	\$155
---	-------

Straight Mounting Bracket

To be used on the extra cushion.

Straight Mounting Bracket for Extra Pad (861106-Kit)	\$125
---	-------

Lateral Side Supports


For use with the Removable Lateral Supports. Select size then state desired quantity in the box below.



Lateral Support Strap, Small: 10-12.5" (Each) (862103-1)	\$75
Lateral Support Strap, Medium: 13-16" (Each) (862103-2)	\$75
Lateral Support Strap, Large: 13.5-9.75" (Each) (862103-3)	\$75
Lateral Support Strap, X-Large: 17.75-23.5" (Each) (862103-4)	\$75

Anatomic Head Support

Cover and Mounting Hardware required. Other headrest options available, please see Head supports order form. Please verify Compatibility.




Anatomic Head Support

Anatomic Head Support, Size Small, A: 5.5" x B: 5.5" x C: 9.5" (5099001)	\$285
Anatomic Head Support, Size Medium, A: 6" x B: 5.5" x C: 10" (5099002)	\$285
Anatomic Head Support, Size Large, A: 6.5" x B: 5.5" x C: 12.5" (5099003)	\$285

Anatomic Head Support Cover

If Anatomic Head Support is selected, you must select a cover. Does not come standard.

Head Support Cover, Small (9923-83)	N/C
Head Support Cover, Medium (9945-83)	N/C
Head Support Cover, Large (9967-83)	N/C




Headrest Mounting Hardware

Required with All Head Supports

Headrest Mounting Hardware, Short (98619-25)	\$250
Headrest Mounting Hardware, Long (98619-26)	\$256

Pommel



Pommel

Pommel, Supine, A: 2", B: 2.75", C: 6", D: 3.5" (861103-1)	\$150
Pommel, Prone, A: 6", B: 2.75", C: 6", D: 6.75" (861103-2)	\$195



Straps for Pommels

Straps for Pommel, Size 1, L: 9.84" - 11.81" (862104-1)	\$175
Straps for Pommel, Size 2, L: 10.63" - 18.5" (862104-2)	\$175


Head Pad



Head Pad, Size 1 and 2 (861116-1)	\$575
Head Pad, Size 3 and 4 (861116-2)	\$575

Knee Supports

Knee Support, Fixed



Knee Support, Standard Pad, Fixed Bracket, Size 1 and 2 (861144-2)	\$575
Knee Support, Hinged/Split Pad, Fixed Bracket, Size 1 and 2 (861144-20)	\$695
Knee Support, Standard Pad, Fixed Bracket, Size 3 and 4 (861144-3)	\$575
Knee Support, Hinged/Split Pad, Fixed Bracket, Size 3 and 4 (861144-30)	\$695

Knee Support, Swing-Away




Knee Support, Standard Pad, Swing Away Bracket, Size 1 and 2 (861154-2)	\$855
Knee Support, Hinged/Split Pad, Swing Away Bracket, Size 1 and 2 (861154-20)	\$1025
Knee Support, Standard Pad, Swing Away Bracket, Size 3 and 4 (861154-3)	\$855
Knee Support, Hinged/Split Pad, Swing Away Bracket, Size 3 and 4 (861154-30)	\$1025

Foot Positioning and Straps



Adjustable Sandals

Adjustable Sandals, Small (869490-RT)	\$233
Adjustable Sandals, Medium (869491-RT)	\$244
Adjustable Sandals, Large (869492-RT)	\$251



Heel Support for Parted Foot Support


Heel Support for Parted Foot Support, Size 1 and 2 (Pair) (861251)	\$61
Heel Support for Parted Foot Support, Size 3 and 4 (Pair) (861252)	\$71

Additional Straps for Adjustable Sandals

Additional Straps, Small/Medium (869497)	\$60
Additional Straps, Large (869498)	\$63


Upholstered Ankle Straps

Upholstered Ankle Straps, XS (85745-XS)	\$202
Upholstered Ankle Straps, Medium (85745-M)	\$202
Upholstered Ankle Straps, XL (85745-XL)	\$202



Foot Straps

Foot Straps, Plastic Lock (Pair) (85096) (K0040)	\$42
Foot Straps, Velcro (Pair) (85102)	\$59



Neoprene Sandals

Neoprene Sandals, Small (Pair) (85696-1)	\$165
Neoprene Sandals, Medium (Pair) (85696-2)	\$165



Phone (877) 227-0097 | Fax (704) 882-0751 | sales.us@etac.com

Caribou

Tray Accessories



Tray, Plastic, with Bowls and Cover

- Tray, Plastic, with Bowls and Cover, Size 1 and 2 (Each) (862130-2) \$775
- Tray, Plastic, with Bowls and Cover, Size 3 and 4 (Each) (862130-3) \$795



Tray Armrests

- Armrest for Plastic Tray (861131) \$225



Mounting Kit

For supine use only

- Low Mount Kit for Plastic Tray (862131-KIT) \$255

Toybar for Plastic Tray

- Toybar for Plastic Tray (861135) \$262

Push Handle

- Attendant Push Handle (861108-U) \$175

Padded Tray Tops

- Square Padded Tray Top (862135-0) \$395
- Small Rounded Tray Top (862135-1) \$395
- Large Rounded Tray Top (862135-2) \$395