



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

Chilla

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 12/16/2025

Chilla Base Frame

Frame size 1 & 2					
		mm (inch), kg (lb)			
		Size 1 		Size 2 	
		Normal	Suspension	Normal	Suspension
	Overall height, to pushbrace	1200 mm (47¼ inch)	1200 mm (47¼ inch)	1200 mm (47¼ inch)	1200 mm (47¼ inch)
	Overall width	600 mm (23½ inch)	620 mm (24½ inch)	660 mm (26 inch)	680 mm (26¾ inch)
	Maximum load	43 kg (94 lb)	43 kg (94 lb)	60 kg (132 lb)	60 kg (132 lb)
	Mass device**	9.2 kg (20 lb)	9.7 kg (21 lb)	9.5 kg (20 lb)	10 kg (22 lb)
	Minimum user mass in transportation	22 kg (48 lb)	22 kg (48 lb)	22 kg (48 lb)	22 kg (48 lb)
	Maximum user mass in transportation	35 kg (77 lb)	35 kg (77 lb)	50 kg (110 lb)	50 kg (110 lb)

*380 mm (15 inch) without seat and rear wheels

**with solid wheels

***frame without rear wheels

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1309.2 Pricing and specifications subject to change without notice.



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

Chilla

Seat size 1 & 2

		mm (inch), kg (lb)	
		Size 1 	Size 2 
 Backrest recline angle		80-160°	80-160°
180° — 0°			
 Backrest height with head support plate		505-665 mm (20-26¼ inch)	620-810 mm (24½-32 inch)
 Backrest height without head support plate		365-525 mm (14½-20¾ inch)	420-610 mm (16½-24 inch)
 Back width, top of backrest (with side supports)		150-225 mm (6-8¾ inch)	170-260 mm (6¾-10¼ inch)
 Back width, top of backrest (without side supports)		315 mm (12½ inch)	365 mm (14½ inch)
 Effective seat depth		180-340 mm (7-13½ inch)	240-420 mm (9½-16½ inch)
 Effective seat width (without hip supports)		330 mm (13 inch)	380 mm (15 inch)
 Effective seat width minimum/maximum (with hip supports)		180-260 mm (7-10¼ inch)	230-310 mm (9-12¼ inch)
 Foot support to seat		130-260 mm (5-10¼ inch)	185-420 mm (7¼-16½ inch)
 Mass device*		5.3 kg (11.7 lb)	6.3 kg (13.8 lb)
 Maximum user mass		35 kg (77 lb)	50 kg (110 lb)
 Maximum user mass in transportation		35 kg (77 lb)	50 kg (110 lb)
 Minimum user mass in transportation		22 kg (48 lb)	22 kg (48 lb)
 Seat recline angle		(-5°)-(40°)	(-5°)-(40°)

*seat only (without accessories)



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

Chilla

Base Frame

Frame, Seat and Back size, must match



Chilla Size 1 Frame Solid Adjustable Suspension (E1234) (881111110)	\$2630
Chilla Size 2 Frame Solid Adjustable Suspension (E1234) (881121110)	\$2630

Adjustable Suspension Option

Omit Adjustable Suspension (\$300)

Storage Basket

Under Storage Basket (881092) N/C

Transit

Transit Fittings (881090) \$224



High-Low:xo Base and Options

High-Low:xo Indoor Base

Manual Lifts

High-Low:xo Manual Lift, Indoor, Size 1 (89010522) (K0108)	\$2196
High-Low:xo Manual Lift, Indoor, Size 2 (89020522) (K0108)	\$2356

Adaptor Mounting Plate

Required selection for High-Low:xo

Chilla Seat Adaptor Plate, For High-Low:xo (890015-15) \$875

High-Low:xo Indoor Base Accessories

Anti-Tip

Anti-Tip for Indoor High-Low:xo Indoor (890007) \$248

Push Brace

Required selection for High-Low:xo

Push Brace, Fixed Angle (890009-04)	\$194
Push Brace, Angle Adjustable (8910555)	\$465

Seat and Back

Seat and Back size must match frame size. Seat and Back are both required selections.

Seat and Back Size 1

Size 1 Seat (E2292) (X881012)	\$435
Size 1 Back (E2291) (XB881012)	\$545

Seat and Back Size 2

Size 2 Seat (E2292) (X881022)	\$455
Size 2 Back (E2291) (XB881022)	\$565

Reclining Back

Selection required

Reclining Back (E1225) (X26000) \$435

Seat and Back Set

Seat and Back Cushion

Required Selection. Pick Standard or Mesh seat and back set.

Standard Seat and Back Set

Seat and Back Set, Standard, Size 1, Pink (881311-10)	\$165
Seat and Back Set, Standard, Size 1, Blue (881311-40)	\$165
Seat and Back Set, Standard, Size 1, Red (881311-50)	\$165
Seat and Back Set, Standard, Size 1, Grey (881311-60)	\$165
Seat and Back Set, Standard, Size 2, Blue (881321-40)	\$165
Seat and Back Set, Standard, Size 2, Red (881321-50)	\$165
Seat and Back Set, Standard, Size 2, Grey (881321-60)	\$165

Mesh Seat and Back Set

Seat and Back Set, Mesh, Size 1, Blue (881312-45)	\$205
Seat and Back Set, Mesh, Size 1, Red (881312-55)	\$205
Seat and Back Set, Mesh, Size 1, Grey (881312-65)	\$205
Seat and Back Set, Mesh, Size 2, Blue (881322-45)	\$205
Seat and Back Set, Mesh, Size 2, Red (881322-55)	\$205
Seat and Back Set, Mesh, Size 2, Grey (881322-65)	\$205

Infant Insert Cushion

Available only with Size 1

Infant Cushion Pink (K0108) (881001-10)	\$375
Infant Cushion Blue (K0108) (881001-20)	\$375
Infant Cushion Black (K0108) (881001-80)	\$375



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

Chilla

Footrests

Mounting Rod and Footplate are required for complete footrest

Footrest Mounting Rod

Footrest Mounting Rod Short 7.5" (E0990) (881050-1)	\$110
Footrest Mounting Rod Long 11" (E0990) (881050-2)	\$110

Angle and Height Adjustable Footplate

Angle and Height Adj Footplate - Short Bracket (K0040) (881055)	\$105
Angle and Height Adj Footplate - Long Bracket (K0040) (881056)	\$105

Headrests

Headrest Pads

Headrest pads in this section do not require mounting hardware



Contoured Head Cushion Size 1 (E0955) (874082-1)	\$214
Contoured Head Cushion Size 2 (E0955) (874082-2)	\$214
Sculpted Head Support Size 1 (E0955) (874084-1)	\$220
Sculpted Head Support Size 2 (E0955) (874084-2)	\$220
Two-Parted Head Support (E0955) (874083)	\$475
Anatomical Head Pad, Small (E0955) (874085-1)	\$265
Anatomical Head Pad, Medium (E0955) (874085-2)	\$275
Curved Head Pad, Small (E0955) (874086)	\$220

Head Supports and Hardware

Head supports in this section require Multi-Adjustable Headrest Hardware and Bracket for Head Support Mount

Pad, Small (E0955) (9923012-80)	\$265
Pad, Medium (E0955) (9945012-80)	\$275
Curved Pad, Small (E0955) (9921012-80)	\$220

Bracket for Head Support

Bracket for Head Support Mount (K0108) (881098)	\$55
---	------

Multi-Adjustable Headrest Hardware

Multi-Adjustable Headrest Hardware (E1033) (X99XX012-80)	\$245
--	-------

Postural Supports

Hip Supports Pad and Hardware

Hardware required when Hip Support Pad is selected

Hip Support Pad Size 1 (E0956) (x8740250)	\$195
Hip Support Pad Size 2 (E0956) (x8740251)	\$225

Hip Support Hardware

Adjustable Removable Mounting Hardware Size 1 (E1033) (881023-1)	\$175
Adjustable Removable Mounting Hardware Size 2 (E1033) (881023-2)	\$175
Adjustable Removable Rotational Mounting Hardware Size 1 (E1033) (881033-1)	\$450
Adjustable Removable Rotational Mounting Hardware Size 2 (E1033) (881033-2)	\$450

Outer Knee Supports and Hardware

Hardware required when Outer Knee Pad is selected

Outer Knee Pads, Size 1 (E0953) (881043-1)	\$195
Outer Knee Pads, Size 2 (E0953) (881043-2)	\$215

Adjustable Outer Knee Hardware

Outer Knee Support Adjustable Removable Hardware (E1033) (X881043-X)	\$175
--	-------

Anti Thrust Wedge

Anti Thrust Wedge (K0108) (8720405)	\$72
-------------------------------------	------

Pommel and Hardware

Hardware required when Pommel is selected

Pommel Size 1 (E0957) (881044-1)	\$145
Pommel Size 2 (E0957) (881044-2)	\$145

Pommel Hardware

Pommel Hardware (E1028) (X881044)	\$175
-----------------------------------	-------

Trunk Supports



Side Supports Size 1 (E0956) (881030-1)	\$220
Side Supports Size 2 (E0956) (881030-2)	\$220
Side Supports Size 3 (E0956) (881030-3)	\$220

Cover

Cover for Knobs on Side Supports (8720360)	\$65
--	------



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

Chilla

Vests and Belts

Torso Vest

Torso Vest Small - S-P: 10.75" - 12.5" (E0960) (85705-S)	\$232
Torso Vest, Medium, S-P: 12.5" - 16" (85705-M)	\$232
Torso Vest Large - S-P: 16" - 18.5" (E0960) (85705-L)	\$232
Torso Vest X-Large - S-P: 18.5" - 20.5" (E0960) (85705-XL)	\$232

Cross Vest

Cross Vest X-Small - S-P: 7.75" - 12.5" (Each) (E0960) (85720-XS)	\$220
Cross Vest Small - S-P: 10.25" - 13.5" (Each) (E0960) (85720-S)	\$220
Cross Vest Medium - S-P: 13.5" - 17.75" (Each) (E0960) (85720-M)	\$220
Cross Vest Large - S-P: 17.75" - 20.75" (Each) (E0960) (85720-L)	\$220
Cross Vest X-Large - S-P: 20.75" - 24.5" (Each) (E0960) (85720-XL)	\$220

2 Point Pelvic Belt

2-Point Pelvic Belt Black - L: 13-14.5" (E0978) (85424-1)	\$98
---	------

4 Point Padded Hip Belt

4-Point Padded Hip Belts 1" X-Small - H-H: 6.75" - 12.5" (E0960) (85725-XS)	\$156
4-Point Padded Hip Belts 1" Small - H-H: 8.25" - 13.75" (E0960) (85725-S)	\$156
4-Point Padded Hip Belts 1" Medium - H-H: 9.75" - 16.25" (E0960) (85725-M)	\$156
4-Point Padded Hip Belts 1" Large - H-H: 12.5" - 17.25" (E0960) (85725-L)	\$156

Pelvic Harness

Pelvic Harness, Small - A:6.5" x B:7.5" (K0108) (89108-0)	\$284
Pelvic Harness, Medium - A:8.5" x B:8.75" (K0108) (89108-1)	\$284
Pelvic Harness, Large - A:10.5" x B:9.75" (K0108) (89108-2)	\$284

Foot Positioning

Foot Straps

Foot Straps, Plastic Lock (Pair) (85096) (K0040)	\$42
Foot Straps, Velcro (Pair) (85102) (K0040)	\$59

Neoprene Sandals

Neoprene Sandals Small (Pair) (85696-1)	\$165
Neoprene Sandals Medium (Pair) (85696-2)	\$165

Ankle Straps

Ankle Straps X-Small - C: 6.25" - 9" (K0108) (85745-XS)	\$202
Ankle Straps Medium - C: 9" - 12.25" (K0108) (85745-M)	\$202
Ankle Straps X-Large - C: 11.5" - 15" (K0108) (85745-XL)	\$202

Accessories

Canopy

For more information on canopy options, please see etac.com

Canopy, Size 1 (K0108) (881080-10)	N/C
Canopy, Size 2 (K0108) (881080-20)	N/C
Premium Canopy, Size 1 (K0108) (881080-30)	\$395
Premium Canopy, Size 2 (K0108) (881080-40)	\$395



Tray / Grip Bar and Hardware

Trays and Grip Bars use the same mounting hardware. Both may be selected; however, only one may be used at a time. Mounting hardware is required for both options.

Tray (E0950) (881071)	\$112
-----------------------	-------

Additional Accessories

Bracket for Communication Aid (K0108) (881099)	\$176
Rain Cover, Size 1 (881093-1)	\$125
<i>Requires Canopy</i>	
Rain Cover, Size 2 (881093-2)	\$125
<i>Requires Canopy</i>	
Hand Brake Kit (881020)	\$375
<i>Available only with the Suspension Frame</i>	
Backpack Size 1 (88043184-1)	\$125
Backpack Size 2 (88043184-2)	\$125



Grip Bar Size 1 (K0108) (881070-1)	\$95
Grip Bar Size 2 (K0108) (881070-2)	\$95

Tray and Grip Bar Hardware

Mounting Hardware for Tray or Grip Bar (E0128) (881060)	\$176
---	-------